**Monoclonal Mouse
Anti-Human PTEN
Clone 6H2.1**

| | |
|---|--|
| English | Code M3627 |
| Intended use | For In Vitro diagnostic use. This antibody is intended for laboratory use to identify qualitatively by light microscopy PTEN expressing cells in normal and neoplastic tissues using immunohistochemical (IHC) test methods. Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. |
| Synonym | MMAC (mutated in multiple advanced cancers), TEP1(TGF- β -regulated and epithelial cell-enriched phosphatase) |
| Summary and explanation | <p>PTEN (phosphatase and tensin homolog deleted on chromosome ten) is a tumor suppressor gene that is mutated in a wide range of cancers.¹⁻⁴ PTEN is located on chromosome sub-band 10q23.3 and encodes a 403 amino acid protein which acts as a dual-specific protein and phospholipid phosphatase.^{2,3,5,6} As a lipid phosphatase, PTEN regulates phosphatidylinositol-3'-kinase by dephosphorylating the D3 position of phosphatidylinositol (3,4,5)-triphosphate and phosphatidylinositol (3,4)-bisphosphate substrates, thereby antagonizing signal transduction downstream of phosphatidylinositol-(PI-3) kinase.⁷ Loss-of-function mutations in the PTEN gene lead to constitutive activation of multiple signaling pathways including the PI3K/Akt pathway which affects cell proliferation, apoptosis and migration.^{8,9} PTEN has also been shown to control p53 protein levels and transcriptional activity.¹⁰</p> <p>Refer to Dako's <i>General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.</p> |
| Reagent provided | Code M3627 Anti-PTEN is provided in liquid form as tissue culture supernatant (containing fetal bovine serum) dialyzed against 0.05 mol/L Tris-HCl, pH 7.2, and 0.015 mol/L sodium azide. Contains stabilizing protein. Clone: 6H2.1 ¹¹ Isotype: IgG _{2a} , kappa Mouse IgG concentration mg/L: See label on vial. M3627 may be used at a dilution of 1:100 when performing IHC using the EnVision™+ detection system. These are guidelines only. Optimal antibody concentrations may vary depending on specimen and preparation method, and should be determined individually in each laboratory. |
| Immunogen | Full-length recombinant human PTEN protein |
| Specificity | Monoclonal mouse anti-PTEN clone 6H2.1, identifies a single band at the 55 kD predicted molecular weight for PTEN in Western blots of MCF-7, T-47D and MDA-MB-435S PTEN homozygous cell lines. PTEN null cell lines BT-549 and MDA-MB-468 cell lines were unreactive with clone 6H2.1 by Western blot analysis. Immunoblots of MCF-7 cells induced to express PTEN demonstrated an increased band intensity at the expected size while a weak band was observed on immunoblots of ZR-75-1 cell lysates containing a hemizygous deletion of PTEN and a mis-sense mutation in the remaining allele. ¹¹ Immunohistochemical staining of paraffin-embedded cell blocks of Ishikawa 3-H-12 cells transfected with PTEN stained positively whereas PTEN deficient 3-H-12 cell block sections were unreactive. ¹² |
| Materials required, but not supplied | Refer to Dako's <i>General Instructions for Immunohistochemical Staining</i> and/or the detection system instructions. Suggested diluent for IHC procedures: Antibody Diluent with Background Reducing Components (code S3022) The following negative control is recommended for IHC procedures: Mouse IgG _{2a} (code X0943) |
| Precautions | <ol style="list-style-type: none">1. For professional users.2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.3. As with any product derived from biological sources, proper handling procedures should be used.4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.5. Unused reagents should be disposed of according to local, State, and Federal regulations. |

Storage

Store at 2–8°C.

Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation*Paraffin Sections*

Anti-PTEN can be used on formalin-fixed, paraffin-embedded tissue sections. Pretreatment of tissue with proteolytic enzymes is not recommended.

The deparaffinized tissue sections must be treated with heat prior to the IHC staining procedure. Heat-induced epitope retrieval (HIER) involves immersion of tissue sections in a pre-heated buffer solution and maintaining heat in a water bath (95–99°C). Alternative heat sources may be used for HIER upon validation against the recommended procedure. Use a 40-minute heating protocol for HIER performed at 95–99°C; after thermal treatment, allow the jar with buffer and slides to cool for 20 minutes at room temperature. Rinse well with buffer or deionized water following HIER. For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended. Target Retrieval Solution pH 9.0 (code S2368) or 10x Concentrate (code S2367) is recommended.

Cryostat Sections And Cell Smears

Anti-PTEN can be used for labeling acetone-fixed cryostat sections or fixed cell smears.

Staining procedure

Follow the procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-PTEN is cytoplasmic and/or nuclear.

Performance characteristics

*Normal Tissues*¹³

| Tissue Type (# tested) | Positively Staining Tissue Elements |
|-------------------------------|--|
| Adrenal (3) | 2/3 zona glomerulosa; 2/3 vascular endothelium |
| Bone Marrow (3) | 1/3 macrophages |
| Breast (3) | 2/3 ductal epithelial cells; 1/3 lobular epithelial cells; 1/3 vascular endothelium |
| Brain/Cerebellum (3) | 3/3 molecular and granular layer, neuropil |
| Brain/Cerebrum (3) | 3/3 neurons, astrocytes; 2/3 neuropil |
| Cervix (3) | 3/3 squamous mucosa; 1/3 vascular smooth muscle, vascular endothelium |
| Colon (2) | 1/2 glandular epithelium, ganglion cells, macrophages; 1/2 germinal center cells; 2/2 nerve, vascular endothelium |
| Esophagus (3) | 1/3 squamous mucosa; 2/3 nerve; 1/3 vascular endothelium |
| Heart (3) | 3/3 vascular endothelium; 1/3 vascular smooth muscle; 1/3 nerve |
| Kidney (3) | 2/3 glomeruli, tubule epithelium, endothelial cells; 1/3 mesangial cells, vessels; 1/3 lining cells, interstitial cells |
| Liver (3) | 3/3 hepatocytes; 1/3 bile ducts |
| Lung (3) | 2/3 bronchial epithelium; 3/3 vascular endothelium; 1/3 vascular smooth muscle, macrophages |
| Mesothelial Cells (3) | 3/3 vascular endothelial cells |
| Nerve, Peripheral (2) | 2/2 nerve cells |
| Ovary (3) | 3/3 stroma; 1/3 surface epithelium |
| Pancreas (3) | 2/3 ductal epithelium, nerve; 1/3 acinar epithelium; 3/3 islets of Langerhans.; 2/3 vascular endothelium |
| Parathyroid (3) | 2/3 chief cells |
| Pituitary (3) | 3/3 pituicytes, neuropil; 2/3 vascular endothelium; 1/3 nerve, vascular smooth muscle |
| Prostate (3) | 3/3 glandular epithelium, fibromuscular stroma; 1/3 vascular endothelium |
| Salivary gland (3) | 1/3 acini and duct epithelium; 2/3 nerve |
| Skeletal muscle (3) | 3/3 vascular endothelium; 1/3 vascular smooth muscle |
| Skin (3) | 3/3 epidermis, hair follicles, sebaceous glands; 1/3 apocrine glands |
| Small intestine (3) | 1/3 glandular epithelium, smooth muscle, macrophages, germinal center cells; 1/3 surface epithelium; 3/3 nerve; 3/3 vascular endothelium |
| Spleen (3) | 1/3 macrophages (rare); 1/3 littoral cells |
| Stomach (3) | 3/3 glandular epithelium; 2/3 vascular endothelium; 1/3 vascular smooth muscle, macrophages; 1/3 nerve; 1/3 germinal center cells |
| Testis (3) | 3/3 seminiferous tubules; 1/3 Leydig cells, vascular endothelium |
| Thymus (3) | 3/3 thymic epithelium |
| Thyroid (2) | 2/2 negative |
| Tonsil (3) | 3/3 squamous epithelium; 2/3 germinal center cells |
| Uterus (3) | 3/3 endometrial glands, endometrial stroma, myometrium |

Abnormal Tissues^{11, 12, 14-16}

| Tumor Type (# tested) | Positively Staining Tumors |
|---|-----------------------------------|
| Breast carcinoma, sporadic primary carcinoma (33) | 28/33 |
| Colorectal carcinoma (90) | 82/90 |
| MSI + HNPCC (45) | 40/45 |
| MSI+ sporadic colorectal carcinomas (22) | 19/22 |
| MSI- sporadic colorectal carcinomas (23) | 23/23 |
| Endometrial cancer (72) | 50/72 |
| Pancreatic endocrine tumors, sporadic (24) | 23/24 |
| Thyroid carcinoma, sporadic nonmedullary (91) | 83/91 |
| Follicular thyroid carcinoma (28) | 27/28 |
| Papillary thyroid carcinoma (37) | 35/37 |
| Undifferentiated thyroid carcinoma (26) | 21/26 |

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Synonyme

MMAC (mutated in multiple advanced cancers), TEP1 (TGF- β -regulated and epithelial cell-enriched phosphatase)

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References
Bibliographie
Literaturangaben

1. Li J, Yen C, Liaw D, Podsypanina K, Bose S, Wang SI, Puc J, Miliareisis C, Rodgers L, McCombie R, Bigner SH, Giovanella BC, Ittmann M, Tycko B, Hibshoosh H, Wigler MH, Parsons R. PTEN, a putative protein tyrosine phosphatase gene mutated in human brain, breast, and prostate cancer. *Science* 1997; 275(5308):1943-1947
2. Steck PA, Pershouse MA, Jasser SA, Yung WK, Lin H, Ligon AH, Langford LA, Baumgard ML, Hattier T, Davis T, Frye C, Hu R, Swedlund B, Teng DH, Tavtigian SV. Identification of a candidate tumour suppressor gene, MMAC1, at chromosome 10q23.3 that is mutated in multiple advanced cancers. *Nat Genet* 1997; 15(4):356-362
3. Li DM, Sun H. TEP1, encoded by a candidate tumor suppressor locus, is a novel protein tyrosine phosphatase regulated by transforming growth factor beta. *Cancer Res* 1997; 57(11):2124-2129
4. Teng DH, Hu R, Lin H, Davis T, Iliev D, Frye C, Swedlund B, Hansen KL, Vinson VL, Gumpfer KL, Ellis L, El-Naggar A, Frazier M, Jasser S, Langford LA, Lee J, Mills GB, Pershouse MA, Pollack RE, Tornos C, Troncoso P, Yung WK, Fujii G, Berson A, Steck PA, et al. MMAC1/PTEN mutations in primary tumor specimens and tumor cell lines. *Cancer Res* 1997; 57(23):5221-5225
5. Maehama T, Dixon JE. The tumor suppressor, PTEN/MMAC1, dephosphorylates the lipid second messenger, phosphatidylinositol 3,4,5-trisphosphate. *J Biol Chem* 1998; 273(22):13375-13378
6. Myers MP, Stolarov JP, Eng C, Li J, Wang SI, Wigler MH, Parsons R, Tonks NK. P-TEN, the tumor suppressor from human chromosome 10q23, is a dual-specificity phosphatase. *Proc Natl Acad Sci USA* 1997; 94(17):9052-9057
7. Cantley LC, Neel BG. New insights into tumor suppression: PTEN suppresses tumor formation by restraining the phosphoinositide 3-kinase/AKT pathway. *Proc Natl Acad Sci USA* 1999; 96(8):4240-4245. Review
8. Downward J. Mechanisms and consequences of activation of protein kinase B/Akt. *Curr Opin Cell Biol* 1998; 10(2):262-267
9. Hill MM, Hemmings BA. Inhibition of protein kinase B/Akt: implications for cancer therapy. *Pharmacol Ther* 2002; 93(2-3):243-251
10. Freeman DJ, Li AG, Wei G, Li HH, Kertesz N, Lesche R, Whale AD, Martinez-Diaz H, Rozengurt N, Cardiff RD, Liu X, Wu H. PTEN tumor suppressor regulates p53 protein levels and activity through phosphatase-dependent and -independent mechanisms. *Cancer Cell* 2003; 3(2):117-130
11. Perren A, Weng LP, Boag AH, Ziebold U, Thakore K, Dahia PL, Komminoth P, Lees JA, Mulligan LM, Mutter GL, Eng C. Immunohistochemical evidence of loss of PTEN expression in primary ductal adenocarcinomas of the breast. *Am J Pathol* 1999; 155(4):1253-1260
12. Pallares J, Bussaglia E, Martinez-Guitarte JL, Dolcet X, Llobet D, Rue M, Sanchez-Verde L, Palacios J, Prat J, Matias-Guiu X. Immunohistochemical analysis of PTEN in endometrial carcinoma: a tissue microarray study with a comparison of four commercial antibodies in correlation with molecular abnormalities. *Mod Pathol* 2005; 18(5):719-727
13. IHC003 Report On File
14. Zhou XP, Loukola A, Salovaara R, Nystrom-Lahti M, Peltomaki P, de la Chapelle A, Aaltonen LA, Eng C. PTEN mutational spectra, expression levels, and subcellular localization in microsatellite stable and unstable colorectal cancers. *Am J Pathol* 2002; 161(2):439-447
15. Perren A, Komminoth P, Saremaslani P, Matter C, Feurer S, Lees JA, Heitz PU, Eng C. Mutation and expression analyses reveal differential subcellular compartmentalization of PTEN in endocrine pancreatic tumors compared to normal islet cells. *Am J Pathol* 2000; 157(4):1097-1103
16. Gimm O, Perren A, Weng LP, Marsh DJ, Yeh JJ, Ziebold U, Gil E, Hinze R, Delbridge L, Lees JA, Mutter GL, Robinson BG, Komminoth P, Dralle H, Eng C. Differential nuclear and cytoplasmic expression of PTEN in normal thyroid tissue, and benign and malignant epithelial thyroid tumors. *Am J Pathol* 2000; 156(5):1693-1700

| | | |
|---|--|---|
| REF Catalogue number Référence du catalogue Bestellnummer |  Temperature limitation Limites de température Zulässiger Temperaturbereich |  Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i> |
|  Manufacturer Fabricant Hersteller | LOT Batch code Code du lot Chargenbezeichnung |  Use by Utiliser jusque Verwendbar bis |
| EC REP Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU | IVD In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum | |

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA
Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

EC REP
Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark
Tel +45 4485 9500
Fax +45 4485 9595
www.dako.com

PT0039/Rev C