**Monoclonal Mouse
Anti-Human
Plasma Cell**
Clone VS38c
Code No./ Code/ Code-Nr. M 7077
Edition/ Ausgabe 16.12.02

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human Plasma Cell, Clone VS38c, is intended for use in immunocytochemistry. The antibody labels normal and neoplastic plasma cells and may be a useful tool in the identification of myeloma and plasmacytoma in bone marrow or other tissues. Additionally it can assist in distinguishing between lymphoplasmacytoid lymphoma and lymphocytic and follicular lymphoma (1). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.</p>
Synonym for antigen	p63 protein (1).
Introduction	<p>Human p63 is a 63 kDa non-glycosylated, reversibly palmitoylated type II transmembrane protein found in plasma cells and epithelium. The p63 protein has a currently unknown function, but because of its homology to both rat and swine proteins, its abundance in secretory cells, and its localization to the rough endoplasmic reticulum, a conserved role in protein processing or secretion is suggested (1).</p>
Reagent provided	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN₃.</p> <p><u>Clone:</u> VS38c (2). <u>Isotype:</u> IgG1, kappa.</p> <p><u>Mouse IgG concentration:</u> see label on vial.</p>
Immunogen	Membrane preparation of MCF-7 AdR cells (breast carcinoma cell line) (2).
Specificity	<p>Anti-Human Plasma Cell, VS38c, was clustered as anti-p63 protein at the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens held in Kobe in 1996. The epitope was localized to the C-terminal 419 residues (1).</p> <p>In Western blotting of cell lysates from plasma and epithelial cell lines under reducing conditions, the antibody labels a band of 64 kDa, corresponding to the p63 protein (1, 2).</p>
Precautions	<ol style="list-style-type: none">1. For professional users.2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	<p>Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.</p>
Specimen preparation	<p><u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with Dako Target Retrieval Solution, High pH, code No. S 3308, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. However, Dako Target Retrieval Solution, code No. S 1700, or 10 mmol/L citrate buffer, pH 6.0 and pre-treatment of tissues with proteinase K were found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.</p> <p><u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections (1, 2).</p>
Staining procedure	<p><u>Dilution:</u> Monoclonal Mouse Anti-Human Plasma Cell, code No. M 7077, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, High pH, code No. S 3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen.</p> <p><u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.</p>

Product-specific limitations	The p63 protein is not plasma cell specific and the antibody has been shown to label a varying proportion of both benign naevi and malignant melanomas in a diffuse cytoplasmic pattern (3). Additionally it labelled 36/50 neuroendocrine tumours of various organs (4). In frozen sections weak endothelial labelling can occasionally be detected (2).
Performance characteristics	Cells labelled by the antibody show staining of the cell membrane and/or cytoplasm. <u>Normal tissues:</u> In frozen and formalin-fixed tissues, plasma cells are strongly labelled, while peripheral blood lymphocytes are negative although some weak labelling of monocytes and neutrophils can occasionally be detected in frozen sections. Epithelial cell labelling observed in frozen sections is considerably reduced or absent in formalin-fixed tissue. Unidentified, strongly labelled spindle cells within the stroma of some tissues were also observed (2). <u>Abnormal tissues:</u> In formalin-fixed tissue, the antibody labelled 7/7 myelomas, 4/4 plasmacytomas, 7/7 lymphoplasmacytoid lymphomas, and weakly 12/38 B-cell large cell non-Hodgkin's lymphomas. None of 12 chronic lymphocytic leukaemias, 6 centroblastic centrocytic lymphomas, 4 hairy cell leukaemias, and 9 T-cell non-Hodgkin's lymphomas were labeled (2).

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