**FLEX**  
**Monoclonal Mouse**  
**Anti-Cytomegalovirus**  
 Clones CCH2 + DDG9  
**Ready-to-Use**  
 (Dako Autostainer/Autostainer Plus)

**Code IS752**

ENGLISH	
<b>Intended use</b>	For in vitro diagnostic use. FLEX Monoclonal Mouse Anti-Human Cytomegalovirus, Clone CCH2 + DDG9, Ready-to-Use (Dako Autostainer/Autostainer Plus), is intended for use in immunohistochemistry together with Dako Autostainer/Autostainer Plus instruments. This antibody is useful for the identification of human cytomegalovirus (HCMV)-infected human tissues (1, 2). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
<b>Summary and explanation</b>	Human cytomegalovirus is a ubiquitous herpes virus. The primary infection usually occurs in childhood where it most often is asymptomatic. After initial exposure, the virus establishes a life-long infection. In neonates, and in immunocompromised individuals, such as AIDS patients and transplant recipients, HCMV can cause severe life-threatening infections, and it is the most common cause of viral birth defects (5-6). Antibody CCH2 (1) reacts with an early nuclear protein identical with the non-structural DNA-binding protein p52 (7). Antibody DDG9 reacts with an immediate early nuclear protein of about 76 kDa (3). For both antibodies, the reactivity persists also at later stages during the HCMV infection where the localization is less distinctly nuclear and appears to be in the cytoplasm. Laser confocal microscopy, however, shows that the reaction is limited to the nuclear membrane (3). The antibodies do not cross-react with adenovirus, herpes simplex virus, and varicella zoster virus (3). Refer to Dako's <i>General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.
<b>Reagent provided</b>	Ready-to-use mixture of two monoclonal mouse antibodies provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, 0.015 mol/L NaN <sub>3</sub> , pH 7.2. <u>Clones</u> : CCH2 + DDG9 (3). <u>Isotypes</u> : IgG1, kappa (CCH2) and IgG2a, kappa (DDG9).
<b>Immunogen</b>	Crude cell lysates from cytomegalovirus strain Ad169-infected human embryonic fibroblasts (3).
<b>Specificity</b>	Antibody CCH2 reacts with an early nuclear protein identical with the non-structural DNA-binding protein p52 (7). More than 200 human CMV isolates have been successfully typed by this antibody which, therefore, seems to recognize an epitope which is highly conserved among HCMV strains (3, 4). Antibody DDG9 reacts with an immediate early nuclear protein of about 76 kDa (3). For both antibodies, the reactivity persists also at later stages during HCMV infection where localization is less distinctly nuclear and appears to be in the cytoplasm. Laser confocal microscopy, however, shows that the reaction is limited to the nuclear membrane (3). The antibodies do not cross-react with adenovirus, herpes simplex virus, and varicella zoster virus (3).
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
<b>Specimen preparation including materials required but not supplied</b>	The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm. Pre-treatment with heat-induced epitope retrieval (HIER) is required using Dako PT Link (Code PT100/PT101). For details, please refer to the PT Link User Guide. Optimal results are obtained by pretreating tissues using EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). <u>Paraffin-embedded sections</u> : Pre-treatment of formalin-fixed, paraffin-embedded tissue sections is recommended using the 3-in-1 specimen preparation procedure for Dako PT Link. Follow the pre-treatment procedure outlined in the package insert for EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). Note: After staining the sections must be dehydrated, cleared and mounted using permanent mounting medium. <u>Deparaffinized sections</u> : Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link and following the same procedure as described for paraffin-embedded sections. After staining the slides should be mounted using aqueous or permanent mounting medium. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended.
<b>Staining procedure including materials required but not supplied</b>	The recommended visualization system is EnVision™ FLEX, High pH, (Dako Autostainer/Autostainer Plus) (Code K8010). The staining steps and incubation times are pre-programmed into the software of Dako Autostainer/Autostainer Plus instruments, using the following protocols: Template protocol: FLEXRTU2 (200 µL dispense volume) or FLEXRTU3 (300 µL dispense volume) Autoprogram: CMV (without counterstaining) or CMVH (with counterstaining) The Auxiliary step should be set to "rinse buffer" in staining runs with ≤10 slides. For staining runs with >10 slides the Auxiliary step should be set to "none". This ascertains comparable wash times. All incubation steps should be performed at room temperature. For details, please refer to the Operator's Manual for the dedicated instrument. If the protocols are not available on the used Dako Autostainer instrument, please contact Dako Technical Services. Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory. If the evaluating pathologist should desire a different staining intensity, a Dako Application Specialist/Technical Service Specialist can be contacted for information on re-programming of the protocol. Verify that the performance of the adjusted protocol is still valid by evaluating that the staining pattern is identical to the staining pattern described in "Performance characteristics". Counterstaining in hematoxylin is recommended using EnVision™ FLEX Hematoxylin, (Dako Autostainer/Autostainer Plus) (Code K8018). Non-aqueous,

(117579-004)

IS752/EFG/MNI/2009.12.04 p. 1/4

Dako Denmark A/S | Produktionsvej 42 | DK-2600 Glostrup | Denmark | Tel. +45 44 85 95 00 | Fax +45 44 85 95 95 | CVR No. 33 21 13 17

permanent mounting medium is recommended.

Positive and negative controls should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include HCMV-infected tissue and the cells/structures should display reaction patterns as described for this tissue in "Performance characteristics" in all positive specimens. The recommended negative control reagent is FLEX Negative Control, Mouse, (Dako Autostainer/Autostainer Plus) (Code IS750).

**Staining interpretation**

The cellular staining pattern is a strong intranuclear reaction. Cytoplasmic staining may be observed.

**Performance characteristics**

Cells labeled by the antibody display a nuclear staining pattern early in the HCMV infection. At a later stage, an apparent cytoplasmic staining pattern might be observed (3).

Normal tissues: Screening of normal human tissues from donors without HCMV infection with the CCH2 antibody showed no reaction with the heart, lung, liver, kidney, spleen, lymph node, tonsil, bone marrow, skeletal muscle, brain, thyroid and parathyroid glands (1).

Abnormal tissues: In HCMV-infected tissues, the CCH2 antibody labeled cells in the lung, liver, kidney, lymph node, brain, thyroid and parathyroid glands (1). The CCH2 and DDG9 antibodies together labeled HCMV-infected tissue of esophagus, stomach, small intestine, colon, gall bladder, lung, adrenal gland, ovary and neural tissue (2).

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
- Niedobitek G, Finn T, Herbst H, Gerdes J, Grillner L, Landqvist M, et al. Detection of cytomegalovirus by in situ hybridization and immunohistochemistry using new monoclonal antibody CCH2: a comparison of methods. J Clin Pathol 1988; 41:1005-9.
- Saiz E, Lubin J, Robinson MJ. The modified Steiner stain: a new use for an old stain? Staining cytomegalovirus-infected cells in gastrointestinal biopsies. Histochem J 1998; 30:549-52.
- Wirgart BZ. Development of rapid methods for early diagnosis of cytomegalovirus infections [dissertation]. Stockholm: Karolinska Hospital; 1991.
- Wirgart BZ, Landqvist M, Hökeberg I, Eriksson B-M, Olding-Stenkvist E, Grillner L. Early detection of cytomegalovirus in cell culture by a new monoclonal antibody, CCH2. J Virol Methods 1990; 27:211-20.
- Doniger J, Muralidhar S, Rosenthal LJ. Human cytomegalovirus and human herpes virus 6 genes that transform and transactivate. Clin Microbiol Rev 1999; 12:367-82.
- Jarvis MA, Nelson JA. Human cytomegalovirus persistence and latency in endothelial cells and macrophages. Curr Opin Microbiol 2002; 5:403-7.
- Plachter B, Nordin M, Wirgart BZ, Mach M, Stein H, Grillner L, et al. The DNA-binding protein P52 of human cytomegalovirus reacts with monoclonal antibody CCH2 and associates with the nuclear membrane at late times after infection. Virus Res 1992;24: 265-76.

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