**Monoclonal Mouse
Anti-Human Podoplanin
Clone D2-40**

Code M3619

Intended use

For In Vitro Diagnostic Use

Monoclonal Mouse Anti-Human Podoplanin, Clone D2-40 is intended for use in immunohistochemistry. This antibody labels the lymphatic endothelium marker podoplanin in normal and neoplastic tissues. The antibody may be useful in the identification of lymphatic invasion of a variety of cancers (1, 2). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and explanation

Podoplanin is a ~38 kDa O-linked transmembrane sialoglycoprotein which is expressed in the endothelium of lymphatic capillaries, but not in the blood vasculature (3). Besides the expression in lymphatic endothelium, podoplanin is also found in a variety of other tissues, including mesothelial cells, reticular cells, follicular dendritic cells, ovarian and testicular germ cells (4). Anti-Podoplanin has also been demonstrated to react with lymphatic endothelium and is unreactive with vascular endothelium (1, 2). In neoplastic tissue, immunostaining of lymphatic endothelium by Anti-Podoplanin has been shown to be useful in identifying lymphatic invasion of primary tumors (1, 2). The use of calretinin and podoplanin markers in conjunction has been shown to improve differentiation of mesothelioma versus metastatic adenocarcinoma (5). Podoplanin has also been shown to be a useful marker for the differential identification of seminomas and follicular dendritic cell tumors (6).

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Monoclonal Mouse antibody provided in liquid form as tissue culture supernatant in 0.05 mol/L Tris-HCl, pH 7.2 and 0.015 mol/L sodium azide. This product contains stabilizing protein.

Clone: D2-40 (7) Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

M3619 may be used at a dilution 1:100–1:200 when performing IHC using the EnVision+™ detection system. These are guidelines only. Optimal antibody concentrations may vary depending on specimen and preparation method, and should be determined individually in each laboratory.

Immunogen

Dysgerminoma tissue (7).

Specificity

Clone D2-40 was tested by Western blotting of cell lysates from a variety of tumor cell lines. D2-40 identified a 40 kD band in immunoblots of HEY ovarian epithelial carcinoma cells and a slightly higher molecular weight form in blots of MGH-U3 cells. D2-40 immunoreactivity has also been observed in Western blots of SK-OV-3 and Caov-3 ovarian epithelial carcinoma cells, PA-1PA-1 ovarian teratocarcinoma cells, Tera-1 and Tera-2 testicular teratocarcinoma cells, 833 embryonal carcinoma cells, RT-4 bladder transitional carcinoma cells and U-2 OS osteosarcoma cells. No reactivity was observed in Western blots of other tumor cell lines derived from a variety of cancers (7).

The effect of enzyme digestion on the immunoreactivity of D2-40 was tested using affinity purified antigen from HEY cells and intact HEY cells. The immunoreactivity of D2-40 with Western blots of purified antigen and with cells by flow cytometry was unaffected by the removal of sialic acid residues through neuraminidase treatment of the cells or protein. Treatment of protein or cells with O-sialoglycoprotease destroyed the immunoreactivity of D2-40 in both Western blotting and flow cytometry, demonstrating that the antigen recognized by D2-40 is an O-linked glycoprotein (7).

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions. Suggested diluent(s) for IHC procedures:

Antibody Diluent (code S0809)

The following negative control is recommended for IHC procedures:

Mouse IgG₁ (code X0931)

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin sections

Anti-Podoplanin can be used on formalin-fixed, paraffin-embedded tissue sections. Pretreatment of tissue with proteolytic enzymes is not recommended.

The deparaffinized tissue sections must be treated with heat prior to the IHC staining procedure. Heat-induced epitope retrieval (HIER) involves immersion of tissue sections in a pre-heated buffer solution and maintaining heat in a water bath or steamer (95–99 °C). Use a 20-minute heating protocol for HIER performed at 95–99 °C; after thermal treatment, allow the jar with buffer and slides to cool for 20 minutes at room temperature. Rinse well with buffer or deionized water following HIER. For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended. Target Retrieval Solution (code S1700) or 10x Concentrate (code S1699) is recommended.

Cryostat sections and cell smears

Anti-Podoplanin can be used for labeling acetone-fixed cryostat sections or fixed cell smears.

Staining procedure

Follow the recommended procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for Anti-Podoplanin is cytoplasmic and sometimes membranous.

Performance characteristics

Normal tissues⁽⁸⁾

Anti-Podoplanin strongly immunostains lymphatic endothelium, ganglion cells and nerve in many normal tissue types. Vascular endothelial cells were negative in all tissues (8).

Tissue Type (# tested)	Positively Staining Tissue Elements
Adrenal (3)	2/3 cytoplasm of capsular fibroblasts; 1/3 cytoplasm of interstitial fibroblasts; 1/3 cytoplasm of pericapsular lymphocytes
Bone marrow (3)	3/3 negative
Brain, cerebellum (3)	3/3 cytoplasm of neuropil
Brain, cerebrum (3)	3/3 nuclei of rare astrocytes
Breast (3)	3/3 cytoplasm of myoepithelial cells
Cervix (3)	3/3 membrane and/or cytoplasm of basilar epithelial cells, 3/3 cytoplasm of stromal fibroblasts
Colon (3)	2/3 lymphocytes
Esophagus (3)	1/3 cytoplasm of ganglion cells in muscle; 1/3 cytoplasm of basilar epithelium
Heart (3)	2/3 cytoplasm of myocytes
Kidney (3)	1/3 cytoplasm of pericapsular fibroblasts
Liver (3)	3/3 negative
Lung (3)	2/3 cytoplasm of bronchial epithelial cells, 3/3 cytoplasmic edge of alveolar cells
Mesothelial cells (3)	1/3 cytoplasm of mesothelial cells
Ovary (3)	3/3 cytoplasm of ovarian stromal cells; 1/3 cytoplasm of surface epithelium of inclusion cysts
Pancreas (3)	2/3 cytoplasm of periductal fibroblasts; 1/3 cytoplasm of ductal myoepithelial cells
Parathyroid (3)	3/3 negative
Peripheral nerve (3)	2/3 cytoplasm of perineurium
Pituitary (3)	1/3 cytoplasm of pericapsular fibroblasts in neurohypophysis; 1/3 cytoplasm of neuropil in neurohypophysis
Prostate(3)	3/3 cytoplasm of basal cells of ducts; 3/3 cytoplasm of stromal fibroblasts
Salivary gland (3)	3/3 cytoplasm of ductal and acinar myoepithelial cells
Skeletal muscle (3)	3/3 cytoplasm of perimysium
Skin (3)	3/3 cytoplasm of myoepithelial cells of adnexal structures
Small intestine (3)	3/3 cytoplasm of myenteric plexus
Spleen (3)	3/3 negative
Stomach (3)	3/3 cytoplasm of periglandular fibroblasts; 1/3 cytoplasm of myenteric plexus; 1/3 cytoplasm and membrane of germinal center lymphocytes
Testes (3/3)	3/3 cytoplasm of peritubular fibroblasts
Thymus (3)	3/3 cytoplasm and membrane of medullary lymphocytes
Thyroid (3)	3/3 negative
Tonsil (3)	3/3 cytoplasm and membrane of basilar epithelial cells and germinal center lymphocytes
Uterus (3)	3/3 cytoplasm of myometrium; 1/3 cytoplasm of endometrial glands

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Bibliographie

Literaturangaben

- Kahn HJ, Marks A. A new monoclonal antibody, D2-40, for detection of lymphatic invasion in primary tumors. Lab Invest 2002;82: 1255-7.
- Kahn HJ, Bailey D, Marks A. Monoclonal antibody D2-40, a new marker of lymphatic endothelium, reacts with Kaposi's sarcoma and a subset of angiosarcomas. Mod Pathol 2002; 15:434-40.
- Breiteneder-Geleff S, Soleiman A, Kowalski H, Horvat R, Amann G, Kriehuber E, et al. Angiosarcomas express mixed endothelial phenotypes of blood and lymphatic capillaries: podoplanin as a specific marker for lymphatic endothelium. Am J Pathol 1999; 154:385-94.
- Kalof AN and Cooper K. D2-40 Immunohistochemistry – So Far! Adv Anat Pathol 2009; 16:62-4.
- Chu AY, Litzky LA, Pasha TL, Acs G, Zhang PJ. Utility of D2-40, a novel mesothelial marker, in the diagnosis of malignant mesothelioma, Mod Pathol 2005, 18:105-110.
- Yu H, Gibson JA, Pinkus GS, Hornick JL. Podoplanin (D2-40) is a novel marker for follicular dendritic cell tumors. Am J Clin Pathol 2007; 128:776-82.
- Marks A, Sutherland DR, Bailey D, Iglesias J, Law J, Lei M, et al. Characterization and distribution of an oncofetal antigen (M2A antigen) expressed on testicular germ cell tumours. Br J Cancer 1999; 80:569-78.
- IHC003 Report On File

REF Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten
 Manufacturer Fabricant Hersteller	LOT Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis
EC REP Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	IVD In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA
Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

EC REP
Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark
Tel +45 4485 9500
Fax +45 4485 9595
www.dako.com

PT0039/Rev C

Edition 07/12