Monoclonal Mouse Anti-Human Desmin
Clone D33
Code No./ Code/ Code-Nr. M 0760
Edition/ Ausgabe 19.12.02

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Desmin, Clone D33, is intended for use in immunocytochemistry. The antibody labels smooth and striated muscle cells as well as mesothelial cells, and is a useful tool for the identification of rhabdomyosarcomas (1, 2), leiomyomas, (3, 4) and mesotheliomas (5). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Introduction	Desmin belongs to the class III of intermediate filaments, constituting part of the cytoskeleton, and is the characteristic intermediate filament of all three types of muscle cells (skeletal, cardiac and smooth muscle). Desmin is a 53-kDa protein encoded by nine exons of a gene located at 2q35. Desmin forms a cytoskeletal network across the muscle fibre bordering at the plasma and nuclear membrane and is particularly localized to the subplasmalemmal region and the Z-band (6).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> D33. <u>IsoType:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial.
Immunogen	Desmin purified from human muscle.
Specificity	In Western blotting of cytoskeletal preparations of MCF-7 cells, SCC-4 cells RT-112 cells, skin epidermis, primary culture of renal carcinoma cells, leiomyoma, fetal heart, and crude extract from human white matter of brain, the antibody labels only a ~53 kDa band in the leiomyoma and fetal heart preparations, proving the reactivity with desmin, and demonstrating the absence of reactivity with other types of intermediate filaments tested (7). As demonstrated by immunocytochemistry, the antibody cross-reacts with the desmin-equivalent protein in mouse and rat.
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin or B5-fixative (1). Pre-treatment of tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, or 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections (7) and fixed cell smears (1, 5).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human Desmin, code No. M 0760, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human colon and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval solution, High pH, code No. S 3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.
Performance characteristics	Cells labelled by the antibody display a cytoplasmic staining pattern. The staining may show a fibrillar aspect (7). <u>Normal tissues:</u> The antibody labels vascular smooth muscle cells and esophageal smooth muscle cells (3). In 8/10 fixed tissues, the antibody labelled mesothelial cells (5). Except for cells of muscle origin, normal bone marrow, kidney, breast, prostate gland, brain (gray matter), jejunum, thyroid gland, lung, colon, stomach, ureter and tonsil are negative (1). In fetal frozen tissues, the antibody labels smooth muscle cells in blood vessels and the small intestine, skeletal muscle cells in the tongue, cardiac muscle cells, stromal cells of the medulla of the kidney, the decidua, placental villi and the umbilical cord, submesothelial stromal cells of the pericardium, as well as mesothelial cells (7). <u>Abnormal tissues:</u> Of human rhabdomyosarcomas, the antibody labelled 4/4 tumours (1). In the human esophagus the antibody labelled 35/35 leiomyomas, 1/3 leiomyosarcomas and 3/17 stromal tumours (3). In addition, the antibody labelled 33/33 highly cellular leiomyomas of the uterus (4). Also, 8/16 malignant, diffuse mesotheliomas of the pleura, were labelled by the antibody (5). In human vaginal and endocervical polyps, the antibody labelled both background stromal cells and the atypical stromal cells within the subepithelial zone (8). The antibody showed no labelling in 10 cases each of neuroblastoma, lymphoma, Wilms' tumour, primitive neuroectodermal tumour, and
(102264-001)	M 0760/EFG/HEW/19.12.02 p. 1/4

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