Monoclonal Mouse
Anti-Epstein-Barr Virus, LMP
Clones CS.1-4
Code No./ Code/ Code-Nr. M 0897
Edition/ Edition/ Ausgabe 24.03.03

ENGLISH**Intended use**

For in vitro diagnostic use.

Monoclonal Mouse Anti-Epstein-Barr Virus, LMP, Clones CS.1-4, is intended for use in immunocytochemistry. The antibody labels Epstein-Barr virus-infected cells expressing the latent membrane protein (LMP), and is a useful tool for the identification of tumour cells in Epstein-Barr virus-associated Hodgkin's disease (1-3) and nasopharyngeal carcinoma (4). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Introduction

Epstein-Barr virus (EBV) is a human herpes virus, normally found as a widespread asymptomatic infection with >95% of the world's adult population being lifelong carriers. The virus genome is a doublestranded DNA molecule of approximately 172 kb encoding at least 80 proteins. Two strains of EBV exist, although the bulk of the genome appears to be identical (5, 6). EBV-infected cells and neoplasms generally show one of three different expressional patterns of latency-associated viral genes. In latency type I, only EBV-encoded-nuclear-antigen 1 (EBNA-1) is expressed together with two small polyadenylated nuclear RNAs (EBER 1 and 2) observed in Burkitt's lymphoma. In latency type II, additional expression of three latent membrane proteins, LMP-1, LMP-2A and LMP-2B are observed, seen in Hodgkin's disease and nasopharyngeal carcinoma. Latency type III is seen in lymphoproliferative diseases arising in immunocompromised individuals and EBV-transformed lymphoblastoid cell lines. In this group all six EBNAs (EBNA-1/2/3A/3B/3C/LP), all three LMPs and the two EBERs are expressed (6).

The integral membrane protein, LMP-1, is encoded by the *BNRF1 gene* and contains 386 amino acids with a molecular mass of 63 kDa. It is regarded as an oncogenic protein as it is implicated in at least four transcription factor signalling pathways, and induces expression of multiple cell surface markers and cell adhesion molecules (4, 5). EBV is associated with certain tumours of lymphoid and epithelial origin including Burkitt's lymphoma (BL), immunoblastic lymphoma (IBL), T-cell lymphoma, Hodgkin's Disease (HD), nasopharyngeal carcinoma (NPC) and gastric carcinoma. The only situation where full virus replication can be reliably observed is in hairy leukoplakia, a benign lesion of the tongue frequently encountered in AIDS patients (5, 6).

Reagent provided

A mixture of four monoclonal mouse antibodies provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L Na₂SO₄.

Clones: CS.1, CS.2, CS.3 and CS. 4 (7). Isotype: IgG1, kappa.

Mouse IgG concentration: See label on vial.

Immunogen

Recombinant fusion protein containing the sequence of bacterial beta-galactosidase and the carboxyl half of EBV-encoded LMP (7).

Specificity

In Western blotting with lysates of EBV-transformed lymphoblastoid cell lines, the antibody labels a major band of 57-66 kDa depending on the virus isolate, sometimes accompanied by a minor band of 50-55 kDa corresponding to full length and truncated form of LMP. The antibody recognizes 20 geographically distinct EBV isolates (7). The antibody cross-reacts with a 43-44 kDa doublet of uncharacterized normal cell proteins in cellular extracts including that from EBV-negative control cell lines (7).

The antibody labels a normal EBV-positive lymphoblastoid cell line established by infection of normal resting B-cells with the B95.8 virus, whereas the BL2 cell line, established from a rare EBV-negative sporadic BL-patient, as well as other unrelated cell isolates, including tonsillar B and T lymphocytes, fetal liver and thymus, cultured fetal fibroblasts, human leukaemic cell lines HSB2 and K562 and 7 EBV-negative B-cell lines show no labelling, except the cellular doublet proteins (7).

Precautions

1. For professional users.
2. This product contains sodium azide (Na₃N), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

Storage

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Specimen preparation

Paraffin sections: The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with proteinase K or heat-induced epitope retrieval is recommended. For heat-induced epitope retrieval, optimal results are obtained with Dako Target Retrieval Solution, code No. S 1700, Dako Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen.

Frozen sections and cell preparations: The antibody can be used for labelling acetone-fixed, frozen sections and cell preparations (1, 7).

Staining procedure

Dilution: Monoclonal Mouse Anti-Epstein-Barr Virus, LMP, code No. M 0897, may be used at a dilution range of 1:100-1:200 when applied on formalin-fixed, paraffin-embedded sections of human EBV-positive colon and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody.

Visualization: DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.

Product-specific limitations

(104986-002)

M 0897/EFG/SUA/24.03.03 p. 1/4

In frozen sections, weak cross-reactivity has been noted in smooth muscle cells and a small fraction of plasma cells (1). Reactivity has also been reported for peripheral nerve tissue (8). In formalin-fixed, paraffin-embedded sections, weak cross-reactivity has been detected in layers three and four of the brain cortex, in peripheral nerves and in the basal cell layer of normal skin. Myoepithelial cells and smooth muscle cells also showed faint cross-reactivity (8). For B and T non-Hodgkin's lymphomas, frozen sections yielded better staining results than formalin-fixed, paraffin-embedded sections (8).

Performance characteristics

Cells labelled by the antibody display a cytoplasmic and membrane staining pattern.

Normal tissues: The antibody does not label normal nasopharyngeal mucosa, neither is small infiltrating lymphocytes, desmoplastic stroma and overlying normal epithelium in primary nasopharyngeal samples labelled (4). The antibody displays no labelling in large intestine, appendix or pancreas (8). Strong labelling of normal early myeloid and erythroid precursors may be seen despite a total absence of evidence of EBV by PCR (9).

Abnormal tissues: In frozen sections of 84 cases of Hodgkin's disease 40/84 was strongly labelled by the antibody in Hodgkin and Reed-Sternberg cells. The portions varied according to histological subtype, where staining was demonstrated in 1/10 with nodular lymphocyte predominance, 16/52 with nodular sclerosis and 23/24 with mixed cellularity type (1). The antibody labelled 22/46 cases of paraffin-embedded tissue of Hodgkin's disease, with 1/12 of lymphocyte predominance, 12/24 of nodular sclerosis, 6/7 of mixed cellularity type and 3/3 of lymphocyte depletion subtypes (2). In another study including 47 Hodgkin's lymphomas, a total of 32/47 were labelled by the antibody, including 28/42 of the nodular sclerosis type and 4/5 of mixed cellularity type (3). Among 51 nasopharyngeal carcinomas, the antibody labelled 40/51, including 8/10 keratinising squamous cell carcinomas, 7/12 non-keratinising squamous cell carcinomas and 25/29 undifferentiated carcinomas (4).

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
1. Pallesen G, Hamilton-Dutoit SJ, Rowe M, Young LS. Expression of Epstein-Barr virus latent gene products in tumour cells of Hodgkin's disease. *Lancet* 1991;337:320-2.
2. Murray PG, Young LS, Rowe M, Crocker J. Immunohistochemical demonstration of the Epstein-Barr virus-encoded latent membrane protein in paraffin sections of Hodgkin's disease. *J Pathol* 1992;166:1-5.
3. Engel M, Essop MF, Close P, Hartley P, Pallesen G, Sinclair-Smith C. Improved prognosis of Epstein-Barr virus associated childhood Hodgkin's lymphoma: study of 47 South African cases. *J Clin Pathol* 2000;53:182-6.
4. Vera-Sempere FJ, Burgos JS, Botella MS, Cordoba J, Gobernado M. Immunohistochemical expression of Epstein-Barr virus-encoded latent membrane protein (LMP-1) in paraffin sections of EBV-associated nasopharyngeal carcinoma in Spanish patients. *Oral Oncol, Eur J Cancer* 1996;32B:163-8.
5. Baumforth KRN, Young LS, Flavell KJ, Constandinou C, Murray PG. The Epstein-Barr virus and its association with human cancers [review]. *J Clin Pathol: Mod Pathol* 1999;52:307-22.
6. Cruchley AT, Williams DM, Niedobitek G, Young LS. Epstein-Barr virus: biology and disease [review]. *Oral Diseases* 1997;3 Suppl 1: S156-63.
7. Rowe M, Evans HS, Young LS, Hennessy K, Kieff E, Rickinson AB. Monoclonal antibodies to the latent membrane protein of Epstein-Barr virus reveal heterogeneity of the protein and inducible expression in virus-transformed cells. *J Gen Virol* 1987;68:1575-86.
8. Jiwa NM, Oudejans JJ, Dukers DF, Vos W, Horstman A, van der Valk P, et al. Immunohistochemical demonstration of different latent membrane protein-1 epitopes of Epstein-Barr virus in lymphoproliferative diseases. *J Clin Pathol* 1995;48:438-42.
9. Leong ASY, Cooper K, Leong FJWM. Manual of diagnostic antibodies for immunohistology. Oxford University Press; 199. p. 162.

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