**Monoclonal Mouse  
Anti-Human Epithelial Membrane Antigen**

Clone E29

**Code No./ Code/ Code-Nr. M 0613**

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**ENGLISH**

<b>Intended use</b>	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Epithelial Membrane Antigen, Clone E29, is intended for use in immunocytochemistry. The antibody labels epithelial cells in a wide variety of tissues and is a useful tool for the identification of neoplastic epithelia (1). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
<b>Introduction</b>	Epithelial membrane antigen (EMA) belongs to a heterogenous population of human milk fat globule (HMFG), proteins. HMFG is a complex secretory product of mammary epithelium and EMA can be recovered from the aqueous phase of skimmed milk following extraction in chloroform and methanol. Besides in milk, these proteins are present in a variety of epithelia of both normal and neoplastic types. A number of monoclonal and polyclonal antisera have been raised against these molecules and anti-EMA has been extensively studied in a large number of neoplastic conditions, most often in conjunction with other antibodies (2). EMA is valuable as a marker in the detection of breast carcinoma metastases in histological sections of liver, lymph node, and bone marrow, and is useful for differentiating anaplastic carcinoma from malignant lymphomas, and for the recognition of spindle cell epithelial malignancies (1, 3).
<b>Reagent provided</b>	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN <sub>3</sub> . <u>Clone:</u> E29 (1). <u>Isotype:</u> IgG2a, kappa. <u>Mouse IgG concentration:</u> see label on vial.
<b>Immunogen</b>	Human milk fat globule membrane preparation (1, 4).
<b>Specificity</b>	In Western blotting of the immunogen the antibody labels bands of 265-400 kDa (1).
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
<b>Specimen preparation</b>	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin, Bouin's fixative, B5 fixative or Zenker's solution (5). Pre-treatment of tissues with proteinase K or heat-induced epitope retrieval is recommended. For heat-induced epitope retrieval of tissues fixed in formalin, optimal results are obtained with Dako Target Retrieval Solution, code No. S 1700, Dako Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone-fixed frozen sections (1).
<b>Staining procedure</b>	<u>Dilution:</u> Monoclonal Mouse Anti-Human Epithelial Membrane Antigen, code No. M 0613, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG2a, code No. X 0943, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+HRP kit, code No. K 0679, and DAKO EnVision™+HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the Dako autostainer.
<b>Product-specific limitations</b>	The antibody labels plasma cells, and EMA appears to be common in plasma cell neoplasms, but is also encountered occasionally among other types of lymphoma. However, it should be emphasized that such cases usually are clearly identifiable on purely morphological grounds as being of lymphoid origin (1).
<b>Performance characteristics</b>	In normal breast and other secretory epithelia, labelling is predominantly localized to apical luminal membranes. In neoplasms, cytoplasmic and apical luminal membrane staining are the most common patterns of immunoreactivity with peripheral membrane staining or other patterns also occurring (5). <u>Normal tissues:</u> The antibody labels epithelial cells in a wide variety of tissues and mesothelial cells. Included are sweat ducts and sebaceous glands of the skin, epithelium of the gastrointestinal tract, acini and ducts of the breast, exocrine cells in pancreas, bladder epithelium, distal tubules of the kidney, cervix, endometrium, respiratory epithelium, thyroid and bile ducts. No labelling was observed in epidermis, endocrine cells in pancreas, glomeruli and proximal tubules of the kidney, central nervous system, peripheral nervous system, connective tissue, hepatocytes, and lymphoid tissue, except for occasional plasma cells (1).

**Abnormal tissues:** The antibody labels a wide variety of neoplastic epithelia, and also neoplastic mesothelial cells (1). In a large study of 2081 epithelial, mesenchymal and hematopoietic neoplasms, the antibody provided a 98.6% specificity and a positive predictive value of 99.2% for neoplastic epithelial differentiation when it was used in combination with anti-leucocyte common antigen (LCA) (2). It was shown that the antibody labelled 105/354 cases of soft tissue and intracranial tumours, with most positive cases among synovial sarcomas, spindle cell malignant mesotheliomas, epithelioid sarcomas, chordomas and choroid plexus tumours, 38/169 cases of small round cell tumours and small cell sarcomas, 13/158 cases of germ cell neoplasms, 23/23 cases of spindle cell "sarcomatoid" carcinomas, 815/918 cases of other epithelial malignancies, including squamous, Merkel cell, breast, gastric adeno-, colonic adeno-, pancreatic, salivary gland, bladder, uterine, ovarian, vaginal, pulmonary, prostatic, thyroid, thymic, hepatic, renal cell, and nasopharyngeal carcinomas (2). In lymphomas, the antibody labelled 10/22 cases of null or T-cell type CD30+ ALCL of small cell, monomorphic, or pleomorphic subtypes (6), 34/35 cases of p80/ALK+, 1/6 cases of CD56/57+ T/NK-cell, 2/7 cases of EBV+ cytotoxic large T-cell, 2/8 cases of low-grade cytotoxic T-cell and 6/10 cases of cytotoxic Hodgkin's-like lymphomas (7). No labelling was observed in 18/18 cases of basal cell carcinoma, 12/12 cases of hepatocellular carcinoma, 43/43 cases of malignant melanoma and 69/69 cases of endocrine neoplasms, with the exception of poorly differentiated lesions in 6 cases of islet cell tumours (2).

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
1. Cordell J, Richardson TC, Pulford KAF, Ghosh AK, Gatter KC, Heyderman E, et al. Production of monoclonal antibodies against human epithelial membrane antigen for use in diagnostic immunocytochemistry. Br J Cancer 1985;52:347-54.
2. Swanson PE. Monoclonal antibodies to human milk fat globule proteins. In: Wick MR, Siegal GP, editors. Monoclonal antibodies in diagnostic immunohistochemistry. New York – Basel: Marcell Dekker Inc; 1988. p. 227-83.
3. Sloane JP, Ormerod MG. Distribution of epithelial membrane antigen in normal and neoplastic tissues and its value in diagnostic tumor pathology. Cancer 1981;47:1786-95.
4. Heyderman E, Strudley I, Powell G, Richardson TC, Cordell JL, Mason DY. A new monoclonal antibody to epithelial membrane antigen (EMA) – E29. A comparison of its immunocytochemical reactivity with polyclonal anti-EMA antibodies and with another monoclonal antibody, HMFG-2. Br J Cancer 1985;52:355-61.
5. Pinkus GS, Kurtin PJ. Epithelial membrane antigen – a diagnostic discriminant in surgical pathology. Immunohistochemical profile in epithelial, mesenchymal, and hematopoietic neoplasms using paraffin sections and monoclonal antibodies. Hum Pathol 1985;16:929-40.
6. Felgar RE, Salhany KE, Macon WR, Pietra GG, Kinney MC. The expression of TIA-1+ cytolytic-type granules and other cytolytic lymphocyte-associated markers in CD30+ anaplastic large cell lymphomas (ALCL): correlation with morphology, immunophenotype, ultrastructure, and clinical features. Hum Pathol 1999;30:228-36.
7. Kagami Y, Suzuki R, Taji H, Yatabe Y, Takeuchi T, Maeda S, et al. Nodal cytotoxic lymphoma spectrum. A clinicopathologic study of 66 patients. Am J Surg Pathol 1999;23:1184-1200.

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