**Monoclonal Mouse
Anti-Human CD61,
Platelet Glycoprotein IIIa**

Clone Y2/51

Code No./ Code/ Code-Nr. M 0753

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ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human CD61, Clone Y2/51, is intended for use in immunocytochemistry. CD61 is a selective marker of platelets and platelet precursors, and antibodies to CD61 may be of value for immunophenotyping of megakaryoblastic leukaemias, and for the identification of Glanzmann's thrombastenia (1). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.</p>
Synonyms for antigen	CD61A, GPIIb/IIIa, β 3 integrin chain, β -3 integrin (1).
Introduction	<p>CD61 has a molecular mass of approximately 110 kDa and is the common β-subunit in the CD41/CD61 and in the CD51/CD61 complexes, both of which are noncovalently associated heterodimers. The complexes require divalent cations for complex integrity and receptor functions (1).</p> <p>The CD41/CD61 complex appears early in megakaryocyte maturation (2). The activated CD41/CD61 complex is a receptor for von Willebrand factor, soluble fibrinogen and fibronectin (1) and plays a central role in platelet activation and aggregation (3).</p> <p>The CD51/CD61 plays a role in tumour metastasis and adenoviral infection (1).</p>
Reagent provided	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN₃.</p> <p><u>Clone:</u> Y2/51. <u>Isotype:</u> IgG1, kappa.</p> <p><u>Mouse IgG concentration:</u> See label on vial.</p>
Immunogen	Phytohemagglutinin (PHA)-stimulated peripheral blood cells (4).
Specificity	<p>Anti-CD61, Y2/51, was included in the Fourth International Workshop and Conference on Human Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD61 (5).</p> <p>Anti-CD61, Y2/51, labels platelets in peripheral blood and bone marrow and it also reacts with megakaryocytes (4).</p>
Precautions	<ol style="list-style-type: none">1. For professional users.2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<p><u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is recommended. Optimal results are obtained with, Dako Target Retrieval Solution, code No. S 1700, or 10 mmol/L citrate buffer, pH 6.0 and pre-treatment of tissues with proteinase K. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.</p> <p><u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections, fixed cell smears and fixed cytopsins.</p>
Staining procedure	<p><u>Dilution:</u> Monoclonal Mouse Anti-CD61, code No. M 0753, may be used at a dilution range of 1:50-1:100, in the DAKO EnVision™+ system when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody.</p> <p>Monoclonal Mouse Anti-CD61, code No. M 0753, may be used at a dilution range of 1:50-1:100, in the APAAP technique, when applied on frozen sections, fixed cell smears and fixed cytopsins, using 30 minutes incubation at room temperature with the primary antibody.</p> <p>Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen.</p> <p><u>Visualization:</u> DAKO LSAB™+HRP kit, code No. K 0679, and DAKO EnVision™+HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.</p> <p><u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the Dako Autostainer.</p>
Performance characteristics	<p>Cells labelled by the antibody display staining confined to the cell surface membrane.</p> <p>In peripheral human lymphoid tissue and in normal human tonsil, the antibody shows labelling pattern corresponding to reactivity restricted to platelets.</p>

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
1. Sun Q-H, Newman PJ. CD Guide. CD61. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 810.
2. Vinci G, Tabilio A, Deschamps JF, van Haeke D, Henri A, Guichard J, et al. Immunological study of *in vitro* maturation of human megakaryocytes. Br J Haematol 1984;56:589-605.
3. Nachman RL, Leung LLK. Complex formation of platelet membrane glycoproteins IIb and IIIa with fibrinogen. J Clin Invest 1982;69:263-9.
4. Gatter KC, Cordell JL, Turley H, Heryet A, Kieffer N, Anstee DJ, et al. The immunohistological detection of platelets, megakaryocytes and thrombi in routinely processed specimens. Histopathology 1988; 13:257 -67.
5. von dem Borne AEG, Modderman PW, Admiraal LG, Nieuwenhuis HK. Joint report of the platelet section. P1: Platelets antibodies, the overall results. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p 951-66.

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