**Monoclonal Mouse
Anti-Thyroid Transcription Factor (TTF-1)
Clone 8G7G3/1**

ENGLISH
Code M3575

Intended use

For In Vitro Diagnostic Use.

This antibody is intended for laboratory use to identify qualitatively by light microscopy thyroid transcription factor-1 (TTF-1) in normal and malignant cells using immunohistochemical (IHC) test methods. Positive results aid in the classification of tumors of the lung and thyroid. The clinical interpretation of any positive staining or its absence should be complemented by morphological and histological studies with proper controls. Evaluations should be made within the context of the patient's clinical history and other diagnostic tests by a qualified individual.

Summary and explanation

Thyroid transcription factor-1 (TTF-1) belongs to a family of homeodomain transcription factors and is selectively expressed in thyroid, lung and diencephalon. TTF-1 has been identified as a transcriptional regulator of thyroid-specific genes and has also been shown to be important in the activation of pulmonary-specific differentiation genes.^{2,3}

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Monoclonal Mouse antibody provided in liquid form as tissue culture supernatant in 0.05 mol/L Tris-HCl, pH 7.2 and 0.015 mol/L sodium azide. This product contains stabilizing protein.

Clone: 8G7G3/1¹ Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

M3575 may be used at a dilution of 1:200 when performing IHC using the LSAB+, HRP System (code K0679). M3575 may be used at a dilution of 1:50 when performing IHC using the LSAB2, HRP System, Liquid DAB (code K0673). These are guidelines only; optimal dilutions should be determined by individual laboratories.

The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

Immunogen

Recombinant rat TTF-1¹

Specificity

Monoclonal anti-rat TTF-1 clone 8G7G3/1 (anti-TTF-1) has been shown to react with rat, mouse and human TTF-1. In Western blotting, anti-TTF-1 specifically identified a 40 kD band in immunoblots of nuclear extracts or whole cell lysates from the TTF-1 positive cell lines MLE 15 (mouse lung epithelial), H441-4 (human lung adenocarcinoma), H 345 (human small cell carcinoma) and rat type II pneumocyte cells. Cell lines lacking TTF-1 mRNA such as HeLa and 3T3 were unreactive. Anti-TTF-1 was also found to function in ELISA assays using recombinant TTF-1.¹

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions. Suggested diluent for IHC procedures:

Antibody Diluent (code S0809)

The following negative control is recommended for IHC procedures:

Mouse IgG₁ (code X0931)

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.^{4,5}
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user.⁶ There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin Sections

Anti-TTF-1 can be used on formalin-fixed, paraffin-embedded tissue sections.

The deparaffinized tissue sections must be treated with heat prior to the IHC staining procedure. For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended. For heat induced epitope retrieval, Target Retrieval Solution (code S1700) is recommended using a 40-minute heating protocol.

Cryostat Sections and Cell Smears

Anti-TTF-1 can be used for labelling acetone-fixed cryostat sections or fixed cell smears.

Staining procedure

Follow the procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-TTF-1 is nuclear.

Product specific limitations

1. For optimal staining performance on paraffin-embedded tissues, the use of DAB as a chromogen with the LSAB2, HRP (Liquid DAB) and LSAB+, HRP (Liquid DAB+) Systems is highly recommended.
2. When using a polyclonal antibody to TTF-1, some focal immunoreactivity was observed in 1/66 gastric adenocarcinomas and 1/8 endometrial adenocarcinomas.¹¹

Performance characteristics

Normal Tissues

Anti-TTF-1 has been shown to be immunoreactive with Type II cells and Clara cells of the lung and follicular cells from thyroid, but was unreactive with all other normal tissues examined. The following tissues were found to lack TTF-1 expression: pituitary, prostate, testes, adrenal gland, skin, breast, kidney, colon, liver, pancreas, small intestine, brain and stomach.^{1,7,8}

Abnormal Tissues

TTF-1 expression has been demonstrated by IHC in neoplastic cells derived from tumors of the lung and thyroid. Anti-TTF-1 was found to react positively with the majority of pulmonary small cell carcinomas and primary and metastatic pulmonary adenocarcinomas, whereas a smaller proportion of large cell undifferentiated lung carcinomas (26%) were immunoreactive.^{1,8} TTF-1 has also been reported to react positively with squamous cell carcinomas of the lung (14%).⁹ TTF-1 immunoreactivity was demonstrated in the majority of pulmonary atypical carcinoids but was rare in typical carcinoids.¹⁰ TTF-1 has also been detected in (3/3) thyroid papillary carcinomas.¹ TTF-1 expression was not found in the majority of other tumors tested such as primary breast carcinomas, breast carcinomas metastatic to the lung, renal cell carcinomas, primary and metastatic colon and prostate adenocarcinomas, and malignant mesotheliomas.^{1,7,8,11,12} When using a polyclonal antibody to TTF-1, some focal immunoreactivity was observed in 1/66 gastric adenocarcinomas and 1/8 endometrial adenocarcinomas.¹¹

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Références

Literatur

1. Holzinger A, Dingle S, Bejarano PA, Miller M-A, Weaver TE, Di Lauro R, Whitsett JA. Monoclonal antibody to thyroid transcription factor-1: Production, characterization and usefulness in tumor diagnosis. *Hybridoma* 1996;15(1):49-53
2. Lazzaro D, Price M, De Felice M, Di Lauro R. The transcription factor TTF-1 is expressed at the onset of thyroid and lung morphogenesis and in restricted regions of the fetal brain. *Development* 1991;113(4):1093-104
3. Bohinski RJ, Di Lauro R, Whitsett JA. The lung-specific surfactant protein B gene promoter is a target for thyroid transcription factor-1 and hepatocyte nuclear factor-3, indicating common factors for organ-specific gene expression along the foregut axis. *Mol Cell Biol* 1994;14(9):5671-81
4. Department of Health, Education and Welfare, National Institute for Occupational Safety and Health, Rockville, MD. "Procedures for the decontamination of plumbing systems containing copper and/or lead azides." 1976
5. Center for Disease Control Manual Guide - Safety Management, No. CDC-22, Atlanta, GA. "Decontamination of laboratory sink drains to remove azide salts." April 30, 1976
6. Clinical Laboratory Improvement Amendments of 1988: Final Rule, 57CFR7163. February 28, 1992
7. Bohinski RJ, Bejarano PA, Balko G, Warnick RE and Whitsett JA. Determination of lung as the primary site of cerebral metastatic adenocarcinomas using monoclonal antibody to thyroid transcription factor-1. *J Neuro-Oncol* 1998;40(3):227-31

8. Khor A, Whitsett JA, Stahlman MT, Stephenson M, Olson SJ, Cagle PT. Utility of surfactant protein B precursor and thyroid transcription factor 1 in differentiating adenocarcinoma of the lung from malignant mesothelioma. Hum Pathol 1999;30(6):695-700
9. Haque AK, Syed S, Lele SM, Freeman DH, Adegboyega PA. Immunohistochemical study of thyroid transcription factor-1 and HER2/neu in non-small cell lung cancer: strong thyroid transcription factor-1 expression predicts better survival. Appl Immunohistochem Mol Morph 2002 June;10(2):103-9
10. Folpe AL, Hansen D, Gown AM, Schmidt RA. Expression of thyroid transcription factor-1 (TTF-1) is common in pulmonary atypical carcinoids and small cell carcinomas, but not typical carcinoids. Lab Invest 1998;78(1):174A
11. Bejarano PA, Baughman RP, Biddinger PW, Miller MA, Fenoglio-Preiser C, Al-Kafaji B, Di Lauro R, Whitsett JA. Surfactant proteins and thyroid transcription factor-1 in pulmonary and breast carcinomas. Mod Pathol 1996;9(4):445-52
12. Harlamert HA, Mira J, Yassin R, Bejarano PA, Baughman R, Miller MA, Whitsett J. Distinguishing primary lung adenocarcinomas from metastatic breast adenocarcinomas in cytology specimens using TTF-1 and cytokeratins 7 and 20. Mod Pathol 1997;10(1):A185

REF	Catalogue number Référence du catalogue Bestellnummer		Temperature limitation Limites de température Zulässiger Temperaturbereich		Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i>
	Manufacturer Fabricant Hersteller	LOT	Batch code Code du lot Chargenbezeichnung		Use by Utiliser jusque Verwendbar bis
EC REP	Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU		IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA

Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

EC REP

Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark

Tel +45 4485 9500
Fax +45 4485 9595

www.dako.com

PT0039/Rev C

Edition 01/09