**Monoclonal Mouse  
Anti-Human  
Cytokeratin 5/6**  
Clone D5/16 B4  
**Code No./ Code/ Code-Nr. M 7237**  
Edition/ Edition/ Ausgabe 17.12.02

ENGLISH	
<b>Intended use</b>	For in vitro diagnostic use.  Monoclonal Mouse Anti-Human Cytokeratin 5/6, Clone D5/16 B4, is intended for use in immunocytochemistry. Antibodies to cytokeratin 5/6 (CK 5/6) have been found valuable for the distinction between low differentiated squamous cell carcinoma and adenocarcinoma. Furthermore, the expression of CK5 by mesothelium makes this antibody valuable in the differentiation between epithelioid mesothelioma and lung carcinoma (1) when used together with other antibodies against mesothelioma markers, such as calretinin and thrombomodulin (2). Anti-CK 5/6 has also been found useful in the differential diagnosis of atypical proliferations of the breast (3). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
<b>Introduction</b>	Cytokeratins are alpha-type fibrous polypeptides with a diameter of 7-11 nm. They are important components of the cytoskeleton in almost all epithelial cells as well as in some non-epithelial cell types. Cytokeratins are, generally, held to be the most ubiquitous markers of epithelial differentiation, and, so far, 20 distinct types numbered by Moll (4, 5) have been revealed. The CK 5 is a high molecular weight, basic type of cytokeratin, with a molecular mass of 58 kDa, expressed in the basal, the intermediate and the superficial cell layers of stratified epithelia as well as in transitional epithelia, complex epithelia, and in mesothelial cells and mesothelioma. CK 5 has not, with few exceptions, been found in simple epithelia and in non-epithelial cells. CK 6 is also a high molecular weight, basic type of cytokeratin, with a molecular mass of 56 kDa, expressed by proliferating squamous epithelium often paired with CK 16 (48 kDa) (4, 5).
<b>Reagent provided</b>	Monoclonal mouse antibody provided in liquid form as purified IgG from ascitic fluid. In 0.05 mol/L Tris/HCl, 15 mmol/L Na <sub>3</sub> , 1% bovine serum albumin, pH 7.2. <u>Clone:</u> D5/16 B4. <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial.
<b>Immunogen</b>	Isolated cytokeratin 5.
<b>Specificity</b>	In Western blotting of cytoskeletal preparations of epidermis and non-keratinizing epithelium, the antibody labels cytokeratin 5. It also labels cytokeratin 6 and weakly labels cytokeratin 4. No cross-reaction with cytokeratins 1, 7, 8, 10, 13/14, 18 and 19 has been found.
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (Na <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
<b>Specimen preparation</b>	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with 10 mmol/L citrate buffer, pH 6.0. However, DakoCytomation Target Retrieval Solution, code No. S 1700 and pre-treatment of tissues with proteinase K were found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone-fixed, frozen sections and cell preparations.
<b>Staining procedure</b>	<u>Dilution:</u> Monoclonal Mouse Anti-Human Cytokeratin 5/6, code No. M 7237, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.

(104963-003)

M 7237/EFG/CE/17.12.02 p. 1/4

**Performance characteristics**

Cells positive for cytokeratin 5/6 show cytoplasmic staining.

Normal tissues: The antibody labels stratified squamous epithelia. In complex epithelia, the basal cells in most cases express CK 5. Simple epithelia and non-epithelial cells, generally, are not labelled by the antibody.Abnormal tissues: In formalin-fixed, paraffin-embedded tissues the antibody labelled 56/61 (92%) of epithelioid pleural mesotheliomas, while 9/63 (14%) of metastatic adenocarcinomas were labelled. Reactive mesothelium was also labelled (2). In another study (1), all of 23 mesotheliomas with epithelioid or acinar differentiation were strongly labelled by the antibody, whilst labelling in the sarcomatoid areas of tumours with a mixed morphology was weak or absent. One case of sarcomatoid mesothelioma gave an equivocal reaction, whereas mesotheliomas of desmoplastic or sarcomatoid type were negative. Of 27 secondary adenocarcinomas of the pleura, 22 were negative, 4 weak or equivocal, and 1 focally positive. High levels of cytokeratin 5/6 positivity in ductal hyperplasia of usual type, and negative staining of most cases of atypical ductal hyperplasia and ductal carcinoma in situ, have been demonstrated with the antibody (3).**FRANÇAIS****Intérêt**

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