**Monoclonal Mouse  
Anti-Human  
Cytokeratin, High Molecular Weight  
Clone 34βE12**

**ENGLISH**  
**Code M0630**

**Intended use**

For In Vitro Diagnostic Use.

Monoclonal mouse anti-human cytokeratin, High Molecular Weight (HMW), clone 34βE12 (anti-CK HMW, 34βE12) is intended for laboratory use to identify qualitatively, by light microscopy, the 66, 57, 51 and 49 kD proteins corresponding to cytokeratins 1, 5, 10 and 14 of the Moll catalog, in normal and neoplastic tissues in frozen or formalin-fixed paraffin-embedded tissue using immunohistochemical (IHC) test methods. Positive results aid in the subclassification of neoplastic tissue as carcinoma or epithelial in origin. Because anti-CK HMW, 34βE12 does not react with all carcinomas, another marker for simple epithelium must be incorporated into a panel of antibodies with anti-CK HMW, 34βE12 to aid in the differential diagnosis of anaplastic tumors of unknown origin. The clinical interpretation of any positive staining or its absence should be complemented by morphological and histological studies with proper controls. Evaluations should be made within the context of the patient's clinical history and other diagnostic tests by a qualified individual.

**Summary and explanation**

Cytokeratins are intermediate filament cytoskeletal proteins essential to development and differentiation of epithelial cells. Approximately twenty different cytokeratins have been identified and are classified and numbered according to molecular weight and isoelectric points.<sup>2</sup> In general, most low molecular weight cytokeratins (40-54 kD) are distributed in nonsquamous epithelium, Moll's catalog numbers 7-8 and/or 17-20.<sup>3</sup> High molecular weight cytokeratins (48-67 kD) are found in the squamous epithelium, Moll's catalog numbers 1-6 and/or 9-16.<sup>3</sup>

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

**Reagent provided**

Monoclonal Mouse antibody provided in liquid form as tissue culture supernatant in 0.05 mol/L Tris-HCl, pH 7.2 and 0.015 mol/L sodium azide. This product contains stabilizing protein.

Clone: 34βE12<sup>1</sup> Isotype: IgG<sub>1</sub>, kappa  
Mouse IgG concentration mg/L: See label on vial.

M0630 may be used at a dilution of 1:50 when performing IHC using the LSAB<sup>TM</sup>2, EnVision<sup>TM</sup> or EnVision<sup>TM</sup> Doublestain detection systems. These are guidelines only. Optimal antibody concentrations may vary depending on specimen and preparation method and should be determined individually in each laboratory.

The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

**Immunogen**

Solubilized immunogen keratin extracted from human stratum corneum

**Specificity**

Anti-CK HMW, 34βE12 has been shown to react with the 66, 57, 51 and 49 kD<sup>1</sup> proteins corresponding to cytokeratins 1, 5, 10 and 14 of the Moll catalog.<sup>1,2,4</sup>

Shah et al. suggested clinical usefulness of anti-CK HMW, 34βE12 as an aid in the differentiation of Paget's disease and Bowen's disease (carcinoma in situ) from Pagetoid superficial spreading melanoma.<sup>5</sup> Anti-CK HMW, 34βE12 positivity was reported in Paget's disease of the breast (5/5) and Bowen's disease (10/10), although 25% (1/4) of Paget's disease of the vulva and 0/6 of the Pagetoid superficial spreading melanomas stained.<sup>5</sup>

**Materials required, but not supplied**

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions.

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN<sub>3</sub> may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.<sup>6</sup>
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

## Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

## Specimen preparation

### Paraffin Sections

Anti-Cytokeratin HMW, 34βE12 can be used on formalin-, Methacarn- or Carnoy's-fixed, paraffin-embedded tissue sections. Tissue pretreatment for epitope enhancement is required.

The following enzymes can be used for pretreatment of formalin-fixed, paraffin-embedded tissue: Proteolytic Enzyme, RTU (code S3007) or Proteinase K, RTU (code S3020), Pepsin (code S3002), Pronase (code S2013) or Trypsin (code S2012) for 5 minutes. Rinse thoroughly with distilled water and continue with the staining procedure of the detection system instructions.

As an alternative to enzyme pretreatment, heat-induced epitope retrieval can be used. Heat-induced epitope retrieval involves immersion of tissue sections in a pre-heated buffer solution and maintaining heat, either in a water bath (95–99 °C), a steamer (95–99 °C) or an autoclave (121 °C). For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended. Target Retrieval Solution (code S1700) or 10x Concentrate (code S1699) is recommended using a 20-minute heating protocol. After thermal treatment, allow the jar with buffer and slides to cool for 20 minutes at room temperature. Rinse well with buffer.

### Cryostat Sections and Cell Smears

Anti-CK HMW, 34βE12 can be used for labelling acetone- or methanol-fixed cryostat sections or fixed cell smears.

## Staining procedure

Follow the procedure for the detection system selected.

## Staining interpretation

The cellular staining pattern for anti-CK HMW, 34βE12 is cytoplasmic.

## Product specific limitations

1. Anti-CK HMW, 34βE12 cannot be used as the sole epithelial marker in the differential diagnosis of carcinomas from non-epithelial tumors because it does not react with simple epithelial tissues. It must be used in a panel of antibodies that includes a marker for simple epithelium, e.g. Cytokeratin clone 35βH11.<sup>7</sup>
2. Anti-CK HMW, 35βE12 displays its unique staining pattern in epithelial tissue. It can be used to differentiate the basal cell layer of the prostate, present in normal prostate and variable or absent in cases of prostate adenocarcinoma as well as benign conditions including inflamed acini, atypical adenomatous hyperplasia and postatrophic hyperplasia.<sup>11</sup>
3. Epithelioid sarcoma and synovial sarcoma have been shown to exhibit positive staining with various cytokeratin antibodies due to their epithelial nature.<sup>9,12</sup>
4. False-positive staining of glial cells and tumors with anti-CK HMW, 34βE12 has been reported on formalin-fixed frozen tissue when proteolytic pretreatment was employed. It was shown by immunocytochemical and biochemical methods that these cells and tumors do not express cytokeratins.<sup>10</sup>
5. An ascites preparation of clone anti-CK HMW, 34βE12 demonstrated unexpected positive staining of smooth muscle cells of normal human myometrium in methanol fixed frozen sections. In the same tissue fixed in methacarn and paraffin embedded, no reactivity was observed.<sup>13</sup>
6. It has been reported that antibodies specific for cytokeratin 8 and possibly 18 and 19 gave unexpected positive staining with leiomyomas and leiomyosarcomas.<sup>13,14</sup> The staining of these neoplasms is more apparent on frozen tissue than on formalin-fixed or paraffin-embedded tissue.<sup>13,14</sup> The anti-CK HMW, 34βE12 antibody is rarely positive with leiomyosarcomas; however, these results underscore the importance of using a panel of monoclonal antibodies in the analysis of poorly differentiated neoplasms. The incorporation of anti-muscle specific antibodies, e.g. desmin, into an antibody panel would aid in the correct interpretation of results.<sup>13</sup>

## Performance characteristics

### Normal Tissues

Normal epithelial tissue reactive with anti-CK HMW, 34βE12 include: squamous epithelium and sweat ducts in skin,<sup>4</sup> all epithelial layers including luminal and basal epithelium and ductal cells in breast,<sup>15</sup> some pneumocytes, mesothelium and bronchial epithelium in lung, collecting duct epithelia in kidney and ductal cells in the pancreas, bile ducts in the liver, and mesothelium and a portion of epithelium (cells with a more basal location) of the gastrointestinal tract.<sup>4</sup> Hepatocytes, pancreatic acinar cells, proximal renal tubules and nonepithelial normal tissues are not labeled by anti-CK HMW, 34βE12.<sup>4</sup>

Anti-CK HMW, 34βE12 stained mammary gland of the breast, squamous epithelia of the cervix, prostate gland, excretory duct of the lingual salivary gland, squamous epithelia of the skin, reticular cells and Hassall's bodies of the thymus, and squamous epithelium of the tonsil. No staining was demonstrated on adrenal, bone marrow, brain cerebellum or cerebrum, colon, heart, kidney, liver, lung, mesothelial cells, ovary, pancreas, parathyroid, pericardium, peripheral nerve, pituitary, skeletal muscle, small intestine, spleen, stomach, testis, thyroid or uterus.

### Abnormal Tissues

The anti-CK HMW, 34βE12 antibody positivity has been reported for squamous cell and ductal or transitional cell carcinomas including: squamous cell carcinoma of the skin, lung and nasopharynx, ductal carcinoma of the breast, pancreas, bile duct and salivary gland as well as transitional cell carcinomas of the bladder and nasopharynx and thymomas.<sup>8,13</sup> Anti-CK HMW, 34βE12 stained epithelial mesotheliomas, but failed to react with sarcomatoid or desmoplastic mesotheliomas.<sup>14</sup> Anti-CK HMW, 34βE12-negative epithelial tumors are either "acinar" type adenomas (e.g., pituitary) or adenocarcinomas of simple epithelia (e.g., endometrial carcinomas, renal and hepatocellular carcinomas) or neuroendocrine tumors.<sup>4,9,13</sup> Thus, anti-CK HMW, 34βE12 may be used as an aid in the subclassification of carcinomas.<sup>13</sup> It has been suggested also that anti-CK HMW, 34βE12 is useful in the differential diagnosis of small acinar lesions of the prostate gland.<sup>16</sup> Its diagnostic value may lie in the identification of basal cells. O'Malley et al.<sup>16</sup> saw positive staining of basal cells in 47 examples of benign prostatic lesions including atypical adenomatous hyperplasia (13), basal cell hyperplasia (11), atrophy (16), post-sclerotic hyperplasia (5) and (105464-004)

fibroepithelial nodule (2), while 21 cases of small-acinar adenocarcinomas showed no reactivity. It was stressed that loss of basal cell layer is not uniform in all prostatic adenocarcinomas. Basal cells are lost only in small-acinar adenocarcinomas. It was recommended that complete absence of staining of prostate lesions with anti-CK HMW, 34βE12 should be regarded as very suggestive, but not diagnostic of malignancy.

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


### Références

### Literatur

1. Gown AM, Vogel AM. Monoclonal antibodies to intermediate filament proteins of human cells: Unique and cross-reacting antibodies. *J Cell Biol* 1982;95(2 Pt 1):414-24
2. Moll R, Franke WW, Schiller DL. The catalog of human cytokeratins: Patterns of expression in normal epithelia, tumors and cultured cells. *Cell* 1982;31(1):11-24
3. Miettinen M. Keratin immunohistochemistry: update of applications and pitfalls. *Pathol Ann* 1993;28 Pt 2:113-43
4. Gown AM, Vogel AM. Monoclonal antibodies to human intermediate filament proteins II. Distribution of filaments in normal human tissues. *Amer J Pathol* 1984;114(2):309-21
5. Shah KD, Tabibzadeh SS, Gergler MA. Immunohistochemical distinction of Paget's Disease from Bowen's Disease and superficial spreading melanoma with the use of monoclonal cytokeratin antibodies. *Amer J Clin Pathol* 1987;88(6):689-95
6. Department of Health, Education and Welfare, National Institute for Occupational Safety and Health, Rockville, MD. "Procedures for the decontamination of plumbing systems containing copper and/or lead azides." DHHS (NIOSH) Publ. No. 78-127, Current 13. August 16, 1976
7. Gown AM, Vogel AM. Monoclonal antibodies to human intermediate filament proteins III. Analysis of tumors. *Amer J Clin Pathol* 1985;84(4):413-24
8. Brawer MK, Peehl DM, Stamey TA, Bostwick DG. Keratin immunoreactivity in the benign and neoplastic human prostate. *Canc Res* 1985;45(8):3663-7
9. Hurlimann J, Gardiol D. Immunohistochemistry in the differential diagnosis of liver carcinomas. *Amer J Surg Pathol* 1991;15(3):280-8
10. Bacchi CA, Zarbo RJ, Jiang JJ, Gown AM. Do glioma cells express cytokeratin? *Appl Immunohistochem* 1995;3(1):45-53
11. Gown AM, Boyd HC, Chang Y, Ferguson M, Reichler B, Tippens D. Smooth muscle cells can express cytokeratins of "simple" epithelium. *Amer J Pathol* 1988;132(2):223-32
12. Brown DC, Theaker JM, Banks PM, Gatter KC, Mason DY. Cytokeratin expression in smooth muscle and smooth muscle tumors. *Histopathol* 1987;11(5):477-86
13. Dairkee SH, Puett L, Hackett AJ. Expression of basal and luminal epithelium-specific keratins in normal, benign and malignant breast tissue. *J Nat Can Inst* 1988;80(9):691-5
14. Bolen JW, Hammar SP, McNutt MA. Reactive and neoplastic serosal tissue: a light-microscopic, ultrastructural, and immunocytochemical study. *Amer J Surg Pathol* 1986;10(1):34-47
15. Varma M, Amin MB, Linden MD, Zarbo RJ. Discriminant staining pattern of small glandular and preneoplastic lesions of the prostate using high molecular weight cytokeratin antibody—A study of 301 consecutive needle biopsies. *Mod Pathol* 1997;10:93A
16. O'Malley FP, Grignon DJ, Shum DT. Usefulness of immunoperoxidase staining with high-molecular-weight cytokeratin in the differential diagnosis of small-acinar lesions of the prostate gland. *Virch Arch Pathol Anat* 1990;417(3):191-6

 <b>REF</b> Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i>
 Manufacturer Fabricant Hersteller	 <b>LOT</b> Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis
 <b>EC REP</b> Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU		 <b>IVD</b> In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum

PT0039/Rev C

Dako North America, Inc.  
 6392 Via Real  
 Carpinteria, California 93013 USA

Tel 805 566 6655  
 Fax 805 566 6688  
 Technical Support 800 424 0021  
 Customer Service 800 235 5763

**EC REP**

Dako Denmark A/S  
 Produktionsvej 42  
 DK-2600 Glostrup Denmark

Tel +45 4485 9500  
 Fax +45 4485 9595

www.dako.com

Edition 01/09