Monoclonal Mouse Anti-Human CD68
Clone KP1
Code No./ Code/ Code-Nr. M 0814
Edition/ Ausgabe 20.01.03

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human CD68, Clone KP1, is intended for use in immunocytochemistry. The antibody labels macrophages and other members of the mononuclear phagocyte lineage and is a useful tool for the identification of neoplasms of myeloid and macrophage/monocyte origin (1). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Introduction	CD68 is a highly glycosylated lysosomal membrane protein with an Mr of 110 000. The CD68 protein belongs to a family of lysosomal glycoprotein (LGP)/plasma membrane shuttling proteins that play a role in endocytosis and/or lysosomal trafficking. CD68 is expressed strongly in cytoplasmic granules, and weakly on the surface of macrophages, monocytes, neutrophils, basophils and NK-cells. Additionally, CD68 is expressed by approximately 40% of peripheral blood B cells and is weakly expressed in 50% of B-cell type acute lymphoblastic leukaemia (B-ALL) cells. CD68 can also be found in the cytoplasm of non-haematopoietic tissues, especially the liver, and renal glomeruli, and tubules (2). Unlike many other CD leucocyte antigens, the CD68 molecule is antigenically very heterogenous, and different antibodies to CD68 show different cellular reactivities (3).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L Na ₃ . <u>Clone:</u> KP1 (4). <u>Iso</u> type: IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial.
Immunogen	Lysosomal fraction of human lung macrophages (4).
Specificity	The antibody was clustered as anti-CD68 at the Fourth International Workshop and Conference on Human Leucocyte Differentiation Antigens held in Vienna in 1989 (5). SDS-PAGE analysis of immunoprecipitates formed between the antibody and ¹²⁵ I-labelled lysates from human spleen with B-cell lymphoma rich in macrophages shows reaction with a 110 kDa polypeptide, corresponding to CD68 (4). In Western blotting of extracts of lung, spleen and U937 cells, diffuse 110, 70 and 40 kDa bands were detected when using reducing conditions. Under non-reducing conditions the spleen extract showed an additional 220 kDa band (4).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formol saline, B5 (4) or formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found less efficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody is not recommended for labelling frozen sections, see "Product-specific limitations".
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human CD68, code No. M 0814, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. Follow the procedure enclosed with the selected visualization kit.
Product-specific limitations	The antibody shows a lack of specificity in frozen sections where it labels a number of normal and neoplastic structures of lymphoid and epithelial components (6). Melanomas are commonly immunoreactive for histiocytic markers and the KP1 antibody labels about 68% of melanomas, typically providing diffuse and weak staining. Accordingly, the antibody should be used cautiously in differentiating melanomas from histiocytic lesions (7).
Performance characteristics	Myelomonocytic cells labelled by the antibody display a diffuse or granular cytoplasmic staining. <u>Normal tissues:</u> In normal peripheral blood smears all monocytes and most granulocytes are labelled by the antibody. Tissue macrophages in a wide range of paraffin-embedded tissues are positive, including lung, germinal centre and bone marrow macrophages, and Kupffer cells in the liver. In the bone marrow also myeloid precursors, and many mature granulocytes are strongly labelled (4). Further, microglial cells in the brain (8), osteoclasts in bone, kidney glomeruli and mast cells are strongly labelled, while moderate reaction is observed with kidney tubules, hepatocytes and occasional lymphoid cells (6). Langerhans' cells, interdigitating reticulum cells, and follicular dendritic cells (FDC) are negative, except for weak positivity observed in a minority of FDCs in dermatopathic lymphadenopathy (4).
(102422-001)	M 0814/EFG/CE/20.01.03 p. 1/4

Abnormal tissues: Acute myeloid leukaemias of M1 to M5 type are strongly labelled by the antibody (6). In another study, 20/20 neoplasms of myeloid, myelomonocytic and presumed macrophage derivation showed strong and extensive cytoplasmic reactivity with the antibody. 14/41 B-lineage lymphomas and leukaemias were also positive, but the labelling was usually confined to small dots. The positive B-cell neoplasms were almost all small-cell proliferations (1). In 23/36 (64%) plasma cell hyperplasias, more than 1% of the plasma cells were labelled by the antibody (9). Of 43 primary and metastatic melanomas, 86% were labelled weakly by the antibody (7). All of 22 T-cell lymphomas were negative, as also 12/12 CD30+ anaplastic large-cell lymphomas (1).

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
- Warnke RA, Pulford KAF, Pallesen G, Ralfkiaer E, Brown DC, Gatter KC, et al. Diagnosis of myelomonocytic and macrophage neoplasms in routinely processed tissue biopsies with monoclonal antibody KP1. Am J Pathol 1989;135:1089-95.
- Goyert SM. MC12. CD68 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 1015-16.
- Falini B, Flenghi L, Pileri S, Gambacorta M, Bigerna B, Durkop H, et al. PG-M1: A new monoclonal antibody directed against a fixative-resistant epitope on the macrophage-restricted form of the CD68 molecule. Am J Pathol 1993;142:1359-72.
- Pulford KAF, Rigney EM, Micklem KJ, Jones M, Stross WP, Gatter KC, et al. KP1: a new monoclonal antibody that detects a monocyte/macrophage associated antigen in routinely processed tissue sections. J Clin Pathol 1989;42:414-21.
- Micklem K, Cordell J, Rigney E, Simmons D, Pulford K, Stross P, et al. M13.1. A macrophage-associated antigen defined by five mAb. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 843-46.
- Cordell JL, Falini B, Flenghi L, Jones DB, Pileri S, Radzun HJ, et al. M15. CD68 cluster workshop report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p.925-7.
- Pernick NL, DaSilva M, Gangi MD, Crissman J, Adsay V. "Histiocytic markers" in melanoma. Mod Pathol 1999;12:1072-7.
- Aoki T, Kobayashi K, Isaki K. Microglial and astrocytic change in brains of Creutzfeldt-Jakob disease: an immunocytochemical and quantitative study. Clin Neuropathol 1999;18:51-60.
- Beschorner R, Horny H-P, Petruch UR, Kaiserling E. Frequent expression of haemopoietic and non-haemopoietic antigens by reactive plasma cells: an immunohistochemical study using formalin-fixed, paraffin-embedded tissue. Histol Histopathol 1999;14:805-12.

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