**Polyclonal Rabbit  
Anti-Laminin  
Code No./ Code/ Code-Nr. Z 0097**

ENGLISH	
<b>Intended use</b>	For in vitro diagnostic use. Polyclonal Rabbit Anti-Laminin is intended for use in immunocytochemistry. The antibody labels laminin, which is present in the basement membrane. The antibody may be useful in the characterization of basement membrane preservation, e.g. in breast cancer (1, 2) and adenocarcinomas of the lung (3). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
<b>Introduction</b>	Laminins are large heterotrimeric basement membrane glycoproteins composed of an $\alpha$ -, a $\beta$ - and a $\gamma$ -chain. At present 5 $\alpha$ -chains, 3 $\beta$ -chains, and 3 $\gamma$ -chains are known to form at least 15 different isoforms (4). The laminin isoforms have different tissue- and developmental-specific localizations, and more than 20 cell surface receptors have been identified for laminin, including integrins, proteoglycans, lectins, gangliosides, sulphatides, amyloid precursor protein and galactosyltransferases (5, 6). Laminins have been found to promote cell adhesion, migration, protease activity, proliferation, tumour growth, angiogenesis and metastasis.
<b>Reagent provided</b>	Purified immunoglobulin fraction of rabbit antiserum provided in liquid form in 0.1 mol/L NaCl, 15 mmol/L NaN <sub>3</sub> . <u>Protein concentration g/L:</u> See label on vial.
<b>Immunogen</b>	Laminin isolated from the rat yolk sac tumour cell line L2.
<b>Specificity</b>	The antibody labels laminin in basement membranes. Traces of contaminating antibodies have been removed by solid-phase absorption with human and rat plasma proteins as well as rat immunoglobulins. In crossed immunoelectrophoresis, using 12.5 $\mu$ L antibody per cm <sup>2</sup> gel area against 2 $\mu$ L human plasma, 2 $\mu$ L rat plasma or 20 $\mu$ g of rat immunoglobulins, only one precipitate corresponding to laminin appears. Staining: Coomassie Brilliant Blue. As demonstrated by immunocytochemistry, the antibody cross-reacts with the laminin-equivalent protein in man, mouse (7) and rat (8).
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed, which cannot be explained by variations in laboratory procedures, and a problem with the antibody is suspected, contact our Technical Services.
<b>Specimen preparation</b>	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with DakoCytomation Proteinase K, code No. S 3020 is recommended. Heat-induced epitope retrieval with DakoCytomation Target Retrieval Solution, pH 6, code No. S 2369, DakoCytomation Target Retrieval Solution, code No. S 1700, DakoCytomation Target Retrieval Solution, pH 9, code No. S 2368, or DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, was found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.
<b>Staining procedure</b>	<u>Dilution:</u> Polyclonal Rabbit Anti-Laminin, code No. Z 0097, should be diluted in DakoCytomation Antibody Diluent, code No. S 0809, and may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of tonsil and using 5 minutes pre-treatment at room temperature with DakoCytomation Proteinase K, code No. S 3020, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), code No. X 0936, diluted to the same protein concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4008 and K 4010 are recommended. Follow the procedure enclosed with the kit.
<b>Performance characteristics</b>	Extracellular matrix labelled by the antibody display staining of the basement membrane. <u>Normal tissues:</u> The antibody was found to label basement membranes in human tonsil and kidney. <u>Abnormal tissues:</u> In a study of 220 cases of ER-positive non-invasive breast lesions, labelling of 20 cases with the antibody revealed no distinct immunoreactivity in the surrounding disrupted myoepithelial cell layers, suggesting disruption of the basement membrane also (1).

In a study of 53 cases of lung carcinomas of the adenocarcinoma mixed bronchioloalveolar subtype, labelling with the antibody was lost in collapsed alveolar regions in 40/53 cases and in scarred regions in 53/53 cases indicating complete destruction of the basement membrane in these areas (3)

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
- Man Y-g, Tai L, Barner R, Vang R, Saenger JS, Shekitka KM, et al. Cell clusters overlying focally disrupted mammary myoepithelial cell layers and adjacent cells within the same duct display different immunohistochemical and genetic features: implications for tumor progression and invasion. Breast Cancer Res 2003;5:R231-41.
- Jing X, Kakudo K, Murakami M, Nakamura Y, Nakamura M, Yokoi T, et al. Intraductal spread of invasive breast carcinomas has a positive correlation with c-erb B-2 overexpression and vascular invasion. Cancer 1999;86:439-48.
- Goto K, Yokose T, Kodama T, Nagai K, Nishiwaki Y, Ando M, et al. Detection of early invasion on the basis of basement membrane destruction in small adenocarcinomas of the lung and its clinical implications. Mod Pathol 2001;14:1237-1245.
- Kutleša S, Siler U, Speiser A, Wessels JT, Virtanen I, Rousselle P, et al. Developmentally regulated interactions of human thymocytes with different laminin isoforms. Immunology 2002;105:407-418.
- Engel J. Laminins and other strange proteins. Biochemistry 1992;31:10643-51.
- Engbring JA, Kleinman HK. The basement membrane matrix in malignancy. J Pathol 2003;200:465-70.
- Kamba T, Higashi S, Kamoto T, Shisa H, Yamada Y, Ogawa O, et al. Animal Model. Failure of ureteric bud invasion. A new model of renal agenesis in mice. Am J Pathol 2001;159:2347-52.
- Paku S, Schnur J, Nagy P, Thorgeirsson SS. Origin and structural evolution of the early proliferating oval cells in rat liver. Am J Pathol 2001;158:1313-23.

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