**Monoclonal Mouse  
Anti-Human  
Mast Cell Tryptase**  
Clone AA1  
**Code No./ Code/ Code-Nr. M 7052**  
Edition/ Ausgabe 18.12.02

**ENGLISH****Intended use**

For in vitro diagnostic use.

Monoclonal Mouse Anti-Human Mast Cell Tryptase, Clone AA1, is intended for use in immunocytochemistry. The antibody labels human mast cells and is a useful tool for the recognition of very atypical or immature mast cells (MC) in mast cell leukaemia, and for the detection of small, even minute, dense focal MC infiltrates in staging procedures in patients with known cutaneous mastocytosis (1). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

**Introduction**

Human mast cell tryptases (EC 3.4.21.59) comprise a family of trypsin-like neutral serine proteases that are predominantly expressed in mast cells (2). In its enzymatically active form, mast cell tryptase exists as a non-covalently linked tetramer of 132 kDa (2, 3). Mast cell tryptase is capable of degrading vasoactive intestinal peptide and activating prekallikrein as well as generating kinins, all important mediators involved in bronchoconstriction and airway hyperresponsiveness, which are major contributors to allergic airway disease (2). Mast cell tryptase also exhibits mitogenic effect on human airway smooth muscle cells, and human lung and dermal fibroblasts. Both smooth muscle hyperplasia and fibrotic changes can lead to thickening of the airway wall and a permanent reduction in airway calibre (2). Mast cells are activated by a number of stimuli, including antigen, superoxides, complement protein, neuropeptides and lipoproteins, resulting in activation and degranulation. Mast cell degranulation and thereby release of mast cell tryptase as well as histamine, leukotrienes and cytokines into the surrounding tissue is a pivotal event in an inflammatory response and seems to play an important role in host defense against pathogens (4). Mast cells play an active role in such diverse diseases as atherosclerosis, asthma, arthritis, bile duct fibrosis, malignancy and pulmonary fibrosis (4, 5).

Identification of mast cells through staining of tissues with antibodies specific for human mast cell tryptase has been useful in identification of focal and diffuse MC infiltrates in primary MC disease and mastocytosis. Demonstration of mast cell tryptase is essential for the identification of highly atypical, hypogranulated or even non-metachromatic MCs, especially in MC leukaemia, and for the detection of small or minute MC infiltrates (1).

**Reagent provided**

Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L Na<sub>3</sub>N.

Clone: AA1 (3). Isotype: IgG1, kappa.

Mouse IgG concentration: see label on vial.

**Immunogen**

Mast cell tryptase isolated from human lung (3).

**Specificity**

In Western blotting of a tryptase-rich extract of human lung under reducing conditions, the antibody labels a band of approximately 32.5 kDa corresponding to human mast cell tryptase (3).

In indirect ELISA the antibody reacts with purified human mast cell tryptase (3).

Antibody-binding to the mast cell tryptase does not affect enzyme activity (3).

Periodate treatment of mast cell tryptase has no effect on antibody labelling, indicating that carbohydrate epitopes are not involved in antibody recognition. Nor is antibody recognition affected by the presence of heparin (3).

**Precautions**

1. For professional users.

2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.

3. As with any product derived from biological sources, proper handling procedures should be used.

**Storage**

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Specimen preparation**

Paraffin sections: The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin or Carnoy's reagent, although the structural integrity of the mast cells is better preserved with formalin fixation (6). Pre-treatment of tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with Dako Target Retrieval Solution, code No. S 1700, Dako Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K or trypsin (6) was found less efficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.

Frozen sections and cell preparations: The antibody can be used for labelling Carnoy's or acetone/methanol fixed, frozen sections and cell preparations (3).

**Staining procedure**

Dilution: Monoclonal Mouse Anti-Human Mast Cell Tryptase, code No. M 7052, may be used at a dilution range of 1:100-1:200 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen.

Visualization: DAKO LSAB™+HRP kit, code No. K 0679, and DAKO EnVision™+HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.

### Product-specific limitations

In specimens fixed in Carnoy's reagent, non-specific labelling of collagen has been reported in one study (7) but was not observed in another (6).

### Performance characteristics

Cells labelled by the antibody display a highly specific granular cytoplasmic staining pattern, corresponding to the secretory granules of mast cells (3).

**Normal tissues:** The antibody specifically labels mast cells (MC) in a number of human tissues, including lung, tonsil, colon, gastric mucosa and skin, whereas only a few MCs are detectable in the pituitary gland (6). The antibody does not label basophils, eosinophils, neutrophils, monocytes or lymphocytes and, likewise, no labelling has been observed in keratinocytes (3).

**Abnormal tissues:** In a study of more than 150 cases of different MC disorders, it was demonstrated that  $<4$  MCs/mm<sup>2</sup> were labelled in bone marrow in most normal and reactive states, whereas  $>5$  but  $<100$  were labelled in myelodysplastic syndromes. More than 100 MCs/mm<sup>2</sup> were detected in MC neoplasias, and cases with a much higher number (maximum 2655 MCs/mm<sup>2</sup>) were often observed (1). In a study of 7 patients with seasonal allergic conjunctivitis, the antibody demonstrated an average MC density within the conjunctival stroma of 106/mm<sup>2</sup> compared to 57.4/mm<sup>2</sup> in 8 normal subjects. Intraepithelial tryptase positive cells were also demonstrated in the allergic subjects, whereas no intraepithelial cells were detected in the normal counterparts (7). In the area around the septal bile ducts of the liver, the number of MCs labelled by the antibody in 45 cirrhotic livers averaged 59.0 MCs/mm<sup>2</sup>, a number that was significantly higher than in 71 normal livers, in which only an average density of 39.4 MCs/mm<sup>2</sup> was found (5).

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
- Horny H-P, Valent P. Diagnosis of mastocytosis: general histopathological aspects, morphological criteria, and immunohistochemical findings. *Leuk Res* 2001;25:543-51.
- Abraham WM. Tryptase: potential role in airway inflammation and remodelling. *Am J Physiol Lung Cell Mol Physiol* 2002;282:L193-L196.
- Walls AF, Bennett AR, McBride HM, Glennie MJ, Holgate ST, Church MK. Production and characterization of monoclonal antibodies specific for human mast cell tryptase. *Clin Exp Allergy* 1990;20:581-9.
- Krishnaswamy G, Kelley J, Johnson D, Youngberg D, Stone W, Huang S-K, et al. The human mast cell: functions in physiology and disease. *Front Biosci* 2001;6:d1109-27.
- Koda W, Harada K, Tsuneyama K, Kono N, Sasaki M, Matsui O, et al. Evidence of the participation of peribiliary mast cells in regulation of the peribiliary vascular plexus along the intrahepatic biliary tree. *Lab Invest* 2000;80:1007-17.
- Walls AF, Jones DB, Williams JH, Church MK, Holgate ST. Immunohistochemical identification of mast cells in formaldehyde-fixed tissue using monoclonal antibodies specific for tryptase. *J Pathol* 1990;162:119-26.
- Morgan SJ, Williams JH, Walls AF, Church MK, Holgate ST, McGill JI. Mast cell numbers and staining characteristics in the normal and allergic human conjunctiva. *J Allergy Clin Immunol* 1991;87:111-6.

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