FLEX
Monoclonal Mouse
Anti-Human CD21
 Clone 1F8
Ready-to-Use
 (Link)

Code IR608

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>FLEX Monoclonal Mouse Anti-Human CD21, Clone 1F8, Ready-to-Use (Dako Autostainer/Autostainer Plus), is intended for use in immunohistochemistry together with Dako Autostainer/Autostainer Plus instruments. This antibody labels CD21-expressing follicular dendritic cells and mature B cells and is a useful aid for the classification of malignant lymphomas (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Synonyms for antigen	C3d-receptor, CR2, EBV-receptor (2).
Summary and explanation	<p>CD21 is a transmembrane glycoprotein belonging to a family of complement regulatory proteins, comprising CD35, C4-binding protein, factor H and CD55. CD21 has a Mr of 145 000 and in its soluble form, sCD21, a Mr of 130 000. The extracellular region consists of multiple short consensus repeats, which contain several distinct ligand-binding sites. Binding of the specific ligands results in differential transmembrane signalling via the cytoplasmic tail that has potential protein kinase C and tyrosine kinase phosphorylation sites (2, 3).</p> <p>CD21 is expressed by follicular dendritic cells (FDCs) and mature B cells (3). FDCs form a three-dimensional meshwork in B-cell follicles, which apparently defines the structure of the follicular compartment (1). On lymphoid tissue B cells, CD21 expression is strong on marginal zone, and moderate on mantle zone B cells. Germinal centre B cells have been found to be negative with most CD21 mAbs, and bone marrow B cells have little or no CD21 expression. The expression of CD21 is gradually lost, together with IgD, after stimulation of resting B cells in vitro. Further, CD21 has been found on several types of epithelial cells, and with low expression on T-cell acute lymphoblastic leukemic cells, subsets of normal thymocytes, and mature T cells (3).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure, Materials Required, Not Supplied, Storage, Specimen Preparation, Staining Procedure, Quality Control, Troubleshooting, Interpretation of Staining, General Limitations.</p>
Reagent provided	<p>Ready-to-use monoclonal mouse antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L sodium azide.</p> <p><u>Clone:</u> 1F8 (4). <u>IsoType:</u> IgG1, kappa.</p>
Immunogen	CD21 preparation purified from human tonsils to 60% purity (4).
Specificity	<p>In Western blotting of the immunogen, the antibody labels a band of 145 kDa, corresponding to CD21. The epitope recognized by the antibody is located within a 72 kDa C3d-binding fragment (4).</p> <p>The antibody labels cells or cell lines known to express CD21 (Raji, NC 37, tonsil cells), whereas no labeling is observed in the CD21-negative Jurkat cells (T-cell line) and human erythrocytes (4).</p>
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	<p>Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.</p>
Specimen preparation	<p>The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm. Pre-treatment with heat-induced epitope retrieval (HIER) is required. Optimal results are obtained by pretreating tissues using EnVision FLEX, Target Retrieval Solution, Low pH (50x) (Link) (Code K8005).</p> <p><u>Deparaffinized sections:</u> Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link. For details, please refer to the PT Link User Guide. The following parameters should be used for PT Link: Pre-heat temperature: 65 °C; epitope retrieval temperature and time: 97 °C for 20 (±1) minutes; cool down to 65 °C. Remove Autostainer slide rack with slides from the PT Link tank and immediately dip slides into a jar/tank (e.g., PT Link Rinse Station, Code PT109) containing diluted room temperature EnVision FLEX Wash Buffer (20x) (Link) (Code K8007). Leave slides in Wash Buffer for 1-5 minutes.</p> <p><u>Paraffin-embedded sections:</u> As alternative specimen preparation, both deparaffinization and epitope retrieval can be performed in the PT Link using a modified procedure. See the PT Link User Guide for instructions. After the staining procedure has been completed, the sections must be dehydrated, cleared and mounted using a permanent mounting method.</p> <p>The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of Dako Silanized Slides (Code S3003) is recommended.</p>
Staining procedure	<p>The recommended visualization system is EnVision FLEX+, Mouse, High pH (Link) (Code K8002), replacing the High pH Target Retrieval Solution from this kit with EnVision FLEX Target Retrieval Solution, Low pH (50x) (Link) (Code K8005). The staining steps and incubation times are pre-programmed into the Autostainer Link software. Please refer to the proper Autostainer Link User Guide for detailed instructions on loading slides and reagents. If the protocols are not available on the used Autostainer instrument, please contact Dako Technical Support. All incubation steps should be performed at room temperature.</p> <p>Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory.</p> <p>Counterstaining in hematoxylin is recommended using EnVision FLEX Hematoxylin (Link) (Code K8008).</p> <p>Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include tonsil and the cells/structures should display reaction patterns as described for this tissue in "Performance characteristics". The recommended negative control reagent is FLEX Negative Control, Mouse (Link) (Code IR750)</p>
Staining interpretation	Cells labeled by the antibody display membrane staining.
Performance characteristics	<p><u>Normal tissues:</u> In reactive lymphoid hyperplasia, the antibody strongly labels germinal centre FDCs in 11/11 cases, with a densely meshed FDC staining of the light zone and a loosely arranged, much less compact FDC staining of the dark zone. The antibody labels plasma cells, and gives a faint and less consistent labeling of lymphocytes in the mantle zone. Sinus-lining cells and monocytoïd B cells are weakly labeled (1). The follicular dendritic cells in the</p>

germinal centers of tonsil show a moderate to strong staining reaction, whereas a subset of activated B cells in the mantle zone show a weak to moderate staining reaction.

Abnormal tissues: In six cases of B-cell chronic lymphocytic leukemia, the antibody showed weak surface labeling of neoplastic cells in four cases, and labeled FDCs in the sparse residual germinal centres in four cases. In 8/8 mantle cell lymphomas, the antibody labeled loosely arranged, ill-defined, and expanded FDC meshworks of either nodular or diffuse patterns, resembling broken-up primary follicles. In 11/11 follicular lymphomas, the abnormal follicles demonstrated dense, sharply defined, expanded and sometimes merging FDC meshworks. However, only loose and patchy FDC labeling was observed in five cases with high-grade transformation in the areas of diffuse large cell lymphoma. In 7/7 low-grade MALT-type B-cell lymphomas, the antibody labeled expanded FDC meshworks, with particularly dense and confluent appearance in cases of primary salivary gland and gastric lymphomas. In 5/5 T-cell and histiocyte-rich B-cell lymphomas, the antibody labeled a few compressed residual follicles at the periphery of the lymphomatous involvement. In 9/9 angioimmunoblastic T-cell lymphomas, clusters of dendritic cells with FDC morphology appeared with frequent incorporation of proliferating postcapillary venules. In 4/4 nodular lymphocyte predominance Hodgkin's lymphomas, the antibody labeled enlarged FDC meshworks overlapping the expanded mantle zones. In 15 cases of Hodgkin's lymphoma of nodular sclerosing subtype, the antibody labeled sharply defined, sometimes irregular FDC meshworks surrounding the negative tumor cells in 11 cases of grade I, and showed sparse or no labeling in 4 cases of grade II disease (1). In follicular dendritic sarcoma, the antibody labeled neoplastic cells in 17/17 cases (5). In another study (6) the antibody labeled Reed-Sternberg and Hodgkin's cells in 7/37 cases of nodular sclerosing, 2/41 cases of mixed cellularity, and 5/12 cases of lymphocyte depletion Hodgkin's lymphoma, whereas no labeling of tumor cells was observed in four cases of lymphocyte predominance type. In nine of the cases where Reed-Sternberg and Hodgkin's cells were labeled by the antibody, the cells expressed no other B- or T-cell associated markers (6).

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
References / Bibliographie / Literaturnachweise

1. Bagdi E, Krenacs L, Krenacs T, Miller K, Isaacson PG. Follicular dendritic cells in reactive and neoplastic lymphoid tissues: a reevaluation of staining patterns of CD21, CD23, and CD35 antibodies in paraffin sections after wet heat-induced epitope retrieval. *Appl Immunohistochem Molecul Morphol* 2001; 9:117-24.
2. Timens W. CD Guide. CD21. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference*; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 1125.
3. Timens W. BC5. CD21 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference*; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 140-2.
4. Petzer AL, Schulz TF, Stauder R, Eigentler A, Myones BL, Dierich MP. Structural and functional analysis of CR2/EBV receptor by means of monoclonal antibodies and limited tryptic digestion. *Immunology* 1988; 63:47-53.
5. Chan JKC, Fletcher CDM, Nayler SJ, Cooper K. Follicular dendritic cell sarcoma. Clinicopathologic analysis of 17 cases suggesting a malignant potential higher than currently recognized. *Cancer* 1997;79:294-313.
6. Nakamura S, Nagahama M, Kagami Y, Yatabe Y, Takeuchi T, Kojima M, et al. Hodgkin's disease expressing follicular dendritic cell marker CD21 without any other B cell marker. A clinicopathologic study of nine cases. *Am J Surg Pathol* 1999; 23:363-76.

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Agilent Technologies Singapore (International) Pte Ltd.