FLEX
Monoclonal Mouse
Anti-Human CD4
 Clone 4B12
Ready-to-Use
 (Link)

Code IR649

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>FLEX Monoclonal Mouse Anti-Human CD4, Clone 4B12, Ready-to-Use (Link), is intended for use in immunohistochemistry (IHC) together with Autostainer Link instruments. This antibody labels thymocytes and T-helper cells (1). Results aid in the classification of anaplastic large cell lymphoma (2). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Summary and explanation	<p>CD4 is a 55 kDa glycoprotein with five immunoglobulin-like external domains, a transmembrane domain and a highly conserved intracellular domain. CD4 acts as a co-receptor in MHC class II-restricted antigen-induced T-cell activation and may increase cell adhesion (1). The CD4 molecule also serves as a receptor for the human immunodeficiency virus on T cells and monocytes/macrophages (4).</p> <p>CD4 is expressed during T-cell development but is absent in immature thymocytes (3). It is expressed in T-helper cells, about 80-90% of mature thymocytes and 55-65% of mature peripheral T cells. It is also present on a subset of suppressor or cytotoxic T cells. CD4 is also expressed on monocytes/macrophages, Langerhans' cells and other dendritic cells, but it is not expressed on B cells (1).</p> <p>CD4 plays an important role in immunophenotyping of reactive lymphocytes and lymphoproliferative disorders. The majority of peripheral T-cell lymphomas are derived from the helper T cell subset (1). Antibodies to CD4 may also aid in the classification of mycosis fungoides and unspecified peripheral T-cell lymphoma (3-5).</p> <p>Refer to Dako <i>General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Ready-to-use monoclonal mouse antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L sodium azide.</p> <p><u>Clone:</u> 4B12. <u>Isotype:</u> IgG1, kappa.</p>
Immunogen	<p>Recombinant protein corresponding to the external domain of the CD4 molecule.</p>
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	<p>Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.</p>
Specimen preparation	<p>The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm.</p> <p>Pre-treatment with heat-induced epitope retrieval (HIER) is required. Optimal results are obtained by pretreating tissues using EnVision FLEX Target Retrieval Solution, High pH (50x) (Code 8004).</p> <p><u>Deparaffinized sections:</u> Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link. For details, please refer to the PT Link User Guide.</p> <p>The following parameters should be used for PT Link: Pre-heat temperature: 65 °C; epitope retrieval temperature and time: 97 °C for 20 (±1) minutes; cool down to 65 °C. Remove Autostainer slide rack with slides from the PT Link tank and immediately dip slides into a jar/tank (e.g., PT Link Rinse Station, Code PT109) containing diluted room temperature EnVision FLEX Wash Buffer (20x) (Link) (Code K8007). Leave slides in Wash Buffer for 1-5 minutes.</p> <p><u>Paraffin-embedded sections:</u> As alternative specimen preparation, both deparaffinization and epitope retrieval can be performed in the PT Link using a modified procedure. See the PT Link User Guide for instructions. After the staining procedure has been completed, the sections must be dehydrated, cleared and mounted using a permanent mounting method.</p> <p>The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of Dako Silanized Slides (Code S3003) is recommended.</p>
Staining procedure	<p>The recommended visualization system is EnVision FLEX+, Mouse, High pH (Link) (Code K8002). The staining steps and incubation times are pre-programmed into the Autostainer Link software. Please refer to the proper Autostainer Link User Guide for detailed instructions on loading slides and reagents. If the protocols are not available on the used Autostainer platform, please contact Dako Technical Support. All incubation steps should be performed at room temperature.</p> <p>Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory.</p>

Counterstaining in hematoxylin is recommended using EnVision FLEX Hematoxylin (Link) (Code K8008).

Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include tonsil and liver and the cells/structures should display reaction patterns as described for this tissue in "Performance characteristics". The recommended negative control reagent is FLEX Negative Control, Mouse, (Link) (Code IR750).

Staining interpretation

Cells labeled by the antibody display membrane staining.

Performance characteristics

Normal tissues: In tonsil, the crowded and isolated T-helper cells show a moderate to strong staining reaction. In liver, the Kupffer and endothelial cells of the sinusoids show a weak to moderate staining reaction.

Abnormal tissues: The antibody labeled 1/1 mucosis fungoides, and 2/2 peripheral T-cell lymphomas. In a study of 75 pediatric patients, the antibody labeled 22/51 anaplastic large cell lymphomas (2).

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3. Dabbs DJ. Ed. *Diagnostic Immunohistochemistry*. 2nd edition, Elsevier Inc; 2006. p. 141-50
Agilent Technologies Singapore (International) Pte Ltd.