FLEX
Monoclonal Mouse
Anti-Human CD56
 Clone 123C3
Ready-to-Use
 (Link)

Code IR628

ENGLISH

Intended use	For in vitro diagnostic use. FLEX Monoclonal Mouse Anti-Human CD56, Clone 123C3, Ready-to-Use (Link), is intended for use in immunohistochemistry (IHC) together with Autostainer Link instruments. This antibody is a useful aid for the identification of natural killer (NK) cells, NK-like T cells, neural/neuroendocrine tissues and classification of related neoplasms (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.
Synonyms for antigen	Leu19, neural cell adhesion molecule (NCAM) and NKH1 (2).
Summary and explanation	CD56 is a membrane glycoprotein which has multiple isoforms generated by alternative splicing from a single gene located on chromosome 11 (3). The core polypeptide is the 140 kDa isoform which can be variably glycosylated to produce mature species. The predominant isoform in NK and T cells is the transmembrane-anchored 140 kDa protein. In tissues, CD56 antigen is found in the cerebellum and cortex of the brain, and at neuromuscular junctions. This molecule is identical to the 140 kDa isoform of the neuronal cell adhesion molecule (NCAM) (2, 4). Additionally, CD56 may also be present in certain large granular lymphocyte leukemias, small cell lung carcinomas, neural-derived tumors, myelomas, and myeloid leukemias (2). Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.
Reagent provided	Ready-to-use monoclonal mouse antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L NaN ₃ . <u>Clone:</u> 123C3 (5). <u>Isotype:</u> IgG1, kappa.
Immunogen	Membrane preparation of a small cell lung carcinoma (SCLC) (5).
Specificity	In Western blotting of homogenates of SCLC and of an undifferentiated spindle cell sarcoma, the antibody labels a 29 kDa band under reducing conditions. Under non-reducing conditions the 29 kDa band shifted to a 150 kDa band suggesting that the antibody recognizes a protein consisting of at least two subunits (5).
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use 2. For professional users. 3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm. Pre-treatment with heat-induced epitope retrieval (HIER) is required. Optimal results are obtained by pretreating tissues using EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code 8004). <u>Deparaffinized sections:</u> Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link. For details, please refer to the PT Link User Guide. The following parameters should be used for PT Link: Pre-heat temperature: 65 °C; epitope retrieval temperature and time: 97 °C for 20 (±1) minutes; cool down to 65 °C. Rinse sections with diluted room temperature EnVision FLEX Wash Buffer (20x) (Code K8002). <u>Paraffin-embedded sections:</u> As alternative specimen preparation, both deparaffinization and epitope retrieval can be performed in the PT Link using a modified procedure. See the PT Link User Guide for instructions. After the staining procedure has been completed, the sections must be dehydrated, cleared and mounted using a permanent mounting method. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of Dako Silanized Slides (Code S3003) is recommended.
Staining procedure	The recommended visualization system is EnVision FLEX+, Mouse, High pH (Link) (Code K8002). The staining steps and incubation times are pre-programmed into the Autostainer Link software. Please refer to the proper Autostainer Link User Guide for detailed instructions on loading slides and reagents. If the protocols are not available on the used Autostainer instrument, please contact Dako Technical Support. All incubation steps should be performed at room temperature. Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory. Counterstaining in hematoxylin is recommended using EnVision FLEX Hematoxylin (Link) (Code K8008). Positive and negative control tissue as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include colon and tonsil and the cells/structures should display reaction patterns as described for these tissues in "Performance characteristics". The recommended negative control reagent is FLEX Negative Control, Mouse (Link) (Code IR750).
Staining interpretation	Cells labeled by the antibody display cytoplasmic and/or membrane staining.

Performance characteristics

Normal tissues: The antibody labels rare lymphoid cells in normal and reactive lymph nodes, nasal/nasopharyngeal mucosa, tonsils, thymus and appendix (1). The antibody also labels the neuroendocrine elements of pituitary and adrenal glands, thyroid follicular epithelium and parafollicular cells, islets of Langerhans, satellite cells and Schwann cells in non-myelinated fibers, Schwann cells of the myelinated nerve fibers usually stained at the Ranvier nodes and Schmidt-Lanterman clefts, peptic cells and chief cells in the stomach, pulmonary veins of the lung, muscularis propria of the esophagus, chief and peptic cells of the stomach, muscularis mucosae of the small and large intestine, some cells in the pancreatic ducts, Henle's loops of the kidney, the ducts of the prostatic gland, the endocrine cells of the testis, rete testis cells at basolateral surface, Leydig cells and boundary tissue of the testis, stroma of the ovary, some cells in the resting acini and ducts of the breast, some ganglion cells and smooth muscle in bronchial and bronchiolar system stained weakly and with varying intensity it also label muscle cells (5). In colon, the ganglion cells and axons of the Auerbach's plexus show a moderate to strong staining reaction. In tonsil, the isolated T cells (NK cells) in the interfollicular areas show a weak to moderate staining reaction.

Abnormal tissues: The antibody labeled 32 cases of CD56+ T/NK-cell lymphomas including 22 cases of nasal/nasopharyngeal lymphomas and 10 cases of extranasal lymphomas. The antibody labeled one case NKH1+ B-lymphoblastic lymphoma. No staining was observed in 24 cases of CD56- T/NK-cell lymphomas and 50 cases of B-cell lymphomas (1).

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

1. Tsang WY, Chan JKC, Ng CS, Pau MY. Utility of a paraffin section-reactive CD56 antibody (123C3) for characterization and diagnosis of lymphomas. Am J Surg Pathol 1996; 20:202-10.
2. Poggi A. CD Guide. CD56. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 1155-6.
3. Costa P. NK4. CD56 workshop panel report In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 271-2.
4. Leong AS-Y, Copper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 1999. p. 101-2.
5. Schol DJ, Mooi WJ, van der Gugten AA, Wagenaar S Sc, Hilgers J. Monoclonal antibody 123C3, identifying small cell carcinoma phenotype in lung tumours, recognizes mainly, but not exclusively, endocrine and neuron-supporting normal tissues. Int J Cancer 1988;Supplement 2:34-40.

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Agilent Technologies Singapore (International) Pte Ltd.