**Monoclonal Mouse Anti-Human Cytokeratin**  
Clone MNF116  
**Code No./ Code/ Code-Nr. M 0821**  
Edition/ Ausgabe 18.12.02

ENGLISH	
<b>Intended use</b>	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Cytokeratin, Clone MNF116, is intended for use in immunocytochemistry. The antibody labels epithelial tissues from simple glandular to stratified squamous epithelium, and is a useful tool for the identification of normal and neoplastic cells of epithelial origin (1, 2). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
<b>Introduction</b>	The cytokeratins (CKs) belong to the intermediate filaments, which create a cytoskeleton in almost all eukaryotic cells. In contrast to other intermediate filaments, CKs are made up of a highly complex multigene family of polypeptides with molecular masses ranging from 40 to 68 kDa. CKs are generally held to belong to the most fundamental markers of epithelial differentiation, and until now, 20 distinct CK polypeptides have been revealed in various human epithelia (3). The CKs can be divided into an acidic type A (class I) and a neutral-basic type B (class II) subfamily.
<b>Reagent provided</b>	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN <sub>3</sub> . <u>Clone:</u> MNF 116. <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial.
<b>Immunogen</b>	Crude extract of splenic cells from a nude mouse engrafted with MCF-7 cells (human breast carcinoma cell line).
<b>Specificity</b>	The antibody is a broad spectrum anti-keratin reagent reacting with intermediate and low-molecular-weight keratins. Thus in immunoblotting, the antibody labels a number of discrete bands ranging from 40 to 58 kDa, corresponding to cytokeratins 5, 6, 8, 17 and probably 19. As demonstrated by immunocytochemistry, the antibody cross-reacts with the cytokeratin-equivalent protein in cow (4), mouse (5), rabbit (6) and rat.
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
<b>Specimen preparation</b>	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with trypsin (1, 7), pepsin (2), proteinase K or heat-induced epitope retrieval is required. For heat-induced epitope retrieval, optimal results are obtained with DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, or 10 mmol/L citrate buffer, pH 6.0. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections as demonstrated with the DakoCytomation EPOS-Conjugate, code No. U 7022, (8).
<b>Staining procedure</b>	<u>Dilution:</u> Monoclonal Mouse Anti-Human Cytokeratin, code No. M 0821, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+HRP kit, code No. K 0679, and DAKO EnVision™+HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.
<b>Performance characteristics</b>	Cells labelled by the antibody display a cytoplasmic staining pattern. <u>Normal tissues:</u> The antibody labels a number of different epithelia, e.g. simple epithelia in biliary and pancreatic ducts, proximal and distal convoluted renal tubules, thyroid-, parathyroid-, gastric-, gall bladder-, small intestinal-, large intestinal-, and endometrial epithelium, as well as hepatocytes (1). The epithelial lining of the rete ovarii, epoophoron, Fallopian tube and ovarian surface is also reactive with the antibody (7). In stratified squamous epithelium of skin, the antibody labels the basal layer, as well as the outer root sheath of follicular epithelium, and the basal layer of sebaceous, eccrine, and apocrine glands and ducts (2). Non-squamous stratified epithelia of esophagus and ectocervix are also reported positive with the MNF116 antibody (1). Transitional epithelium (urothelium) and different complex epithelia, e.g. ducts and acini of breast, endocervical glands, prostate epithelium, epididymis, bronchial epithelium as well as ducts and acini of salivary glands are also labelled by the antibody (1). Mesenchymal cells are generally not reactive with the antibody, however a weak and
(102466-001)	M 0821/EFG/HEW/18.12.02 p. 1/4

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focal labelling of smooth muscle cells is reported. Endothelium, skeletal muscle, cartilage, and lymphoid tissues are not labelled by the antibody (1).

Abnormal tissues: Of epithelial tumours, the antibody labelled 4/4 colorectal carcinomas, 3/3 gastric carcinomas, 4/4 breast carcinomas, 3/3 prostatic carcinomas, 3/3 renal cell carcinomas, 3/4 hepatocellular carcinomas, 3/3 transitional cell carcinomas, 2/3 carcinoids of appendix, and, moreover, 2/2 teratomas and 2/2 pleomorphic adenomas (epithelial elements), and 14/14 squamous cell carcinomas, i.e. 6 of epidermal-, 3 of cervical- and 5 of bronchial origin (1). In addition, 16/16 cases of mixed tumours and myoepitheliomas arising in soft tissues, displayed positive reactivity with the antibody (9). Positive labelling of mesothelioma has also been observed. The antibody, further, labels micrometastatic cells of non-small cell lung carcinomas in lymph nodes (10). Non-epithelial tumours, e.g. melanomas, lymphomas, benign and malignant fibrous histiocytomas, leiomyomas, liposarcomas, chondrosarcomas, Ewings sarcomas, angiofibromas, angiosarcomas, atypical fibroxanthomas, juvenile xanthogranulomas, hemangiomas, granular cell tumours, hair follicle myxomas, epithelioid histiocytomas and dermatofibrosarcomas are not labelled by the antibody. However, in 1/6 cases (1), and 1/2 cases (2) of leiomyosarcomas, the antibody displayed a weak and focal reactivity.

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
<b>References/ Références/ Literatur</b>	<ol style="list-style-type: none"><li>Goddard MJ, Wilson B, Grant JW. Comparison of commercially available cytokeratin antibodies in normal and neoplastic adult epithelial and non-epithelial tissues. J Clin Pathol 1991;44:660-3.</li> <li>Prieto VG, Lugo J, McNutt NS. Intermediate- and low-molecular-weight keratin detection with the monoclonal antibody MNF116. An immunohistochemical study on 232 paraffin-embedded cutaneous lesions. J Cutan Pathol 1996;23:234-41.</li> <li>Moll R. Cytokeratins as markers of differentiation in the diagnosis of epithelial tumors. In: Hermann, Harris, editors. Subcellular Biochemistry. Volume 31, New York: Plenum Press; 1998. p. 205-62.</li> <li>Berkovitz BKB, Whatling R, Barrett AW, Omar SS. The structure of bovine periodontal ligament with special reference to the epithelial cell rests. J Periodontol 1997;68:905-13.</li> <li>Lee I, Yu E, Ikehara S. Thymic epithelial cells of severe combined immunodeficiency (SCID) mice. J Korean Med Sci 1994;9:35-41.</li> <li>Millar TJ, Herok G, Koutavas H, Martin DK, Anderton PJ. Immunohistochemical and histochemical characterisation of epithelial cells of rabbit lacrimal glands in tissue sections and cell cultures. Tissue Cell 1996;28:301-12.</li> <li>Woolnough E, Russo L, Khan MS, Heatley MK. An immunohistochemical study of the rete ovarii and epoophoron. Pathology 2000;32:77-83.</li> <li>Richter T, Nährig J, Komminoth P, Kowolik J, Werner M. Protocol for ultrarapid immunostaining of frozen sections. J Clin Pathol 1999;52:461-3.</li> <li>Kilpatrick SE, Hitchcock MG, Kraus MD, Calonje E, Fletcher CDM. Mixed tumors and myoepitheliomas of soft tissue. A clinicopathologic study of 19 cases with a unifying concept. Am J Surg Pathol 1997;21:13-22.</li> <li>Nicholson AG, Graham ANJ, Pezzella F, Agnetta G, Goldstraw P, Pastorino U. Does the use of immunohistochemistry to identify micrometastases provide useful information in the staging of node-negative non-small cell lung carcinomas? Lung Cancer 1997;18:231-40.</li></ol>
--	--

<b>Explanation of symbols/ Légende des symboles/ Erläuterung der Symbole</b>			
<b>REF</b>	Catalogue number Référence du catalogue Bestellnummer	 <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p>	 <p>Manufacturer Fabricant Hersteller</p>
<b>IVD</b>	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	<b>LOT</b>	Batch code Code du Lot Chargenbezeichnung
	Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten		Use by Utiliser jusque Verwendbar bis

(102466-001)	M 0821/EFG/HEW/18.12.02 p. 3/4
DakoCytomation Denmark A/S · Produktionsvej 42 · DK-2600 Glostrup · Denmark · Tel. +45 44 85 95 00 · Fax +45 44 85 95 95 · CVR No. 33 21 13 17	

DakoCytomation Denmark A/S · Produktionsvej 42 · DK-2600 Glostrup · Denmark · Tel. +45 44 85 95 00 · Fax +45 44 85 95 95 · CVR No. 33 21 13 17

M 0821/EFG/HEW/18.12.02 p. 4/4