FLEX
Polyclonal Rabbit
Anti-Human
Lambda Light Chains
Ready-to-Use
 (Dako Omnis)

Code GA507

ENGLISH

Intended use For in vitro diagnostic use.
 FLEX Polyclonal Rabbit Anti-Human Lambda Light Chains, Ready-to-Use (Dako Omnis), is intended for use in immunohistochemistry together with the Dako Omnis instrument. This antibody is useful for the identification of plasma cells and related lymphoid cells containing lambda light chains, and it is a useful aid for the classification of monoclonal gammopathies (1) and neoplastic monoclonal proliferation (2-4). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.

Summary and explanation The human immunoglobulins basically consist of two identical heavy chains (Mr about 50 000) and two identical light chains (Mr about 20 000). The four chains are covalently linked together by disulphide bonds. The light chains are either kappa or lambda. The five immunoglobulins, IgA, IgD, IgE, IgG and IgM, differ regarding the heavy chains, and IgA and IgM also occur in polymeric forms (5).
 Individual B cells express either kappa or lambda light chains, but never both. In a polyclonal population the ratio of kappa to lambda bearing B cells is 2:1, and the occurrence of a mixture of kappa and lambda light chain-bearing cell types suggests polyclonality and a reactive or non-neoplastic proliferation of B cells (6). Multiple myeloma and B-cell lymphoma are characterized by the proliferation of monoclonal neoplastic plasma cells producing only one type of light chain. Thus, the demonstration of a kappa or lambda light chain-restricted cell population is very useful in the classification of these (7).
 Refer to *Dako General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.

Reagent provided Ready-to-use polyclonal rabbit antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L sodium azide.

Immunogen Polyclonal light chains of type lambda isolated from a pool of normal human sera

Specificity The antibody reacts with free lambda chains as well as lambda chains in intact immunoglobulin molecules. Traces of contaminating antibodies have been removed by solid-phase absorption with human plasma proteins.
 Specificity tests; staining with Coomassie Brilliant Blue:
Double immunodiffusion: No reaction with kappa or gamma chains is seen when using 15 µL of concentrated antibody against 15 µL of kappa IgG.
Crossed immunoelectrophoresis: Only lambda-related precipitates appear when using 12.5 µL concentrated antibody per square cm gel area against 2 µL human plasma or 2 µL of a pool of concentrated urine from patients with tubular proteinuria.

- Precautions**
1. For in vitro diagnostic use.
 2. For professional users.
 3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
 4. As with any product derived from biological sources, proper handling procedures should be used.
 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
 6. Unused solution should be disposed of according to local, State and Federal regulations.

Storage Store at 2-8 °C. During storage the cap should be closed. Do not use after expiration date stamped on vial. Onboard stability is 375 hours. Onboard time is tracked by the Dako Omnis software. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Step		Comments
Fixation/embedding	Formalin-fixed, paraffin-embedded	Onboard deparaffinization
Pre-treatment	EnVision FLEX, Low pH (Code GV805)	30 min HIER
Antibody	Ready-to-use	12.5 min incubation
Negative Control	FLEX Negative Control, Rabbit (Code GA600)	12.5 min incubation
Visualization	EnVision FLEX (Code GV800)	Block: 3 min; Polymer: 20 min; Chromogen: 5 min
Counterstain	Hematoxylin (Code GC808)	3 min incubation
Control Tissue	Tonsil	Cytoplasmic and/or membranous staining
Slides	FLEX IHC Microscope Slides (Code K8020)	Recommended for greater adherence of tissue sections to glass slides
Mounting	Non-aqueous, permanent mounting required	After staining, the sections must be dehydrated, cleared and mounted using a permanent mounting method
Instrumentation	Dako Omnis	Reagents are provided in instrument-specific vials

*The user must always read the package insert for detailed instructions of the staining procedure and handling of the product.

Specimen preparation Paraffin sections: The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of 4 µm.

Pre-treatment: Pre-treatment of formalin-fixed, paraffin-embedded tissue sections with heat-induced epitope retrieval (HIER) is required. Pretreating tissues with HIER using diluted EnVision FLEX Target Retrieval Solution, Low pH (50x) (Dako Omnis), Code GV805, is recommended. Deparaffinization, rehydration and target retrieval are performed onboard Dako Omnis. Please refer to Dako Omnis Basic User Guide.

Dako Omnis ensures that the tissue sections do not dry out during the pre-treatment process or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides, Code K8020 is recommended.

Staining procedure Program: The staining steps and incubation times are pre-programmed into the Dako Omnis software. Please refer to the Dako Omnis Basic User Guide for detailed instructions on loading slides and reagents. If the protocols are not available in the Dako Omnis system, please contact Dako Technical Support. All incubation steps are performed at 32 °C onboard Dako Omnis.

Visualization: The recommended visualization system is EnVision FLEX, High pH (Dako Omnis), Code GV800. The visualization is performed onboard Dako Omnis.

Note: Use EnVision FLEX Target Retrieval Solution, **Low pH** (50x) (Dako Omnis), Code GV805, for HIER.

Counterstaining: The recommended counterstain is Hematoxylin (Dako Omnis), Code GC808. The counterstaining is performed onboard Dako Omnis.

Mounting: After staining onboard Dako Omnis the sections must be dehydrated, cleared and mounted using permanent mounting method.

Quality control: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include tonsil and the cells/structures should display reaction patterns as described for this tissue in the "Performance characteristics" section. The recommended negative control reagent is FLEX Universal Negative Control, Rabbit (Dako Omnis), Code GA600.

Staining interpretation

The cellular staining pattern for lambda light chains may be granular cytoplasmic and/or perinuclear, or may be observed as membranous.

Performance characteristics

Normal tissues: In normal tonsil, the antibody labels plasma cells, follicle centre cells and mantle zone cells providing a polyclonal pattern with a clear distinction between labeled and non-labeled cells. In tonsil, approximately 40-50% of the mantle zone B cells show weak to moderate membranous staining, while immunoblasts and plasma cells in the germinal center show moderate to strong cytoplasmic staining. The B lymphocytes should not be labeled or only weakly labeled in the germinal center. Some background reaction in serum, connective tissue and epithelial cells may be observed.

Abnormal tissues: In reactive lymph nodes, a characteristic staining pattern was demonstrated with Dako Anti-Lambda, and an indistinguishable pattern was obtained with Dako Anti-Kappa clearly demonstrating the polyclonality of the B cells (3). In plasmacytomas, the abnormal plasma cells were monotypic and either kappa or lambda-labeled (2). In 113 cases of formalin-fixed B-cell non-Hodgkin lymphoma, including several small core biopsy specimens with extremely limited tissue, monotypic light chain expression was demonstrated in 91 (81%) of the cases with a 100% specificity using Dako Anti-Kappa and Anti-Lambda (4).

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
References/ Bibliographie/ Literaturnachweis

- Peterson LC, Brown BA, Crosson JT, Mladenovic J. Application of the immunoperoxidase technic to bone marrow trephine biopsies in the classification of patients with monoclonal gammopathies. Am J Clin Pathol 1986;85:688-93.
- Taylor CR, Burns J. The demonstration of plasma cells and other immunoglobulin-containing cells in formalin-fixed, paraffin-embedded tissues using peroxidase-labelled antibody. J Clin Path 1974;27:14-20.
- Harris NL, Poppema S, Data RE. Demonstration of immunoglobulin in malignant lymphomas. Use of an immunoperoxidase technic on frozen sections. Am J Clin Pathol 1982;78:14-21.
- Marshall-Taylor CE, Cartun RW, Mandich D, DiGuseppe JA. Immunohistochemical detection of immunoglobulin light chain expression in B-cell non-Hodgkin lymphomas using formalin-fixed, paraffin-embedded tissues and a heat-induced epitope retrieval technique. Appl Immunohistochem Mol Morphol 2002;10:258-62.
- Klein J, Hořejší V. Immunology. 2nd ed. Abingdon (UK): Blackwell Science Ltd; 1999. p. 226-7.
- Leong ASY, Cooper K, Leong FJWA. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 1999. p. 218.
- Jessup E. Antigen retrieval techniques for the demonstration of immunoglobulin light chains in formalin-fixed, paraffin-wax embedded sections. UK NEQAS immunocytochemistry news. Winter 1994/95; Issue 4:12-6.

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Agilent Technologies Singapore (International) Pte Ltd.