FLEX
Monoclonal Mouse
Anti-Human
Epithelial Antigen
 Clone Ber-EP4
Ready-to-Use
 (Dako Omnis)
Code GA637

ENGLISH

Intended use

For in vitro diagnostic use.

FLEX Monoclonal Mouse Anti-Human Epithelial Antigen, Clone Ber-EP4, Ready-to-Use (Dako Omnis), is intended for use in immunohistochemistry (IHC) together with the Dako Omnis instrument. This antibody labels most epithelial cells and is a useful aid in the classification of adenocarcinoma (1). The antibody may also aid in the classification of esophageal carcinoma (2) and basal and squamous cell carcinoma of the skin (3). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.

Synonyms for antigen

Ep-CAM

Summary and explanation

Epithelial antigen is a transmembrane glycoprotein functioning as a cellular adhesion molecule. This epithelium-specific antigen is broadly distributed in epithelial cells, and displays a highly conserved expression in carcinomas (4, 5). As exceptions to the general expression in normal epithelia, adult hepatocytes, in contrast to fetal hepatocytes, parietal cells in gastric glands, and apical cells in squamous epithelia are negative. Epithelial antigen may rarely be present in mesotheliomas (1, 4). It has been reported that epithelial antigen may play an important role as tumor-cell marker in lymph nodes from patients with esophageal carcinoma otherwise classified as node-negative (2).

Refer to *Dako General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.

Reagent provided

Ready-to-use monoclonal mouse antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L sodium azide.

Clone: Ber-EP4 (4). Isotype: IgG1, kappa.

Immunogen

MCF-7 cells (human breast carcinoma cell line) (4).

Specificity

Analysis of immunocomplexes between the antibody and lysate of ¹²⁵I surface-labeled MCF-7 cells in SDS-PAGE under reducing conditions shows that the antibody labels two polypeptides of 34 kDa and 39 kDa, respectively, corresponding to epithelial antigen. Under non-reducing conditions, the polypeptides appear as 39 kDa and 41 kDa, while deglycosylation reduces the size to 31 kDa, and 36 kDa.

In immunoprecipitation experiments, the Ber-EP4 antibody blocks the reaction of the HEA125 antibody with MCF-7 cell lysate and vice versa, showing that the two antibodies react with the same antigen. The two antibodies also produce identical staining results in cells and tissues (4).

Of 37 cell lines tested, the antibody homogeneously labels all (10/10) carcinoma cell lines, whereas all non-epithelial cell lines (26/27) are not labeled except for the erythromyeloid cell line K562 (4).

Precautions

1. For in vitro diagnostic use.
2. For professional users.
3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
4. As with any product derived from biological sources, proper handling procedures should be used.
5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
6. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store at 2-8 °C. During storage the cap should be closed. Do not use after expiration date stamped on vial. Onboard stability is 375 hours. Onboard time is tracked by the Dako Omnis software. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Staining protocol overview*

Step		Comments
Fixation/embedding	Formalin-fixed, paraffin-embedded	Onboard deparaffinization
Pre-treatment	EnVision FLEX, Low pH (Code GV805)	30 min HIER
Antibody	Ready-to-use	20 min incubation
Negative Control	FLEX Negative Control, Mouse (Code GA750)	20 min incubation
Visualization	EnVision FLEX (Code GV800) + EnVision FLEX+ Mouse LINKER (Code GV821)	Block: 3 min; Link: 10 min; Polymer: 20 min; Chromogen: 5 min
Counterstain	Hematoxylin (Code GC808)	3 min incubation
Control Tissue	Colon, kidney	Cytoplasmic and membrane staining
Slides	FLEX IHC Microscope Slides (Code K8020)	Recommended for greater adherence of tissue sections to glass slides
Mounting	Non-aqueous, permanent mounting required	After staining, the sections must be dehydrated, cleared and mounted using a permanent mounting method
Instrumentation	Dako Omnis	Reagents are provided in instrument-specific vials

*The user must always read the package insert for detailed instructions of the staining procedure and handling of the product.

Specimen preparation

Paraffin sections: The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of 4 µm.

Pre-treatment: Pre-treatment of formalin-fixed, paraffin-embedded tissue sections with heat-induced epitope retrieval (HIER) is required. Pretreating tissues with HIER using diluted EnVision FLEX Target Retrieval Solution, Low pH (50x) (Dako Omnis), Code GV805, is recommended. Deparaffinization, rehydration and target retrieval are performed onboard Dako Omnis. Please refer to Dako Omnis Basic User Guide.

Dako Omnis ensures that the tissue sections do not dry out during the pre-treatment process or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides, Code K8020 is recommended.

Staining procedure

Program: The staining steps and incubation times are pre-programmed into the Dako Omnis software. Please refer to the Dako Omnis Basic User Guide for detailed instructions on loading slides and reagents. If the protocols are not available in the Dako Omnis system, please contact Dako Technical Support. All incubation steps are performed at 32 °C onboard Dako Omnis.

Visualization: The recommended visualization system is EnVision FLEX, High pH (Dako Omnis), Code GV800 in combination with EnVision FLEX+ Mouse (LINKER) (Dako Omnis), Code GV821. The visualization is performed onboard Dako Omnis.

Note: Use EnVision FLEX Target Retrieval Solution, **Low pH** (50x) (Dako Omnis), Code GV805, for HIER.

Counterstaining: The recommended counterstain is Hematoxylin (Dako Omnis), Code GC808. The counterstaining is performed onboard Dako Omnis.

Mounting: After staining onboard Dako Omnis the sections must be dehydrated, cleared and mounted using a permanent mounting method.

Quality control: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include colon and kidney and the cells/structures should display reaction patterns as described for this tissue in the "Performance characteristics" section. The recommended negative control reagent is FLEX Universal Negative Control, Mouse (Dako Omnis), Code GA750.

Staining interpretation

Cells labeled by the antibody display cytoplasmic and membranous staining. The membrane staining is preferentially basolateral (4).

Performance characteristics

Normal tissues: All normal epithelial tissues are labeled by the antibody. Epithelial cells of different origin display varying levels of staining, but most epithelia are strongly positive. Only parietal cells in gastric glands, apical cell layers in squamous epithelia, and adult hepatocytes are negative (4). The antibody does not label non-epithelial tissues, including spleen, peripheral blood, bone marrow, brain, connective tissue, smooth and striated muscle, heart, endothelia, and myoepithelia. Additionally pleura and peritoneum-lining cells are negative, but cells covering the ovary display a slight staining (4). In colon, the column epithelial cells show a moderate to strong staining reaction. In kidney, the epithelial cells lining the Bowman capsule show a weak to moderate staining reaction.

Abnormal tissues: The antibody labeled 142 of 144 epithelial tumor specimens, irrespective of their differentiation, derived from breast, esophagus, stomach, colon, rectum, pancreas, kidney, liver, lung, thyroid and salivary glands, vagina, ovary, cervix uteri and nasopharynx, reflecting the staining pattern in their non-malignant counterparts. Hepatocellular carcinomas displayed heterogeneous staining and included the two non-labeled cases. In this study (4), 2 of 2 squamous cell carcinomas of the lung and cervix uteri, respectively, were labeled in formalin-fixed, paraffin-embedded specimens, even though only the basal cell layers were labeled in normal tissues. In some carcinomas, such as gastric carcinomas, the antibody demonstrated a stronger labeling than in normal tissues, especially on the membrane. None of 88 non-epithelial tumors and 20 cases of leukemia were labeled by the antibody (4).

In a study of 83 adenocarcinomas and 115 malignant mesotheliomas, 72/83 adenocarcinomas were labeled by the antibody whereas only 1/115 malignant mesotheliomas were labeled (1). In another study, 20/20 lung adenocarcinomas and 4/46 mesotheliomas were labeled. Of the 4 labeled mesotheliomas, the 2 showed a strictly focal labeling (6). In lymph nodes classified as tumor free by conventional histopathological staging, the antibody labeled micrometastatic tumor cells in 89 of 126 patients with completely resected esophageal carcinomas (2). In a study of 75 skin tumors, the antibody labeled 39/39 basal cell carcinomas, 0/23 squamous cell carcinomas, and showed some areas of staining in 13/13 basosquamous carcinomas (3).

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

References/ Bibliographie/ Literatur

1. Sheibani K, Shin SS, Kezirian J, Weiss LM. Ber-EP4 antibody as a discriminant in the differential diagnosis of malignant mesothelioma versus adenocarcinoma. Am J Surg Pathol 1991;15:779-84.
2. Hosch S, Kraus J, Scheunemann P, Izbicki JR, Schneider C, Schumacher U, et al. Malignant potential and cytogenetic characteristics of occult disseminated tumor cells in esophageal cancer. Cancer Res 2000;60:6836-40.
3. Beer TW, Shepherd P, Theaker JM. Ber EP4 and epithelial membrane antigen aid distinction of basal cell, squamous cell and basosquamous carcinomas of the skin. Histopathology 2000;37:218-23.
4. Latza U, Niedobitek G, Schwarting R, Nekarda H, Stein H. Ber-EP4: new monoclonal antibody which distinguishes epithelia from mesothelia. J Clin Pathol 1990;43:213-9.
5. Momburg F, Moldenhauer G, Hämmerling GJ, Möller P. Immunohistochemical study of the expression of a Mr 34,000 human epithelium-specific surface glycoprotein in normal and malignant tissues. Cancer Res 1987;47:2883-91.
6. Carella R, Deleonardi G, D'Errico A, Salerno A, Egarter-Vigl E, Seebacher C, et al. Immunohistochemical panels for differentiating epithelial malignant mesothelioma from lung adenocarcinoma. Am J Surg Pathol 2001;25:43-50.

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Agilent Technologies Singapore (International) Pte Ltd.