**Monoclonal Mouse
Anti-Human MITF
(microphthalmia transcription factor)
Clone D5**

**English
Code M3621**

Intended use

For In Vitro Diagnostic Use.

This antibody is intended for laboratory use to identify qualitatively by light microscopy MITF expressing cells in normal and neoplastic tissues using immunohistochemical (IHC) test methods. Positive results aid in the classification of melanoma.³ Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Synonym

Microphthalmia Transcription Factor

Summary and explanation

MITF (microphthalmia transcription factor) is a basic-helix-loop-helix-leucine-zipper (bHLHzip) transcription factor which regulates the development and survival of melanocytes and retinal pigment epithelium, and is also involved in transcription of pigmentation enzyme genes such as tyrosinase, TRP1 and TRP2.^{2,4} MITF has been shown to be phosphorylated by MAP kinase in response to c-Kit activation, resulting in upregulation of MITF transcriptional activity.¹ Mutations of the MITF gene are associated with the autosomal dominant hereditary deafness and pigmentation condition, Waardenburg Syndrome type 2A.⁵ Multiple isoforms of MITF exist including MITF-A, MITF-B, MITF-C, MITF-H, and MITF-M, which differ in the amino-terminal domain and in their expression patterns. The MITF-M isoform is restricted to the melanocyte cell lineage.⁴

Anti-MITF recognizes a nuclear protein which is expressed in the majority of primary and metastatic epithelioid malignant melanomas as well as in normal melanocytes, benign nevi and dysplastic nevi. A wide range of non-melanoma malignancies have been tested with anti-MITF, and only a small subset of tumors were found to be immunoreactive.^{6,7,8}

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Anti-MITF is provided in liquid form as tissue culture supernatant (containing fetal bovine serum) dialyzed against 0.05 mol/L Tris-HCl, pH 7.2, and 0.015 mol/L sodium azide. Contains stabilizing protein.

Clone: D5 Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

M3621 may be used at a dilution of 1:100 when performing IHC using the EnVision+™ detection system. These are guidelines only. Optimal antibody concentrations may vary depending on specimen and preparation method, and should be determined individually in each laboratory.

Immunogen

Histidine fusion protein of the amino-terminal *Taq-Sac* fragment of human *MITF* cDNA^{1,2}

Specificity

In Western blot analysis of human melanoma cell line, 501mel, anti-MITF demonstrated 54 and 60 kD bands corresponding to MITF variants differing in phosphorylation status. Further specificity was demonstrated in gel-mobility assays where anti-MITF produced a specific gel-mobility supershift with MITF but not with related bHLHzip factors TFEB, TFEC and TFE3.¹

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions. Suggested diluent(s) for IHC procedures:

Antibody Diluent (code S0809)

The following negative control is recommended for IHC procedures:

Mouse IgG₁ (code X0931)

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.^{9,10}
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin Sections

Anti-MITF can be used on formalin-fixed, paraffin-embedded tissue sections. Pretreatment of tissue with proteolytic enzymes is not recommended.

The deparaffinized tissue sections must be treated with heat prior to the IHC staining procedure. Heat-induced epitope retrieval (HIER) for this primary antibody requires immersion of tissue sections in a buffer solution and heating in a pressure cooker (121 °C). Use a 5 minute heating protocol in a pressure cooker set at 121 °C. Rinse well with buffer or deionized water following HIER. For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended. Target Retrieval Solution (code S1700) or 10x Concentrate (code S1699) is recommended.

Cryostat Sections And Cell Smears

Anti-MITF can be used for labelling acetone-fixed cryostat sections or fixed cell smears.

Staining procedure

Follow the recommended procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-MITF is nuclear.

Performance characteristics

Normal Tissues¹¹

Tissue Type (# tested)	Positively Staining Tissue Elements
Skin ^b	nuclei of melanocytes
Nevi ^b	nuclei of melanocytes
Dysplastic nevi ^b	nuclei of melanocytes
Adrenal (3) ¹¹	0/3
Bone Marrow (3) ¹¹	0/3
Brain, cerebellum (3) ¹¹	0/3
Brain, cerebrum (3) ¹¹	0/3
Breast (3) ¹¹	0/3
Cervix (3) ¹¹	0/3
Colon (3) ¹¹	0/3
Esophagus (3) ¹¹	0/3
Heart (3) ¹¹	1/3 nuclei of myocytes
Kidney (3) ¹¹	0/3
Liver (3) ¹¹	0/3
Lung (3) ¹¹	1/3 nuclei and cytoplasm of alveolar macrophages (weak)
Mesothelial Cells (3) ¹¹	0/3
Nerve, peripheral (3) ¹¹	0/3
Ovary (3) ¹¹	0/3
Pancreas (3) ¹¹	0/3
Parathyroid (3) ¹¹	0/3
Pituitary (3) ¹¹	1/3 cytoplasm (granular) of rare cells in neurohypophysis*
Prostate (3) ¹¹	0/3
Salivary Gland (3) ¹¹	0/3
Skeletal Muscle (3) ¹¹	0/3
Skin (3) ¹¹	3/3 nuclei of melanocytes
Small Intestine (3) ¹¹	0/3
Spleen (3) ¹¹	0/3
Stomach (3) ¹¹	0/3
Testis (3) ¹¹	0/3
Thymus (3) ¹¹	0/3
Thyroid (3)	0/3
Tonsil (3) ¹¹	0/3
Uterus (3) ¹¹	1/3 nuclei and cytoplasm of endometrial stromal cells

*only nuclear immunoreactivity was considered as positive

Abnormal Tissues^{5,6}

Tissue Type (# tested)	Positively Staining Tissue Elements
Malignant melanoma (76) ⁶	76/76 nuclei of tumor cells
In situ (19) ⁶	19/19 nuclei of tumor cells
Primary (50) ⁶	50/50 nuclei of tumor cells
Metastatic (7) ⁶	7/7 nuclei of tumor cells
Epithelioid melanoma (32) ⁷	29/32 nuclei of tumor cells
Melanoma (44) ⁸	44/44 nuclei of tumor cells
Ameloblastoma (1) ⁷	0/1
Angiosarcoma (2) ⁷	0/2
Breast carcinoma, invasive ductal cell (10) ⁶	0/10 nuclei of tumor cells 2/10 cytoplasm of tumor cells*
Breast carcinoma (8) ⁷	1/8 nuclei of tumor cells (multi-focal and weak)
Carcinoid (2) ⁷	0/2
Colon carcinoma (18) ⁷	0/18
Endometrial adenocarcinoma (10) ⁶	0/10
Endometrial stroma sarcoma (2) ⁷	0/2
Esophagus carcinoma (2) ⁷	0/2
Gastric carcinoma (3) ⁷	0/3
Leiomyosarcoma (5) ⁷	2/5 nuclei of tumor cells (focal)
Lung adenocarcinoma (6) ⁷	0/6
Lung small-cell carcinoma (5) ⁷	0/5
Lung squamous cell carcinoma (19) ^{6,7}	0/19
Malignant mixed mullerian tumor (1) ⁷	1/1 nuclei of tumor cells (weak)
Malignant peripheral nerve sheath tumor (3) ⁷	0/3
Nasopharyngeal carcinoma (1) ⁷	0/1
Osteosarcoma (1) ⁷	0/1
Ovarian carcinoma (1) ⁷	0/1
Pancreatic carcinoma (1) ⁷	0/1
Papillary thyroid (2) ⁷	0/2
Pheochromocytoma (1) ⁷	0/1
Prostate carcinoma (2) ⁷	0/2
Renal cell carcinoma (17) ⁷	2/17 nuclei of tumor cells
Rhabdomyosarcoma (1) ⁷	0/1
Squamous carcinoma of the vulva (10) ⁶	0/10
Synovial sarcoma (1) ⁷	0/1
Testicular cancer (10) ⁶	0/10
Thymoma (1) ⁷	0/1
Thyroid carcinoma (10) ⁶	0/10
Wilm's tumor (1) ⁷	0/1

*only nuclear immunoreactivity was considered as positive

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Bibliographie

Literaturangaben

- Hemesath TJ, Price ER, Takemoto C, Badalian T, Fisher DE. MAP kinase links the transcription factor Microphthalmia to c-Kit signalling in melanocytes. *Nature*. 1998; 391(6664):298
- Tachibana M, Perez-Jurado LA, Nakayama A, Hodgkinson CA, Li X, Schneider M, Miki T, Fex J, Francke U, Arnheiter H. Cloning of *MITF*, the human homolog of the mouse *microphthalmia* gene and assignment to chromosome 3p14.1-p12.3. *Hum Mol Genet*. 1994; 3(4):553
- Chang KL, Folpe AL. Diagnostic utility of microphthalmia transcription factor in malignant melanoma and other tumors. *Adv Anat Pathol*. 2001; 8(5):273
- Widlund HR, Fisher DE. Microphthalmia-associated transcription factor: a critical regulator of pigment cell development and survival. *Oncogene*. 2003; 22(20):3035
- Tassabehji M, Newton VE, Read AP. Waardenburg syndrome type 2 caused by mutations in the human microphthalmia (*MITF*) gene. *Nat Genet*. 1994; 8(3):251
- King R, Weilbaecher KN, McGill G, Cooley E, Mihm M, Fisher DE. Microphthalmia transcription factor. A sensitive and specific melanocyte marker for MelanomaDiagnosis. *Am J Pathol*. 1999; 155(3):731
- Granter SR, Weilbaecher KN, Quigley C, Fisher DE. Role for microphthalmia transcription factor in the diagnosis of metastatic malignant melanoma. *Appl Immunohistochem Mol Morphol*. 2002; 10(1):47

8. Dorvault CC, Weilbaeher KN, Yee H, Fisher DE, Chiriboga LA, Xu Y, Chhieng DC. Microphthalmia transcription factor: a sensitive and specific marker for malignant melanoma in cytologic specimens. *Cancer*. 2001; 93(5):337
9. Department of Health, Education and Welfare, National Institute for Occupational Safety and Health, Rockville, MD. "Procedures for the decontamination of plumbing systems containing copper and/or lead azides." DHHS (NIOSH) Publ. No. 78-127, Current 13. August 16, 1976
10. Centers for Disease Control Manual Guide – Safety Management, No. CDC-22, Atlanta, GA. "Decontamination of laboratory sink drains to remove azide salts." April 30, 1976
11. IHC003 Report On File

REF	Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i>
	Manufacturer Fabricant Hersteller	LOT	Batch code Code du lot Chargenbezeichnung
EC REP	Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU		IVD
		In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA

Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

EC **REP**

Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark

Tel +45 4485 9500
Fax +45 4485 9595

www.dako.com

PT0039/Rev C