## Monoclonal Mouse Anti-Human MUM1 Protein

Clone MUM1p

**Code No./ Code/ Code-Nr. M 7259**

Edition/ Edition/ Ausgabe 30.09.03

ENGLISH
---------

**Intended use**  
For in vitro diagnostic use.

Monoclonal Mouse Anti-Human MUM1 Protein, Clone MUM1p, is intended for use in immunocytochemistry. The antibody labels the MUM1 protein, which is expressed in a subset of B cells in the light zone of the germinal centre (representing late stages of B-cell differentiation), plasma cells, activated T cells, and a wide spectrum of related haematolymphoid neoplasms. Of non-haematolymphoid neoplasms only a proportion of melanomas are labelled. The antibody is of value, for example, for the subclassification of lymphoid malignancies (1, 2). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

**Synonyms for antigen**  
IRF4 (interferon regulatory factor 4), ICSAT (interferon consensus sequence binding protein for activated T cells), and PIP (PU.1 interaction partner) (3).

**Introduction**  
MUM1 protein (multiple myeloma oncogene-1) is a 50 kDa protein encoded by the *MUM1* gene that was originally identified because of its involvement in the rare t(6;14) (p25;q32) translocation observed in multiple myeloma (MM), causing juxtaposition of the *MUM1* gene to the Ig heavy chain locus (4, 5). MUM1-protein-deficient mice are unable to form germinal centres, lack plasma cells in the spleen and lamina propria, and exhibit a profound reduction of serum immunoglobulin levels and an inability to mount detectable antibody responses or to generate T-lymphocyte cytotoxic or antitumour responses (1). Expression of MUM1 protein is independent from the t(6;14) (p25;q32) translocation and has been detected in multiple myelomas, lymphoplasmacytic lymphomas (LPL), diffuse large B-cell lymphomas (DLBCL), and activated T cells. MUM1 protein is not specific for plasmacytic differentiation, but it may be added to the panel of phenotypic markers available for characterization of B-cell lymphoma histogenesis, in addition to being a useful marker in the subclassification of lymphoid malignancies, e.g. in B-CLL (1-3, 6).

MUM1 protein expression is restricted to cells in the lymphocytic and melanocytic lineages (5).

**Reagent provided**  
Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN<sub>3</sub>.  
Clone: MUM1p (1). Iso­type: IgG1, kappa.  
Mouse IgG concentration: See label on vial.

**Immunogen**  
Recombinant GST-MUM1 fusion protein (1).

**Specificity**  
In Western blotting of pHeBo-CMV-MUM1-HA-transfected HeLa cell lysate, the antibody labels a band of 52 kDa corresponding to MUM1-HA protein, whereas no labelling is observed in non-transfected HeLa cell lysate. In Western blotting of lysates of IM9 (myeloma cell line), unfractionated tonsil cell suspensions, and PHA-stimulated T cells, the antibody labels a band of 50 kDa corresponding to MUM1 protein. No labelling is observed with lysates of U937 cells (myeloid cell line), and unstimulated T cells (1).

SDS-PAGE analysis of immunoprecipitates formed between lysate of the IM9 myeloma cell line and the antibody shows reaction with a 50 kDa protein corresponding to MUM1 protein (1).

**Precautions**  
1. For professional users.  
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.  
3. As with any product derived from biological sources, proper handling procedures should be used.

**Storage**  
Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Specimen preparation**  
Paraffin sections: The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin or B5 (1). Pre-treatment of tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, DakoCytomation Target Retrieval Solution, pH 9.0, code No. S 2368, or DakoCytomation Target Retrieval Solution, High pH, code No. S 3308. Less optimal results are obtained with 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with DakoCytomation Proteinase K, code No. S 3020, was found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.

Frozen sections and cell preparations: The antibody can be used for labelling acetone-fixed, frozen sections. However, the labelling is usually weaker than in paraffin sections (1).

**Staining procedure**  
Dilution: Mono­clonal Mouse Anti-Human MUM1 Protein, code No. M 7259, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, pH 9.0, code No. S 2368, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen.

Visualization: DAKO LSAB™+/ HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.

Automation: The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.

(104293-002)	M 7259/EFG/CE/30.09.03 p. 1/4
DakoCytomation Denmark A/S · Produktionsvej 42 · DK-2600 Glostrup · Denmark · Tel. +45 44 85 95 00 · Fax +45 44 85 95 95 · CVR No. 33 21 13 17	

### Performance characteristics

Cells labelled by the antibody display staining predominantly confined to the nucleus, although weak to moderate staining of the cytoplasm is also present in most cases with nuclear positivity (2).

Normal tissues: In tonsil and spleen, the antibody strongly labels plasma cells and between 3 and 10% of the germinal centre (GC) B cells, mainly located in the light zone. Rare MUM1 protein-positive cells in the mantle zone represent B cells exiting GC and transiting through the follicle mantle zone. The antibody also labels 1-5% of the T cells present in the GC and the interfollicular area. No labelling is observed in tonsil macrophages, follicular dendritic cells, endothelia and epithelia. Thymus cortex, thymus medulla, Hassal's corpuscles, and bone marrow erythroid and myeloid precursors, megakaryocytes, osteoblasts and osteoclasts are also MUM1 protein-negative (1). In an additional study of 142 normal, non-lymphoid tissues representing 33 tissue types, all tissues were found to be negative with the antibody (2).

Abnormal tissues: In a study of a large number (N=210) and variety (18 types) of human haematolymphoid malignancies, the antibody labelled a wide spectrum of B-, T-, and natural killer (NK)-cell lymphomas. All types of haematolymphoid neoplasms tested were labelled, except for Burkitt's lymphoma, mast cell, and histiocytic disorders. A total of 56% of T- and NK-cell lymphomas, and a total of 35% B-cell neoplasms were labelled, with a strong labelling of MM (100%), DLBCL (51%), marginal-zone lymphoma (58%), LPL (50%), post-transplant lymphoproliferative disorders (67%), and Hodgkin's disease (100%). Of 122 follicle centre lymphomas, 23% were labelled by the antibody, i.e. 79% of grade III and 21% of grade I+II. Lesional cells labelled by the antibody did not necessarily show plasmacytic differentiation, suggesting that MUM1 protein expression may precede the expression of CD138 and VS38 proteins that are more specific markers for plasmacytic differentiation. Of 731 non-lymphoid neoplasms representing 20 different types, only 5/22 melanomas were labelled by the antibody (2).

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