**Monoclonal Mouse
Anti-Human
Serotonin**
Clone 5HT-H209
Code M0758

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Serotonin, Clone 5HT-H209, is intended for use in immunocytochemistry. The antibody labels cells expressing serotonin. The antibody may be a useful tool for characterization of carcinoids (1-3). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Synonym for antigen	5-hydroxytryptamine, 5HT (4).
Summary and explanation	Carcinoid tumours are the most commonly occurring endocrine tumours. Although generally found in the gut wall, they also occur in the pancreas, rectum, ovary and lung (5). Carcinoid tumours typically contain numerous membrane-bound neurosecretory granules, which are composed of a variety of hormones and biogenic amines. One of the best characterized of these substances is serotonin. Serotonin is synthesized from its precursor, 5-hydroxytryptophan, by the enzyme aromatic acid decarboxylase. Serotonin is subsequently metabolised by monoamine oxidase to 5-hydroxyindoleacetic acid, which is excreted in the urine (6). In addition to serotonin, carcinoid tumours have been found to secrete corticotropin, histamine, dopamine, substance P, neurotensin, prostaglandins, and kallikrein (6). Among endocrine tumours of the gastrointestinal tract, carcinoid tumours arising from the midgut secrete serotonin, while carcinoid tumours arising from the hindgut rarely secrete serotonin (5).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> 5HT-H209. <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> See label on vial.
Immunogen	5-hydroxytryptamine hydrochloride.
Specificity	In immunocytochemistry the antibody labels entero-chromaffine (EC) cells in colon. As demonstrated by immunocytochemistry, the antibody cross-reacts with the serotonin-equivalent molecule in guinea pig (4), mouse (7) and rat (8).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with Dako Proteinase K, code S3020, or heat-induced epitope retrieval is required. For heat-induced epitope retrieval, optimal results are obtained with Dako Target Retrieval Solution, code S1700. Less optimal results are obtained with Dako Target Retrieval Solution, High pH, code S3308, Dako Target Retrieval Solution, pH 9.0, code S2368, or 10 mmol/L citrate buffer, pH 6.0. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections (4) or cell smears (2).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human Serotonin, code M0758, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human colon and using 5 minutes proteolytic epitope retrieval with Dako Proteinase K, code S3020, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG1, code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code S0809. Positive and negative controls should be run simultaneously with patient specimen.

Visualization: Dako LSAB™+/HRP kit, code K0679, and Dako EnVision™+/HRP kits, codes K4004 and K4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code K0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.

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
1. Torske KR, Thompson LDR. Adenoma versus carcinoid tumor of the middle ear: a study of 48 cases and review of the literature. Mod Pathol 2002;15:543-55.
2. Menezes G, Wakely P. Aspiration cytopathology of middle-ear neuroendocrine carcinoma. Diagn Cytopathol 2001;25:168-71.
3. Granberg D, Wilander E, Öberg K, Skogseid B. Prognostic markers in patients with typical bronchial carcinoid tumors. J Clin Endocrinol Metab 2000;85:3425-30.
4. Suzuki A, Naruse S, Kitagawa M, Ishiguro H, Yoshikawa T, Ko SBH, et al. 5-hydroxytryptamine strongly inhibits fluid secretion in guinea pig pancreatic duct cells. J Clin Invest 2001;108:749-56.
5. Perry RR, Vinik AI. Endocrine tumors of the gastrointestinal tract. Annu Rev Med 1996;47:57-68.
6. Kulke MH, Mayer RJ. Carcinoid tumors. N Engl J Med 1999;340:858-68.
7. Batt J, Asa S, Fladd C, Rotin D. Pituitary, pancreatic and gut neuroendocrine defects in protein tyrosine phosphatase-sigma-deficient mice. Mol Endocrinol 2002;16:155-69.
8. Date Y, Kojima M, Hosoda H, Sawaguchi A, Mondal MS, Suganuma T, et al. Ghrelin, a novel growth hormone-releasing acylated peptide, is synthesized in a distinct endocrine cell type in the gastrointestinal tracts of rats and humans. Endocrinology 2000;141:4255-61.

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