Monoclonal Mouse Anti-Human CD45, Leucocyte Common Antigen
Clones 2B11 + PD7/26
Code M0701

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human CD45, Leucocyte Common Antigen, Clones 2B11 + PD7/26, is intended for use in immunohistochemistry. The antibody labels CD45 in both normal and neoplastic cells, and is a useful tool for identifying tumour cells of lymphoid origin (1-3). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Synonyms for antigen	T200, Ly-5.
Summary and explanation	CD45 is a transmembrane glycoprotein expressed on most nucleated cells of haematopoietic origin. CD45, encoded by a single gene mapped to chromosome 1, has various isoforms based on differential splicing of exons 4, 5 and 6. On human leucocytes, five different isoforms of CD45, named ABC, AB, BC, B and 0, have been identified. These isoforms are recognized by CD45RA, CD45RB, CD45RC and CD45R0 antibodies. The isoforms range in Mr from 180000 to 220000. All the CD45 isoforms share the same intracellular segment, which has been shown to have tyrosine phosphatase activity. Various leucocytes express characteristic CD45 isoforms, thus T cells express CD45 isoforms corresponding to their development and activation, B cells predominantly express the ABC isoform, and monocytes and dendritic cells predominantly express the B and 0 isoforms. Granulocytes principally express only the B and 0 isoforms (4).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> 2B11 (1) and PD7/26 (1). <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial. The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.
Immunogen	2B11: Isolated neoplastic cells from a case of T-cell lymphoma/leukaemia (1). PD7/26: Human peripheral blood lymphocytes maintained in T-cell growth factor (1).
Specificity	Anti-CD45 is a mixture of two monoclonal antibodies, clones 2B11 and PD7/26, directed against different epitopes. Clone 2B11 was clustered as anti-CD45 at the Third International Workshop and Conference on Human Leucocyte Differentiation Antigens, held in Oxford in 1986 and reacts with all the known isotypes of the CD45 family (5). Clone PD7/26 was clustered as anti-CD45RB at the Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens, held in Boston in 1993 (6).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin, B5 (2), or Bouin's (1). Pretreatment of tissues with heat-induced epitope retrieval is recommended. Optimal results are obtained with Dako Target Retrieval Solution, code S1700, Dako Target Retrieval Solution, High pH, code S3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone-fixed, frozen sections (1), and cell smears.
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human CD45, Leucocyte Common Antigen, code M0701, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval solution, code S1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG1, code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code S0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> Dako LSAB™+HRP kit, code K0679, and Dako EnVision™+HRP kits, codes K4004 and K4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code K0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunohistochemical staining using automated platforms, such as the Dako Autostainer.
Product-specific limitations	Labelling of the surface membrane of mammary ductal cells in specimens of fibrocystic disease and fibroadenoma has been reported for the antibody. This expression, however, is not considered to present a differential diagnostic problem that might lead to misidentification (7).
Performance characteristics	Cells labelled by the antibody predominantly display staining of the cell membrane, but cytoplasmic staining may also occur.

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
- Warnke RA, Gatter KC, Falini B, Hildreth P, Woolston R-E, Pulford K, et al. Diagnosis of human lymphoma with monoclonal antileukocyte antibodies. N Engl J Med 1983;309:1275-81.
- Kurtin PJ, Pinkus GS. Leukocyte common antigen - a diagnostic discriminant between hematopoietic and nonhematopoietic neoplasms in paraffin sections using monoclonal antibodies: Correlation with immunologic studies and ultrastructural localization. Hum Pathol 1985;16:353.
- Michie SA, Spagnolo DV, Dunn KA, Warnke RA, Rouse RV. A panel approach to the evaluation of the sensitivity and specificity of antibodies for the diagnosis of routinely processed histologically undifferentiated human neoplasms. Am J Clin Pathol 1987;88:457-62.
- Sewell WA, Cooley MA, Hegen M. NL6. CD45 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 499-502.
- Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1, and leucocyte common antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leukocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 788-803.
- Morimoto C. T18. CD45 cluster report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leukocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 386-9.
- Herman GE and Elfont E. Aberrant CD45 (leukocyte common antigen) staining of non-malignant breast lesions in zinc formalin fixed tissue. J Histotechnol 1993;16:151-3.
- Hall PA, d’Ardenne AJ, Stansfeld AG. Paraffin section immunohistochemistry. I. Non-Hodgkin’s lymphoma. Histopathol 1988;13:149-60.

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