Monoclonal Mouse Anti-Human CD45R0
Clone UCHL1
Code No./ Code/ Code-Nr. M 0742
Edition/ Ausgabe 20.01.03

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human CD45R0, Clone UCHL1, is intended for use in immunocytochemistry. The antibody labels CD45R0 in both normal and neoplastic cells, and is a useful tool for identifying T-cell lymphomas, and for the differentiation of low-grade B-cell from T-cell lymphomas (1, 2). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Introduction	CD45 is a transmembrane glycoprotein expressed on most nucleated cells of haematopoietic origin. CD45, encoded by a single gene mapped to chromosome 1, has various isoforms based on differential splicing of exons 4, 5 and 6. On human leucocytes, five different isoforms of CD45, named ABC, AB, BC, B and 0, have been identified. These isoforms are recognized by CD45RA, CD45RB, CD45RC and CD45R0 antibodies. The Mr of the CD45R0 isoform is 180 000, and the gene corresponding to the CD45R0 protein contains none (hence the 0) of the differentially spliced exons 4, 5 and 6. All the CD45 isoforms share the same intracellular segment, which has been shown to have tyrosine phosphatase activity. Various leucocytes express characteristic CD45 isoforms, thus T cells express CD45 isoforms corresponding to their development and activation. B cells predominantly express the ABC isoform, and monocytes and dendritic cells predominantly express the B and 0 isoforms. Granulocytes principally express only the B and 0 isoforms (3).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> Clone UCHL1 (4). <u>Isotype:</u> IgG2a, kappa. <u>Mouse IgG concentration:</u> see label on vial.
Immunogen	IL-2 dependent T-cell line, CA1 (4).
Specificity	The antibody was clustered as anti-CD45R0 at the Fourth (5) and Fifth (6) International Workshops and Conferences on Human Leucocyte Differentiation Antigens. In Western blotting of lysate of peripheral blood mononuclear cells, the antibody was found to react specifically with a 180 kDa band corresponding to CD45R0. In leucocyte common antigen (LCA) transfected cell lines, expressing the ABC, the AB, the BC, the B or the 0 CD45 isoform only, the antibody reacted solely with the transfectant expressing CD45R0 (7).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Bouin's and acidic formalin (pH < 3.1) are unsuitable (4). Pre-treatment of tissues with heat-induced epitope retrieval is recommended. Optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone-fixed, frozen sections (1, 7), cell smears (1), and cytospins (7).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human CD45R0, code No. M 0742, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval solution, code No. S 1700, pH 6.0, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG2a, code No. X 0943, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.
Performance characteristics	Cells labelled by the antibody generally display a staining confined to the cell membrane. <u>Normal tissues:</u> The antibody labels most thymocytes, a subpopulation of resting T cells within both the CD4 and CD8 subsets, and mature activated T cells. Additionally, granulocytes and monocytes are labelled, whereas normal B cells and NK cells are negative (4). In a study covering a large panel of normal tissues (1), the antibody labelled lymphocytes in the T-cell areas of normal tonsil, spleen, and normal and reactive lymph nodes. In the thymus 90% of cortical and about 50% of medullary thymocytes were labelled. In contrast, the large cortical thymic blasts were negative. Membrane labelling in mature myeloid cells and about 20% of macrophages was a constant finding. Diffuse, and rarely
(102420-001)	M 0742/EFG/CE/20.01.03 p. 1/4

strong, cytoplasmic labelling was observed in glandular epithelia, squamous epithelium transitional epithelium, hepatocytes, syncytiotrophoblasts and all smooth muscle.

Abnormal tissues: Of T-cell tumours, the antibody labelled 100% of cases of *Mycosis fungoides* (n=10), 83% of cases of peripheral T-cell lymphoma (n=25), 78% of cases of T acute lymphoblastic lymphoma (n=9) and 100% of malignant histiocytosis of the intestine (n=13). 2/2 cases of true histiocytic lymphoma were also positive. No labelling was observed in a large range of B-cell lymphomas (n=62), but some large B-cell lymphomas of centroblastic and immunoblastic types showed diffuse cytoplasmic labelling. In Hodgkin's lymphoma (n=16) Reed-Stenberg cells showed positive cytoplasmic and negative membrane labelling. In 4/4 cases of granulocytic sarcoma, the antibody labelled mature myeloid cells only (1). With few exceptions, the antibody labelled tumour cells in 28 cases of cutaneous malignant T-cell lymphoma. In 85 biopsy specimens from inflammatory skin diseases, most reactive T cells were positive with the antibody. Cells in cutaneous B-cell lymphomas were consistently negative, as were tumour cells in 125 cases of various non-lymphatic skin tumours (8).

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
References/ Références/ Literatur												
<ol style="list-style-type: none">Norton AJ, Ramsay AD, Smith SH, Beverley PCL, Isaacson PG. Monoclonal antibody (UCHL1) that recognises normal and neoplastic T cells in routinely fixed tissues. J Clin Pathol 1986;39:399-405. Davey FR, Elghetany MT, Kurec AS. Immunophenotyping of hematologic neoplasms in paraffin-embedded tissue sections. Am J Clin Pathol 1990;93. Suppl 1:S17-S26. Sewell WA, Cooley MA, Hegen M. NL6. CD45 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 499-502. Smith SH, Brown MH, Rowe D, Callard RE, Beverley PCL. Functional subsets of human helper-inducer cells defined by a new monoclonal antibody, UCHL1. Immunology 1986;58:63-70. Schmidt RE. Non-lineage/natural killer section report: new and previously defined clusters. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leukocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 517-42. Morimoto C. T18. CD45 cluster report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leukocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 386-9. Poppema S, Lai R, Visser L. Monoclonal Antibody OPD4 is reactive with CD45R0, but differs from UCHL1 by the absence of monocyte reactivity. Am J Pathol 1991;139:725-9. Sterry W, Hauschild A. Use of monoclonal antibodies (UCHL1, Ki-B3) against T and B cell antigens in routine paraffin-embedded skin biopsy specimens. J Am Acad Dermatol 1989;21:98-107.												
Explanation of symbols/ Légende des symboles/ Erläuterung der Symbole												
<table> <tbody><tr> <td>REF</td> <td>Catalogue number Référence du catalogue Bestellnummer</td> <td> <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p></td> <td> <p>Manufacturer Fabricant Hersteller</p></td></tr> <tr> <td>IVD</td> <td>In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum</td> <td>LOT</td> <td>Batch code Code du Lot Chargenbezeichnung</td></tr> <tr> <td></td> <td>Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten</td> <td></td> <td>Use by Utiliser jusque Verwendbar bis</td></tr> </tbody></table>	REF	Catalogue number Référence du catalogue Bestellnummer	 <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p>	 <p>Manufacturer Fabricant Hersteller</p>	IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	LOT	Batch code Code du Lot Chargenbezeichnung		Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten		Use by Utiliser jusque Verwendbar bis
REF	Catalogue number Référence du catalogue Bestellnummer	 <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p>	 <p>Manufacturer Fabricant Hersteller</p>									
IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	LOT	Batch code Code du Lot Chargenbezeichnung									
	Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten		Use by Utiliser jusque Verwendbar bis									

(102420-001)	M 0742/EFG/CE/20.01.03 p. 4/4
DakoCytomation Denmark A/S · Produktionsvej 42 · DK-2600 Glostrup · Denmark · Tel. +45 44 85 95 00 · Fax +45 44 85 95 95 · CVR No. 33 21 13 17	