**Monoclonal Mouse  
Anti-Human  
Wilms' Tumor 1 (WT1) Protein  
Clone 6F-H2**

**ENGLISH**  
**Code M3561**

**Intended use**

For In Vitro Diagnostic Use.

**Summary and explanation**

*WT1* is a gene involved in the induction of Wilms' tumor, a pediatric renal malignancy. The *WT1* gene, located on chromosome 11p13, is inactivated in 5 to 10% of sporadic Wilms' tumors and nearly 100% of Denys-Drash patients, a syndrome associated with genitourinary abnormalities and Wilms' tumor. *WT1* encodes a zinc finger transcription factor which recognizes the early growth response (EGR-1) consensus sequence found in promoters of growth factor genes. The protein encoded by *WT1* regulates transcription of other genes and can function both as a transcriptional activator and repressor. *WT1* has been demonstrated to repress transcription of various growth-related genes such as platelet-derived growth factor A chain (*PDGF-A*) and insulin-like growth factor (*IGF*).<sup>3-5</sup> The transcriptional activity of this protein, however, can be modulated through physical interactions between *WT1* and p53. In the absence of wild type p53, *WT1* has been demonstrated to function as a transcriptional activator of the *EGR-1* promoter in transfection assays.<sup>6</sup>

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

**Reagent provided**

Anti-human *WT1*, 6F-H2 is a mouse monoclonal antibody supplied in liquid form as tissue culture supernatant (containing fetal bovine serum) dialyzed against 0.05 mol/L Tris-HCl, pH 7.2, and 0.015 mol/L sodium azide.

Clone: 6F-H2<sup>1,2</sup> Isotype: IgG<sub>1</sub>, kappa  
Mouse IgG concentration mg/L: See label on vial.

Anti-human *WT1* may be used at a dilution of 1:50 in the LSAB<sup>TM</sup>+ method determined on formalin-fixed, paraffin-embedded tissue. These are guidelines only; optimal dilutions should be determined by the individual laboratory.

**Immunogen**

Truncated human *WT1* protein corresponding to N-terminal amino acids 1–181<sup>1,2</sup>

**Specificity**

Alternate splicing of the *WT1* transcript gene results in four distinct mRNA species which are present in different amounts. The predicted molecular weight of the *WT1* protein ranges from 46 to 49 kD depending on the mRNA species however, the *WT1* proteins derived from cells migrate at 50 to 55 kD in SDS-PAGE.<sup>4</sup> Monoclonal mouse anti-human *WT1* (anti-*WT1*) recognizes an epitope found in the amino terminal 84 amino acids of *WT1*. Anti-*WT1* reacts with all isoforms of the full-length *WT1* and also identifies *WT1* lacking exon 2, frequently found in subsets of sporadic Wilms' tumors. In immunoprecipitation assays, anti-*WT1* has been shown to recognize full-length and *in vitro* translated *WT1*. Anti-*WT1* also detects full-length denatured *WT1* (55 kD) in Western blots of *WT1* baculovirus-infected cell lysates.<sup>1</sup>

**Materials required, but not supplied**

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions.

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN<sub>3</sub> may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

**Storage**

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

## Specimen preparation

### Paraffin Sections

Anti-human WT1 can be used on formalin-fixed, paraffin-embedded tissue sections. Tissue sections should be pretreated with a 0.4% pepsin solution in 0.2 N HCl for 30 minutes at 37°C. For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended.

### Cryostat Sections and Cell Smears

Anti-human WT1 can also be used to label cryostat sections or cell smears.

## Staining procedure

Follow the recommended procedure for the detection system selected.

## Staining interpretation

The cellular staining pattern for anti-WT1 is nuclear and/or cytoplasmic.

## Performance characteristics

### Normal Cells

WT1 mRNA is strongly expressed in mesenchymally derived tissue. During embryonic development, WT1 expression has been observed in human kidney, spleen, and gonadal ridge mesoderm and the mesothelial lining of the coelomic cavity and the organs contained within.<sup>7</sup> WT1 gene expression in adult human tissues has been observed in sertoli cells of testes, decidual cells of the uterus and in granulosa cells of the ovary. A number of other adult tissues were shown to contain WT1 transcripts including pleura, spleen, heart, foreskin, endometrium, and glomerular epithelium of the kidney.<sup>7,8</sup> Antibodies to WT1, including clone 6F-H2, have been demonstrated to stain the cytoplasm of the following normal cells in some tissue specimens. Cytoplasmic staining of blood vessels and connective tissue of lung alveoli were observed to stain positively.<sup>8</sup> Granular cytoplasmic staining has also been reported in mature myelocytic cells.<sup>11</sup>

### Tumor Cells

High levels of WT1 expression have been demonstrated in the epithelial and blastemal components of Wilms' tumors, whereas stromal elements were found to be expressed at very low levels.<sup>10</sup> The majority of malignant mesotheliomas were found to be strongly express WT1 message and also the gene product. The WT1 protein was immunolocalized, using anti-WT1 (clone 6F-H2), to the nuclei of malignant mesothelioma cells.<sup>8,9</sup> WT1 gene and protein expression have also been demonstrated in the majority of acute leukemias but not in cells from chronic myelogenous leukemia, except when in blast crisis. Anti-WT1 immunostained the nuclei of leukemia blast cells but was unreactive with normal mononuclear cells and normal peripheral CD34+ hematopoietic progenitor cells.<sup>11,11</sup> WT1 mRNA expression has been reported in a number of other tumors including melanomas, ovarian cancer and sex cord-stromal tumors.<sup>11-14</sup> In addition to specific nuclear immunoreactivity, antibodies to WT1 (including clone 6F-H2) have been reported to stain the cytoplasm of tumor cells in some cases of adenocarcinoma and the desmoplastic stroma and basement membrane of some carcinoma specimens. Perinuclear staining of tumor cells in malignant mesothelioma has also been observed.<sup>9</sup> This cytoplasmic staining may represent cross-reactivity with an epitope unrelated to WT1.<sup>9,11</sup>

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
### Références

### Literatur

1. Rauscher JF, Morris JF, Fredericks WJ, Lopez-Guisa J, Balakrishnan C, Jost M, Herlyn M, Rodeck U. Characterization of monoclonal antibodies directed to the amino-terminus of the WT1, Wilms' tumor suppressor protein. *Hybridoma* 1998; 17:191
2. Menssen HD, et al. Presence of Wilms' tumor (*wt1*) transcripts and the WT1 nuclear protein in the majority of human acute leukemias. *Leukemia* 1995; 9:1060
3. Coppes MJ, et al. The role of *WT1* in Wilms tumorigenesis. *FASEB J* 1993; 7:886
4. Haber DA and Buckler AJ. *WT1*: A novel tumor suppressor gene inactivated in Wilms' tumor. *New Biol* 1992; 4(2):97
5. Rauscher FJ. The WT1 Wilms tumor gene product: A developmentally regulated transcription factor in the kidney that functions as a tumor suppressor. *FASEB J* 1993; 7:896
6. Maheswaran S, et al. Physical and functional interaction between WT1 and p53 proteins. *Proc Natl Acad Sci USA* 1993; 90:5100
7. Pritchard-Jones K, et al. The candidate Wilms' tumour gene is involved in genitourinary development. *Nature* 1990; 346:194
8. Amin KM, et al. Wilms' tumor 1 susceptibility (*WT1*) gene products are selectively expressed in malignant mesothelioma. *Amer J Pathol* 1995; 146(2):344
9. Kumar-Singh S, et al. *WT1* mutation in malignant mesothelioma and WT1 immunoreactivity in relation to *p53* and growth factor receptor expression, cell-type transition, and prognosis. *J Pathol* 1997; 181:67
10. Varansi R, et al. Fine structure analysis of the *WT1* gene in sporadic Wilms tumors. *Proc Natl Acad Sci USA* 1994; 91:3554
11. Menssen HD, et al. Detection by monoclonal antibodies of the Wilms' tumor (*WT1*) nuclear protein in patients with acute leukemia. *Intl J Canc* 1997; 70:518
12. Bruening W, et al. Analysis of the 11p13 Wilms' tumor suppressor gene (*WT1*) in ovarian tumors. *Canc Invest* 1993; 11(4):393
13. Rodeck U, et al. Expression of the *WT1* Wilms' tumor gene by normal and malignant human melanocytes. *Intl J Canc* 1994; 59:78
14. Coppes MJ, et al. Analysis of *WT1* in granulosa cell and other sex cord-stromal tumors. *Canc Res* 1993; 53:2712

 Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten
 Manufacturer Fabricant Hersteller	 Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis
 Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	 In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	

  
Dako North America, Inc.  
6392 Via Real  
Carpinteria, California 93013 USA

  
Dako Denmark A/S  
Produktionsvej 42  
DK-2600 Glostrup Denmark

Tel 805 566 6655  
Fax 805 566 6688  
Technical Support 800 424 0021  
Customer Service 800 235 5763

Tel +45 4485 9500  
Fax +45 4485 9595

www.dako.com

PT0039/Rev C