**Monoclonal Mouse
Anti-Human
Amyloid A**
Clone mc1
Code No./ Code/ Code-Nr. M 0759
Edition/ Ausgabe 20.01.03

ENGLISH

Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Amyloid A, Clone mc1, is intended for use in immunocytochemistry. The antibody labels amyloid A (AA) in tissues, and aids in the identification and classification of AA-amyloidosis. Immunocytochemical staining using Dako Monoclonal Mouse Anti-Human Amyloid A in combination with Congo Red staining (1), or particularly Congo Red fluorescence (2), is much more sensitive than Congo Red staining alone. Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Introduction	Amyloidosis is a group of diseases that have in common the extracellular deposition of fibrillar proteins with a specific biochemical conformation known as β -pleated sheets. Approximately 20 different precursor proteins that may be deposited as amyloid fibrils have been identified (3). Amyloid proteins deposited in the tissues can be identified by Congo Red staining, and chemically classified via amino acid sequence studies, by immunochemistry, or via immunocytochemistry (1). Amyloid A (AA) is an extracellular deposited insoluble fibrillar protein, highly resistant to proteolytic degradation, and produced from the precursor protein, serum amyloid A (SAA) (3). Before the therapy of a suspected amyloid disease can be planned, both the presence of amyloid and its chemical origin must be known (2).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCL, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> mc1 (4, 5). <u>Isotype:</u> IgG2a, kappa. <u>Mouse IgG concentration:</u> see label on vial.
Immunogen	An equal mixture of human amyloid A coupled to horseradish peroxidase and human amyloid A coupled to high molecular weight kininogen (5).
Specificity	In micro-ELISA, the antibody reacts with amyloid A and the serum precursor of amyloid A indicating cross-reactivity among amyloid A protein and the serum precursor of amyloid A. In contrast, no reactivity to other non-AA amyloid fibril proteins or human serum proteins, such as albumin, transferrin, and IgG is seen (4, 5). In immunocytochemistry, the antibody labels tissues from AA patients, but shows no reactivity with a host of unknown antigens in tissue sections of various organs (5), nor with amyloid types A κ , A λ , and amyloid fibril proteins in familial amyloid polyneuropathy and senile cardiovascular amyloidosis, respectively. Neither are Alzheimer's amyloid plaques, amyloid in microangiopathy, lichen amyloidosis, and senile islets of Langerhans labelled by the antibody (5-7).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with proteinase K or heat-induced epitope retrieval is recommended. For heat-induced epitope retrieval, optimal results are obtained with Dako Target Retrieval Solution, code No. S 1700, Dako Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections (5).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human Amyloid A, code No. M 0759, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human kidney from a patient with amyloidosis A and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG2a, code No. X 0943, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the Dako Autostainer.
Product-specific limitations	Cross-reactivity with the serum precursor of amyloid A has been observed (5, 6).
Performance characteristics	Amyloid A labelled by the antibody has an extracellular localization in the majority of cases.

Normal tissues: The antibody does not label normal tissue (5).

Abnormal tissues: The antibody labelled amyloid A in renal and rectal biopsies of patients with inflammatory pediatric disease (1). Tissues from patients with a clinical diagnosis of rheumatoid arthritis, sporadic Muckle-Wells syndrome, idiopathic polyneuritis, idiopathic amyloidosis, familial Mediterranean fever, and Still's syndrome were also labelled (5).

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
1. Linke RP, Gärtner HV, Michels H. High-sensitivity diagnosis of AA amyloidosis using Congo Red and immunohistochemistry detects missed amyloid deposits. *J Histochem Cytochem* 1995;43:863-9.
2. Linke RP. Highly sensitive diagnosis of amyloid and various amyloid syndromes using Congo red fluorescence. *Virchows Arch* 2000;436:439-48.
3. Sipe JD, Cohen AS. History of the amyloid fibril [review]. *J Struct Biol* 2000;130:88-98.
4. Linke RP. Identification of amyloid protein AA with a monoclonal antibody. *Blut* 1982;45:407-9.
5. Linke RP. Monoclonal antibodies against amyloid fibril protein AA. Production, specificity, and use for immunohistochemical localization and classification of AA-type amyloidosis. *J Histochem Cytochem* 1984;32:322-8.
6. Linke RP. Identification of AA-type amyloid in tissue sections using monoclonal antibodies. In: Peeters H, editor. *Protides Biol Fluids*. Oxford: Pergamon Press; 1983. p. 835-8.
7. Linke RP, Nathrath WBJ, Eulitz M. Classification of amyloid syndromes from tissue sections using antibodies against various amyloid fibril proteins: report of 142 cases. In: Glenner CG, Osseman EF, Benditt EP, Calkins E, Cohen AS, Zucker-Franklin D, editors. *Amyloidosis*. Amsterdam: Plenum Publishing Corporation; 1986. p. 599-605.

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