INSULIN

CE

Polyclonal Guinea Pig Anti-Insulin

ENGLISH
Code A0564

Intended use

For In Vitro Diagnostic Use.

Polyclonal guinea pig anti-insulin is intended for laboratory use to identify qualitatively by light microscopy insulin positive cells in normal and neoplastic tissues using immunohistochemical (IHC) test methods. The clinical interpretation of any positive staining or its absence should be complemented by morphological and histological studies with proper controls. Evaluations should be made within the context of the patient's clinical history and other diagnostic tests by a qualified individual.

Summary and explanation

Insulin is one of seven known polypeptide hormones produced in the pancreas. Insulin, secreted by B-cells of the islets of Langerhans, participates in glucose utilization, protein synthesis and in the formation and storage of neutral lipids.¹

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Guinea pig anti-insulin is an unfractionated guinea pig antiserum in 0.05 mol/L Tris-HCl, pH 7.6, and 0.015 mol/L sodium azide.

Protein concentration g/L: See label on vial.

Guinea pig anti-insulin may be used at a dilution of 1:50 to 1:150 in the LSAB method, determined on formalin-fixed, paraffin-embedded tissue. These are guidelines only; optimal dilutions should be determined by the individual laboratory.

Immunogen

Porcine pancreatic insulin

Specificity

The antibody cross-reacts with insulin from several mammalian species. Specificity as determined by radioimmunoassay was 100% for human insulin, 100% for porcine insulin and less than 0.05% for glucagon and human growth hormone. This product has been optimized for use on human tissues.

For more information on insulin detection and pancreatic tumors, please see references 3–8.

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions.

Precautions

1. For professional users.
2. This product contains sodium azide (NaN_3), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN_3 may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin Sections

Guinea pig anti-insulin can be used on formalin-fixed, paraffin-embedded tissue sections. Pretreatment of tissue with proteolytic enzymes is not required.

Staining procedure

Follow the recommended procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-insulin is cytoplasmic.

Performance characteristics

Normal Tissues

This antibody labels the cytoplasm of insulin-producing B-cells.

Abnormal Tissues

Insulin producing islet cell tumors, hyperplastic islet cells and islet cells originating in pancreatic ductules are positively labelled.² Mesenchymal tumors are not reactive.

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Literatur

1. Sternberger LA. In: Immunocytochemistry (2nd ed.). New York: John Wiley and Sons 1979;8
2. Weidenheim KM, Hinchey WW, Campbell WG Jr. Hyperinsulinemic hypoglycemia in adults with islet-cell hyperplasia and degranulation of exocrine cells of the pancreas. Amer J Clin Pathol 1983;79(1):14-24
3. Heitz PU, Kasper M, Polak JM, Kloppel G. Pancreatic endocrine tumors: Immunocytochemical analysis of 125 tumors. Hum Pathol 1982;13(3):263-71
4. Creutzfeldt W. Endocrine tumors of the pancreas. In: Bruno WE, et al. (eds). Diabetic Pancreas. New York: Plenum Press 1977:551
5. Creutzfeldt W, et al. Insulinomas and gastrinomas. In: Bloom SR and Grossmann MI (eds). Gut Hormones. Edinburgh: Churchill Livingstone 1978;96:589
6. Creutzfeldt W. Endocrine tumors of the pancreas. Clinical, chemical and morphological findings. In: Fitzgerald PJ and Morrison AB (eds.). The Pancreas. Baltimore: Williams & Wilkins 1980:208
7. Larsson L. Endocrine pancreatic tumors. Hum Pathol 1978;9(4):401-16
8. Friesen SR. Tumors of the endocrine pancreas. N Eng J Med 1982;306(10):580-90

REF Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten
 Manufacturer Fabricant Hersteller	LOT Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis
EC REP Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	IVD In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA

EC REP

Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark

Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

Tel +45 4485 9500
Fax +45 4485 9595

www.dako.com

P10039/Rev C

Edition 05/07