**Polyclonal Rabbit  
Anti-Human  
CD117, c-kit**

---

**English  
Code A4502**

**Intended use**

For In Vitro Diagnostic Use.

Polyclonal Rabbit Anti-Human CD117, c-kit is intended for use in immunocytochemistry. The antibody labels the transmembrane tyrosine kinase receptor CD117/c-kit, located in hematopoietic stem cells, melanocytes, mast cells, Cajal cells, germ cells, basal cells of skin, and mammary ductal epithelia.<sup>1-5</sup> The antibody may be a useful tool for the identification of several cancers expressing c-kit, including mast cell diseases, acute myeloid leukaemia (AML), small cell lung carcinoma (SCLC), and Ewing sarcoma, and it may aid in the differentiation between gastrointestinal stromal tumors (GISTs) and other intraabdominal mesenchymal tumors.<sup>1-4</sup> Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

**Synonyms**

KIT; stem cell factor receptor; mast cell growth factor receptor; steel factor receptor; p145(c-kit)

**Summary and explanation**

The protooncogene *c-kit*, localized to human chromosome 4,<sup>6</sup> encodes a transmembrane receptor, CD117/c-kit, belonging to the class III receptor tyrosine kinase family, which includes the receptor for colony-stimulating factor 1, and the platelet-derived growth factor receptors type A and B. The extracellular region of these receptors consists of five immunoglobulin-like domains where the second and third loops are thought to be involved in ligand binding. The intracytoplasmic tyrosine kinase domain is split by a long hydrophilic insert between the ATP-binding region and the phosphotransferase active site. Receptor activation is accompanied by receptor dimerization, substrate phosphorylation and autophosphorylation, receptor internalization, activation of protein kinases and phospholipases, and transcription of different protooncogenes.<sup>5</sup> Generally, the c-kit tyrosine kinase receptor pathway has been shown to be important for tumour growth and progression in several cancers<sup>4</sup> and mutations in the *c-kit* gene leading to ligand-independent phosphorylation (activation) of the c-kit tyrosine kinase, and are believed to have a central pathogenetic role in e.g. gastrointestinal stromal tumours.<sup>7</sup>

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

**Reagent provided**

Affinity-isolated rabbit antibody purified through antigen-bound, activated thiol AvidGel F affinity chromatography and provided in liquid form, in 0.05 mol/L Tris/HCl, 0.1 mol/L NaCl, 0.015 mol/L NaN<sub>3</sub>, pH 7.2, 1% bovine serum albumin. Package size is 1 mL.

Protein concentration g/L: See label on vial.

**Immunogen**

Peptide corresponding to amino acids 963 to 976 at the cytoplasmic C-terminal part of c-kit<sup>1,5</sup>

**Specificity**

In Western blotting of an extract of the small cell lung carcinoma cell line SY, that over-expresses c-kit mRNA, the antibody labels a band of 145 kD corresponding to the c-kit protein. The labelled band is rather broad, from 120 to 155 kD. An additional band of 100 kD is also labelled.<sup>1</sup> In another study applying a different antibody to c-kit, a 100 kD protein was, likewise, labelled. This labelling was abolished when the antibody was incubated with the synthetic c-kit peptide antigen used for immunization.<sup>6</sup>

In Western blotting, the antibody did not react with an extract of the adenocarcinoma cell line HS, that is without *c-kit* gene expression.<sup>1</sup>

**Materials required, but not supplied**

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions.

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.<sup>9,10</sup>
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused solution should be disposed of according to local, State and Federal regulations.

## Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

## Specimen preparation

### *Paraffin Sections*

The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is recommended. For heat-induced epitope retrieval, the following solutions were found efficient: Target Retrieval Solution (code S1700) and Target Retrieval Solution, pH 9 (code S2368). The following solution was found inefficient: Target Retrieval Solution, High pH (code S3308). The staining intensity was found to be reduced by a pretreatment using 0.01 mol/L citrate buffer, pH 6.0. Pre-treatment of tissues with Proteinase K destroyed the epitope. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.

## Staining procedure

### *Dilution*

Polyclonal Rabbit Anti-Human CD117, c-kit (code A4502) may be used at a dilution range of 1:400 to 1:600 when applied on formalin-fixed, paraffin-embedded sections of gastrointestinal stromal tumor and using 20 minutes heat-induced epitope retrieval in Target Retrieval Solution (code S1700) and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed, code X0936), diluted to a protein concentration identical to the concentration of the primary antibody.

### *Visualization*

LSAB+ HRP Kit (code K0679) and EnVision+ HRP Kits (codes K4008 and K4010) are recommended. Follow the procedure enclosed with the selected visualization kit.

### *Automation*

The antibody is well-suited for immunocytochemical staining using automated platforms, such as the Dako Autostainer.

## Staining interpretation

The cellular staining pattern for anti-CD117, c-kit is membranous and/or cytoplasmic.

## Performance characteristics

### *Normal Tissues*

The antibody labels breast epithelium, skin basal cells (seemingly melanocytes), spermatocytes, oocytes, tissue mast cells,<sup>1</sup> fibroblast-like bone lining cells<sup>8</sup> and intestinal Cajal cells.

### *Abnormal Tissues*

The antibody labelled 50/56 seminomas/dysgerminomas, 45/123 small cell carcinomas of the lung, 3/7 immature teratomas, 2/11 serous papillary ovarian adenocarcinomas, 2/14 malignant melanomas of skin, 2/18 neuroblastomas, 7/86 pulmonary squamous cell carcinomas, 1/18 basal cell carcinomas, 9/227 non-small cell cervical carcinomas, 1/35 pulmonary adenocarcinomas and 1/92 breast carcinomas.<sup>1</sup> The antibody was also tested on 365 soft tissue sarcomas. Tumors showing occasional immunoreactivity, in most cases focal labelling, were 3/5 melanotic schwannomas, 5/20 angiosarcomas, 4/20 metastatic melanomas, 4/20 Ewing sarcomas/malignant primitive peripheral neuroectodermal tumors, 2/10 perineuriomas, 2/20 extraskelatal myxoid chondrosarcomas. Rare tumor cells were labelled by the antibody in 1/10 low-grade fibromyxoid sarcomas and 1/20 desmoid fibromatosis. More than 35 consecutive GISTs were found to overexpress c-kit when examined with the antibody, results that were confirmed by immunoblotting. The labelling was typically diffuse and generally strong in these tumors.<sup>3</sup> The antibody also labelled activated bone marrow stem cells in bone marrow biopsies from cases of infiltrative malignant lymphomas and osteomyelitis. Labelling of activated bone marrow stem cells were more heterogeneous in cases of aplastic anaemia, myelodysplasia and metastasis from prostate or breast cancer.<sup>8</sup>

No labelling was observed in 17 large cell pulmonary carcinomas, 5 adenosquamous pulmonary carcinomas, 1 mucoepidermoid pulmonary carcinoma, 2 small cell type cervical carcinomas, 9 ovarian clear cell carcinomas, 4 ovarian mucinous cystadenocarcinomas, 4 endometrial carcinomas, 17 germ cell tumors excluding seminomas/dysgerminomas and immature teratomas, 355 gastrointestinal tract carcinomas, 59 liver carcinomas, 31 pancreas carcinomas, 7 gall bladder carcinomas, 16 bile duct carcinomas, 29 kidney carcinomas, 49 urinary bladder carcinomas, 159 prostate carcinomas, 14 adrenal carcinomas, 91 thyroid carcinomas, 45 skin carcinomas, 7 ganglioneuromas, 23 pheochromocytomas, 12 paragangliomas, 7 thyroid medullary carcinomas, 1 retinoblastoma, 30 carcinoids, 10 pulmonary neuroendocrine tumors, 20 gastrointestinal tract neuroendocrine tumors, 3 Ewing tumors and 8 pancreatic islet cell tumors.<sup>1</sup> In soft tissue tumors, no labelling was observed in 40 leiomyosarcomas, 25 rhabdomyosarcomas, 20 myxofibrosarcomas, 10 myxoid/round cell liposarcomas, 10 dedifferentiated liposarcomas, 20 solitary fibrous tumors, 20 synovial sarcomas, 30 dermatofibrosarcoma protuberans, 20 schwannomas, 20 malignant peripheral nerve sheath tumors, 10 clear cell sarcomas, 10 low-grade endometrial stromal sarcomas and 5 follicular dendritic cell sarcomas.<sup>3</sup>

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**Références**

**Literatur**

1. Tsuura Y, Hiraki H, Watanabe K, Igarashi S, Shimamura K, Fukuda T, et al. Preferential localization of *c-kit* product in tissue mast cells, basal cells of skin, epithelial cells of breast, small cell lung carcinoma and seminoma/dysgerminoma in human: immunohistochemical study on formalin-fixed, paraffin-embedded tissues. *Virchows Arch* 1994; 424:135-41
2. van Oosterom AT, Judson I, Verweij J, Stroobants S, di Paola ED, Dimitrijevic S, et al. Safety and efficacy of imatinib (STI571) in metastatic gastrointestinal stromal tumours: a phase I study. *Lancet* 2001; 358: 1421-3
3. Hornick JL, Fletcher CDM. Immunohistochemical staining for KIT (CD117) in soft tissue sarcomas is very limited in distribution. *Am J Clin Pathol* 2002; 117:188-93
4. Smithey BE, Pappo AS, Hill DA. C-kit expression in pediatric solid tumors: a comparative immuno-histochemical study. *Am J Surg Pathol* 2002; 26:486-92
5. Bühring H-J, Ashman LK, Gattei V, Kniep B, Larregina A, Pinto A, et al. CR2.7. Stem-cell factor receptor (p145(*c-kit*)) summary report (CD117). In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. *Leukocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 1882-8*
6. Yarden Y, Kuang W-J, Yang-Feng T, Coussens L, Munemitsu S, Dull TJ, et al. Human proto-oncogene *c-kit*: a new cell surface receptor tyrosine kinase for an unidentified ligand. *EMBO J* 1987; 6:3341-51
7. Miettinen M, Lasota J. Gastrointestinal stromal tumors – definition, clinical, histological, immunohisto-chemical, and molecular genetic features and differential diagnosis (review). *Virchows Arch* 2001; 438:1-12
8. Huss R, Moosmann S. The co-expression of CD117 (*c-kit*) and osteocalcin in activated bone marrow stem cells in different diseases. *Br J Haematol* 2002; 118:305-12
9. Center for Disease Control Manual Guide–Safety Management, No. CDC-22, Atlanta, GA. “Decontamination of laboratory sink drains to remove azide salts.” April 30, 1976
10. Department of Health, Education and Welfare, National Institute for Occupational Safety and Health, Rockville, MD. “Procedures for the decontamination of plumbing systems containing copper and/or lead azides.” DHHS (NIOSH) Publ. No. 78-127, Current 13. August 16, 1976

 Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i>
 Manufacturer Fabricant Hersteller	 Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis
 Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	 In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	

PT0039/Rev C

Dako North America, Inc.  
6392 Via Real  
Carpinteria, California 93013 USA

Tel 805 566 6655  
Fax 805 566 6688  
Technical Support 800 424 0021  
Customer Service 800 235 5763

**EC REP**

Dako Denmark A/S  
Produktionsvej 42  
DK-2600 Glostrup Denmark

Tel +45 4485 9500  
Fax +45 4485 9595

www.dako.com