**Polyclonal Rabbit
Anti-Human
Gastrin
Code No./ Code/ Code-Nr. A 0568
Edition/ Edition/ Ausgabe 01.12.02**

ENGLISH	
Intended use	For in vitro diagnostic use. Polyclonal Rabbit Anti-Human Gastrin, code No. A 0568, is intended for use in immunocytochemistry (1) and gel immunoprecipitation techniques. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
Reagent provided	Purified immunoglobulin fraction of rabbit antiserum provided in liquid form. In 0.1 mol/L NaCl, 15 mmol/L NaN ₃ . <u>Protein concentration g/L:</u> See label on vial.
Immunogen	Synthetic human gastrin-17 non-sulphated form, coupled to bovine serum albumin.
Specificity	The antibody reacts with human non-sulphated and sulphated forms of gastrin-17 as well as gastrin-34, and cross-reacts about 20% with cholecystokinin octapeptide (CCK-8). Traces of contaminating antibodies have been removed by solid-phase absorption with bovine plasma proteins. The specificity of the antibody has been ascertained as follows: <u>Crossed immunoelectrophoresis:</u> No precipitation arch appears when using 12.5 µL A 0568 per cm ² gel area against 2 µL of human plasma or 2 µL bovine serum. Staining: Coomassie Brilliant Blue. As demonstrated by immunocytochemistry, the antibody cross-reacts with the gastrin-equivalent protein in cow and mouse.
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. The product may be used in different techniques and in combination with different sample types and materials, therefore each individual laboratory should validate the test system applied.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, relevant controls should be run simultaneously with patient specimens. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pretreatment of tissues with heat-induced epitope retrieval is recommended. Optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.
Staining procedure	<u>Dilution:</u> Polyclonal Rabbit Anti-Human Gastrin, code No. A 0568, may be used at a dilution range of 1:400-1:800 when applied on formalin-fixed, paraffin-embedded sections of stomach and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), code No. X 0936, diluted to the same protein concentration as the primary antibody. Unless the stability in the actual test system has been established, it is recommended to dilute the product immediately before use or dilute in DakoCytomation Antibody Diluent, code No. S 0809. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4008 and K 4010, are recommended. Follow the procedure enclosed with the selected visualization kit.

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