**Monoclonal Mouse
Anti-Human
Glial Fibrillary Acidic Protein**
Clone 6F2
Code No./ Code/ Code-Nr. M 0761
Edition/ Edition/ Ausgabe 13.08.03

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Glial Fibrillary Acidic Protein, Clone 6F2, is intended for use in immunocytochemistry. The antibody labels glial fibrillary acidic protein (GFAP) and may be a useful tool for the identification of astrocytes and astrocytic cells under normal and pathological conditions (1). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Introduction	GFAP is a 50 kDa intracytoplasmic filamentous protein that constitutes a portion of the cytoskeleton in astro-cytes, and it has proved to be the most specific marker for cells of astrocytic origin. In the central rod domain of the molecule, GFAP shares considerable structural homology with the other intermediate filaments (2). Functionally, GFAP is thought to be important in astrocyte motility and shape by providing structural stability to astrocytic processes (3). Following injury to the human CNS caused by trauma, genetic disorders, or chemicals, astrocytes proliferate and show extensive hypertrophy of the cell body and processes, and GFAP is markedly upregulated. In contrast, with increasing astrocyte malignancy, there is a progressive loss of GFAP production. Thus, malignant astrocytomas have fewer tumour cells that stain positively and intensely for GFAP than do less malignant astrocytomas and normal brain specimens (3). Outside the CNS, sensitive detection methods may demonstrate GFAP in Schwann cells, enteric glia cells, salivary gland neoplasms, and metastasizing renal carcinomas. Additionally, GFAP has been demonstrated in epiglottic cartilage, pituicytes, immature oligodendrocytes, papillary meningiomas (3), and myoepithelial cells of the breast (4).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> 6F2. <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> See label on vial.
Immunogen	GFAP from human brain.
Specificity	As demonstrated by immunocytochemistry, the antibody labels GFAP in astrocytes and cells of astrocytic origin (1, 5, 6).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, DakoCytomation Target Retrieval Solution, pH 9.0, code No. S 2368, DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, or 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with DakoCytomation Proteinase K, code No. S 3020, was found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone fixed, frozen sections.
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human Glial Fibrillary Acidic Protein, code No. M 0761, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human brain and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, code No. S 1700 and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No.

K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.

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