**Polyclonal Rabbit
Anti-Human
Chorionic Gonadotropin
Code No./ Code/ Code-Nr. A 0231
Edition/ Ausgabe 01.12.02**

ENGLISH	
Intended use	For in vitro diagnostic use. Polyclonal Rabbit Anti-Human Chorionic Gonadotropin (hCG), code No. A 0231, is intended for use in immunocytochemistry. The antibody labels hCG-containing cells (1) and may be used for the demonstration of trophoblastic elements, e.g. in germ cell tumours (2). Differential identification is aided by the results from a panel of antibodies. The antibody is also well-suited for pregnancy tests provided their sensitivities are standardized to a level avoiding reaction with luteinizing hormone (LH). This level must be established for each individual test system. An example of using A 0231 in ELISA is given in (3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
Introduction	hCG is a 40 kDa glycoprotein composed of two subunits, α and β . The α subunit is common for hCG, LH, follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH), while the β subunit is unique for hCG although its amino acid sequence partially is identical with that of the β subunit of LH (3, 4). hCG is the most important marker of gestational trophoblastic cells, being present in syncytiotrophoblastic cells and cells of the intermediate trophoblast but absent in cytotrophoblasts. hCG expression has also been demonstrated in some non-trophoblastic tumours, e.g. in 14% of patients with hepatocellular carcinoma (4). In pregnancy, the serum hCG concentration increases at an exponential rate initially, then peaks at around 100 kIU/L at about the tenth week after conception. The concentration decreases about tenfold over the next several weeks followed by a small increase near term (5).
Reagent provided	Purified immunoglobulin fraction of rabbit antiserum provided in liquid form. In 0.1 mol/L NaCl, 15 mmol/L NaN ₃ . <u>Protein concentration g/L:</u> See label on vial. The titre variation between different lots of A 0231 as determined by single radial immunodiffusion is less than 10%. This is achieved by adjusting the titre of each individual lot to match the titre of an antibody reference preparation kept at -80 °C.
Immunogen	The isolated beta-chain of hCG.
Specificity	The antibody reacts with hCG, and like most other polyclonal antibodies to beta-hCG it also shows a slight cross-reaction (about 4%) with LH (3). Traces of contaminating antibodies to human plasma proteins have been removed by solid-phase absorption. The specificity of the antibody has been tested as follows: <u>Crossed immunoelectrophoresis:</u> Only the hCG precipitation arch appears when using 12.5 μ L A 0231 per cm ² gel area against 10 μ L of concentrated urine from pregnant women (15th-16th week of gestation). No precipitation is seen when 12.5 μ L antibody per cm ² gel area is tested against 10 μ L of normal human plasma. Staining: Coomassie Brilliant Blue. As demonstrated by rocket immunoelectrophoresis, the antibody cross-reacts with the hCG-equivalent protein in hamster, horse, mouse, musk ox, polecat, swine, and wildebeest (6).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. The product may be used in different techniques and in combination with different sample types and materials, therefore each individual laboratory should validate the test system applied.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, relevant controls should be run simultaneously with patient specimens. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	IMMUNOCYTOCHEMISTRY <u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, or 10 mmol/L citrate buffer, pH 6.0 and pre-treatment of tissues with proteinase K. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.
Staining procedure	<u>Dilution:</u> Polyclonal Rabbit Anti-Human Chorionic Gonadotropin (hCG), code No. A 0231, may be used at a dilution range of 1:300-1:600 when applied on formalin-fixed, paraffin-embedded sections of human placenta and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), code No. X 0936, diluted to the same protein concentration as the primary antibody. Unless the stability in the actual test system has been established, it is recommended to dilute the product immediately before use or dilute in DakoCytomation Antibody Diluent, code No. S 0809. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4008 and K 4010, are recommended. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.
Performance characteristics	Cells labelled by the antibody display a cytoplasmic staining pattern (1). <u>Normal tissues:</u> In first trimester placenta, the antibody labels plentiful syncytiotrophoblast cytoplasmic granules, some undergoing exocytosis (1).
(102357-001)	A 0231/EFG/KSS/01.12.02 p. 1/4

Abnormal tissues: 189 orchidectomy specimens with germ cell tumours (95 pure seminomas and 94 non-seminomas) were studied with the antibody. hCG was demonstrated in all choriocarcinoma components in the syncytiotrophoblasts and syncytiotrophoblast-like cells in 8% of seminomas and 30% of non-seminomas. Carcinoma-in-situ present in seminiferous tubules adjacent to the tumours was positive for hCG in a few cases (2).

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
- Morrish DW, Marusyk H, Siy O. Demonstration of specific secretory granules for human chorionic gonadotropin in placenta. J Histochem Cytochem 1987;35:93-101.
- Jacobsen GK, Jacobsen M. Alpha-fetoprotein (AFP) and human chorionic gonadotropin (hCG) in testicular germ cell tumours. A prospective immunohistochemical study. Acta Path Microbiol Immunol Scand Sect A 1983;91:165-76.
- Mehta HC, MacDonald DJ. A sensitive enzyme immunoassay specific for human chorionic gonadotrophin. Clin Chim Acta 1982;121:245-50.
- Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology; London: Oxford University Press; 1999. p. 25-6.
- The National Committee for Clinical Laboratory Standards (NCCLS). Choriongonadotropin testing; Nomenclature, reference preparations, assay performance, and clinical application; Approved guideline. NCCLS document I/LA10-A (1-56238-314-0) NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA; 1996.
- Hau J, Nilsson M, Skovgaard-Jensen H-J, de Souza A, Eriksen E, Wandall LT. Analysis of animal serum proteins using antisera against human analogous proteins. Scand J Lab Anim Sci 1990;17:3-7.

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