**Monoclonal Mouse  
Anti-Human  
Octamer-Binding Transcription Factor 3/4  
Clone N1NK**

**English  
Code M3649**

<b>Intended use</b>	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human Octamer-Binding Transcription Factor 3/4, Clone N1NK, is intended for use in immunohistochemistry. Antibodies to octamer-binding transcription factor 3/4 (OCT3/4) may be useful for the identification of specific subtypes of germ cell tumors including seminoma, embryonal carcinoma, and intratubular germ cell neoplasia of unclassified type (IGCNU) (1-4). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.</p>
<b>Synonyms for antigen</b>	<p>Octamer binding protein-3 (OCT3), octamer transcription factor-3 (OTF3), octamer binding protein-4 (OCT4), octamer transcription factor-4 (OTF4), POU5F1 (1-3, 5, 6).</p>
<b>Summary and explanation</b>	<p>Alternative splicing gives rise to two mRNA variants from the human OCT4 gene on chromosome 6, originally named OCT3A and OCT3B (7), thus the adapted OCT3/4 name. The OCT3A and OCT3B isoforms are composed of 360 and 265 amino acids, respectively, of which the 225 amino acids at the C-termini are identical. OCT3/4 belongs to the POU (Pit-Oct-Unc) transcription factor family which can activate the expression of their target genes through binding an octameric sequence motif of an AGTCAAAT consensus sequence (6). OCT3/4 is expressed in early embryonic cells and germ cells (1) and is central to the gene regulatory network responsible for self-renewal, pluripotency, and lineage commitment in embryonic stem cells and induced pluripotent stem cells (5). OCT3/4 expression is also detected in specific subtypes of germ cell tumors, such as seminomas, germinomas, dysgerminoma, embryonal carcinoma, the germ cell component of gonadoblastomas, and IGCNU (1, 2, 4).</p> <p>Refer to Dako's <a href="#">General Instructions for Immunohistochemical Staining</a> or the detection system instructions of IHC procedures.</p>
<b>Reagent provided</b>	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant in 0.05 mol/L Tris-HCl, 0.015 mol/L sodium azide, pH 7.2. This product contains stabilizing protein.</p> <p><u>Clone:</u> N1NK. <u>Isotype:</u> IgG1.</p> <p><u>Mouse IgG concentration mg/L:</u> See label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
<b>Immunogen</b>	<p>Recombinant protein corresponding to 147 amino acids at N-terminus of human OCT3/4.</p>
<b>Specificity</b>	<p>In Western blotting of hES-TW1 cell lysates, Anti-OCT3/4, clone N1NK recognizes a major band of ~40 kDa corresponding to the expected molecular weight of OCT3/4 (8).</p>
<b>Precautions</b>	<ol style="list-style-type: none"><li>1. For professional users.</li><li>2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li><li>3. As with any product derived from biological sources, proper handling procedures should be used.</li><li>4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.</li><li>5. Unused solution should be disposed of according to local, State and Federal regulations.</li></ol>
<b>Storage</b>	<p>Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.</p>

**Quick Guide\***

Step		Comments
Fixation	Formalin	N/A
Pre-treatment	EnVision FLEX™, High pH (Code K8004)	20 min HIER, 3-in-1 using PT Link and PT Link Rinse Station
Dilution	1:50	20 min incubation
Dilution Buffer	Dako Antibody Diluent (Code S0809)	Dilute immediately prior to use
Negative Control	Negative Control, Mouse IgG1 (Code X0931)	20 min incubation
Visualization	EnVision™ FLEX+, High pH (Code K8002/K8012)	20 min incubation, 15 min Mouse (LINKER) incubation, 2x5 min DAB+ incubation
Counterstain	EnVision™ FLEX Hematoxylin (Code K8008/K8018)	5 min incubation
Control Tissue	Seminoma	Nuclear staining. Cytoplasmic staining in the presence of strong nuclear staining is acceptable.
Slides	FLEX IHC Microscope Slides (Code K8020)	Recommended for greater adherence of tissue sections to glass slides.
Mounting	Non-aqueous, permanent mounting required	After staining, the sections must be dehydrated, cleared and mounted using permanent mounting medium.
Instrumentation	Autostainer Link 48 and Autostainer Plus	Use instrument-specific vials (Code SK200-SK203 and Code S3425)

*\*The user must always read the package insert for detailed instructions of the staining procedure and handling of the product.*

**Specimen preparation**

**Paraffin sections:** The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Tissue specimens should be cut into sections of approximately 4 µm.

**Pre-treatment:** Pre-treatment of formalin-fixed, paraffin-embedded tissue sections with heat-induced epitope retrieval (HIER) is required. Optimal results are obtained by pretreating tissues with HIER using diluted EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8004). Deparaffinization, rehydration and epitope retrieval can be performed in Dako PT Link (Code PT100/PT101). For details, please refer to PT Link User Guide. The following parameters should be used for PT Link: Pre-heat temperature: 65 °C; epitope retrieval temperature and time: 97 °C for 20 (±1) minutes; cool down to 65 °C. Remove slide rack from PT tank and immediately dip slides in jar/tank (e.g., PT Link Rinse Station (Code PT109)) containing diluted room temperature EnVision™ FLEX Wash Buffer (20x) (Code K8007). Leave slides in Wash Buffer for 1-5 minutes.

The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended. After staining, the sections must be dehydrated, cleared and mounted using permanent mounting medium.

**Staining procedure**

**Dilution:** The recommended dilution of Monoclonal Mouse Anti-Human Octamer-Binding Transcription Factor 3/4, Clone N1NK, Code M3649, is 1:50. Dilute the antibody in Dako Antibody Diluent (Code S0809). Incubate pretreated tissue sections for 20 minutes at room temperature. These are guidelines only. Optimal conditions may vary depending on specimen and preparation method, and should be validated individually by each laboratory.

**Negative control:** The recommended negative control reagent is Dako Negative Control, Mouse IgG1 (Code X0931), diluted to the same Ig concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, dilute these reagents immediately prior to use. Positive and negative controls should be run simultaneously with patient specimens.

**Visualization:** The recommended visualization system is EnVision™ FLEX+, High pH (Code K8002/K8012) in combination with EnVision™ FLEX+ Mouse (LINKER) (Code K8021/K8022) using a 20 minute incubation at room temperature. The incubation time for EnVision™ FLEX+ Mouse (LINKER) (Code K8021/K8022) is 15 minutes at room temperature. Follow the procedure enclosed with the selected visualization system(s).

**Automation:** The antibody is well-suited for immunohistochemical staining using automated platforms, such as Dako Autostainer, Autostainer Plus and Autostainer Link as well as PT Link for pre-treatment.

**Counterstaining:** The recommended counterstain is EnVision™ FLEX Hematoxylin (Code K8008/K8018). For optimal results, non-aqueous, permanent mounting medium is recommended.

**Controls:** Positive and negative control tissues should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include seminoma or embryonal carcinoma and the cells/structures should display reaction patterns as described for this tissue in the "Performance characteristics" section.

**Staining interpretation**

The cellular staining pattern is nuclear. Diffuse cytoplasmic staining in the presence of strong nuclear staining is acceptable.

**Performance characteristics****Normal tissues:**

Anti-OCT3/4 was non-reactive in the cell nuclei in all normal tissue types tested (9).

Tissue Type (# tested)	Positive Tissue Elements	Tissue Type (# tested)	Positive Tissue Elements
Adrenal (3)	0/3	Pancreas (3)	0/3
Bone marrow (3)	0/3	Parathyroid (3)	0/3
Breast (3)	0/3	Pituitary (3)	0/3
Cerebellum (3)	0/3	Prostate (3)	0/3
Cerebrum (3)	0/3	Salivary gland (3)	0/3
Cervix (3)	0/3	Skeletal Muscle	0/3
Colon (3)	0/3	Skin (3)	0/3
Esophagus (3)	0/3	Small intestine (3)	0/3
Heart (3)	0/3	Spleen (3)	0/3
Kidney (3)	0/3	Stomach (3)	0/3
Liver (3)	0/3	Testis (3)	0/3
Lung (3)	0/3	Thymus (3)	0/3
Mesothelial cells (3)	0/3	Thyroid (3)	0/3
Nerve, peripheral (3)	0/3	Tonsil (3)	0/3
Ovary (3)	0/3	Uterus (3)	0/3

Note: weakly positive immunostaining in the cytoplasm of epithelial cells was observed in 3 different normal tissue types. The reactive elements include epithelia of the testis (2/3), stomach (1/3), and small intestine (2/3). The weak positive staining was present in a percentage of epithelial cells in the testis (10-20%), stomach (50%), and small intestine (20%).

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**References  
Bibliographie  
Literaturnachweis**

1. Cheng L. Establishing a germ cell origin for metastatic tumors using OCT4 immunohistochemistry. *Cancer* 2004;101:2006-10.
2. Jones TD, Ulbright TM, Eble JN, Cheng L. OCT4: A sensitive and specific biomarker for intratubular germ cell neoplasia of the testis. *Clin Cancer Res* 2004;10:8544-47.
3. Ulbright TM. Germ cell tumors of the gonads: a selective review emphasizing problems in differential diagnosis, newly appreciated, and controversial issues. *Modern Path* 2005;18:S61-79
4. Høi-Hansen CE, Kraggerud SM, Abeler VM, Karen J, Meyts ERD, Lothe RA. Ovarian dysgerminoma are characterised by frequent KIT mutations and abundant expression of pluripotency markers. *Molecular Cancer* 2007;6:12
5. Bhartiya D, Kasiviswanathan S, Unni SK, Pethe P, Dhabalia J, Patwardhan S, et al. Newer insights into premeiotic development of germ cells in adult human testis using Oct-4 as a stem cell marker. *J Histochem Cytochem* 2010;58:1093-1106.
6. Pan GJ, Chang ZY, Pei D. Stem cell pluripotency and transcription factor Oct4. *Cell Res* 2002;12:321-9.
7. Takeda J, Seino S, Bell GI. Human Oct3 gene family: cDNA sequences, alternative splicing, gene organization, chromosomal location, and expression at low levels in adult tissues. *Nucleic Acids Res* 1992;20:4613-20.
8. Yu J. Mouse Monoclonal Anti-Oct3/4: A test of the specificity of Anti-Oct3/4 clone N1NK by Western blotting. 2013;D15899.
9. Yu J. Report of normal tissue immunohistochemical testing using mouse anti-human Octamer-binding transcription factor (Oct3/4), Clone N1NK. 2013;D13564.

REF Catalogue number	Temperature limitation	Consult instructions for use
Manufacturer	LOT Batch code	Use by
EC REP Authorized representative in the European Community		IVD In vitro diagnostic medical device

Dako North America, Inc.  
6392 Via Real  
Carpinteria, California 93013 USA  
Tel 805 566 6655  
Fax 805 566 6688  
Technical Support 800 424 0021  
Customer Service 800 235 5763

EC REP  
Dako Denmark A/S  
Produktionsvej 42  
DK-2600 Glostrup Denmark  
Tel +45 4485 9500  
Fax +45 4485 9595  
www.dako.com

PT0039/Rev C

Edition 08/13