**Polyclonal Rabbit
Anti-Human
Prostate-Specific Antigen
Code A0562**

ENGLISH	
Intended use	For in vitro diagnostic use. Polyclonal Rabbit Anti-Human Prostate-Specific Antigen, Code A0562, is useful for the identification of human prostate-specific antigen (PSA) in gel immunoprecipitation techniques, immunohistochemistry (1, 2) and immunoblotting (2, 3). The antibody is useful for the identification of prostate cancer especially poorly differential neoplasms (4). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Synonym for antigen	Human kallikrein 3 (hK3) (5).
Summary and explanation	Prostate-Specific Antigen (PSA) is a 33 kDa protein belonging to the kallikrein family of proteases (6). It is primarily produced by the prostatic epithelium and the epithelial lining of the periurethral glands (5), and is secreted at concentrations of 0.5-2 g/L into the seminal fluid, where it liquefies the gel forming after ejaculation by digesting the major gel-forming proteins. Five to six isoenzymes of PSA can be separated from seminal fluid by isoelectric focusing and ion exchange chromatography (6). From the normal and hyperplastic prostate, PSA leaks into the extracellular fluid and diffuses into the blood, where the concentration is approximately one million-fold lower than in seminal fluid. In prostate cancer, the architecture and polarization of the epithelial cells are deranged, leading to active secretion into the extracellular space and into the circulation, and causing increased serum PSA concentrations (6). PSA is strongly expressed in both normal and neoplastic prostatic tissue. Although PSA can be considered prostate-specific, PSA and/or PSA gene expression has been detected at low levels in some extraprostatic tissues such as normal breast tissue (5, 6), breast tumours, breast milk, female serum, endometrium, adrenal neoplasms and renal cell carcinomas (5).
Reagent provided	Purified immunoglobulin fraction of rabbit antiserum provided in liquid form. In 0.1 mol/L NaCl, 15 mmol/L NaN ₃ . <u>Protein concentration g/L:</u> See label on vial. The titre variation between different lots of A0562 is less than 10%. This is achieved by adjusting the titre of each individual lot to match the titre of an antibody reference preparation kept at -80 °C.
Immunogen	PSA isolated from pooled human seminal fluid.
Specificity	The antibody reacts with human prostate-specific antigen. Traces of contaminating antibodies have been removed by solid-phase absorption with human plasma proteins. The specificity of the antibody has been ascertained as follows: <u>Crossed immunoelectrophoresis:</u> Only the PSA precipitation arch appears when using 12.5 µL antibody per cm ² gel area against 2 µL human seminal fluid. When using 12.5 µL antibody per cm ² gel area against 2 µL human plasma no precipitate is seen. Staining: Coomassie Brilliant Blue.
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. The product may be used in different techniques and in combination with different sample types and materials, therefore each individual laboratory should validate the test system applied. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, relevant controls should be run simultaneously with patient specimens. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Services.
Specimen preparation	IMMUNOHISTOCHEMISTRY <u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. The tissue sections should not dry out during the immunohistochemical staining procedure.
Staining procedure	<u>Dilution:</u> Polyclonal Rabbit Anti-Human Prostate-Specific Antigen, Code A0562, may be used at a dilution range of 1:150-1:300 when applied on formalin-fixed, paraffin-embedded sections of human prostate and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), Code X0936, diluted to the same protein concentration as the primary antibody. Unless the stability in the actual test system has been established, it is recommended to dilute the reagent immediately before use. <u>Visualization:</u> Dako LSAB™+/HRP kit, Code K0679, and Dako EnVision™+/HRP kits, Codes K4008 and K4010, are recommended. Follow the procedure enclosed with the selected visualization kit.
Product-specific limitations	The antibody has been reported to weakly label some collecting ducts and tubular cells in the kidney, ductal cells in the parotid gland and some pancreatic cells (7).
Staining interpretation	The cellular staining pattern is cytoplasmic.
Performance characteristics	<u>Normal tissues:</u> The antibody strongly labels prostatic epithelium. No labelling was observed in 34 other human tissues, including brain, skin, muscle, nerve, breast, lymph node, spleen, thymus, different parts of the respiratory, digestive, urinary and genital systems (7).
(104258-003)	A0562/EFG/MER/21.08.08 p. 1/4

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Abnormal tissues: All 60 studied prostate carcinomas (Gleason pattern 3-5) showed diffuse positive labelling by the antibody. The percentage of PSA-positive cells was >95% in 57/60 of the carcinomas (4). The antibody labelled 3/15 (20%) ductal breast carcinomas and 2/7 (28%) apocrine breast carcinomas showed patchy areas of immunoreactive cells (7). 310 nonprostatic tumours including 40 breast adenocarcinomas, 40 lung adenocarcinomas, 40 colonic adenocarcinomas, 40 renal cell carcinomas, 40 urothelial carcinomas, 30 glandular lesions of the urinary bladder, 40 carcinoid tumours and 40 salivary gland tumours were negative with the antibody. However, diffuse strong labelling was present along the luminal aspect of non-neoplastic salivary ducts in 27/30 (90%) of the salivary gland tumours. Also 13/40 (32%) benign seminal vesicles showed cytoplasmic labelling (4).

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
- Gillespie JW, Best CJM, Bichsel VE, Cole KA, Greenhut SF, Hewitt SM, et al. Technical advance. Evaluation of non-formalin tissue fixation for molecular profiling studies. Am J Pathol 2002;160:449-57.
- Meng FJ, Shan A, Jin L, Young CYF. Short communication: The expression of a variant prostate-specific antigen in human prostate. Cancer Epidemiol Biomarkers Prev 2002;11:305-9.
- Qian Y, Sensibar JA, Zelner DJ, Schaeffer AJ, Finlay JA, Rittenhouse HG, Lee C. Two-dimensional gel electrophoresis detects prostate-specific antigen-α1-antichymotrypsin complex in serum but not in prostatic fluid. Clin Chem 1999;43:352-9.
- Varma M, Morgan M, Jasani B, Tamboli P, Amin MB. Polyclonal anti-psa is more sensitive but less specific than monoclonal anti-psa. Implications for diagnostic prostaic pathology. Am J Clin Pathol 2002;118:202-7.
- Polascik TJ, Oesterling JE, Partin AW. Prostate specific antigen: a decade of discovery – what we have learned and where we are going [Review]. J Urol 1999;162:293-306.
- Stenman U-H, Leinonen J, Zhang W-M, Finne P. Prostate-specific antigen. Cancer Biol 1999;9:83-93.
- Alanen KA, Kuopio T, Koskinen PJ, Nevalainen TJ. Immunohistochemical labelling for prostate specific antigen in non-prostatic tissues. Path Res Pract 1996;192:233-7.

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