**Monoclonal Mouse
Anti-Human
CD105, Endoglin
Clone SN6h**

ENGLISH
Code M3527

Intended use

For In Vitro Diagnostic Use.

Synonym

GP160

Summary and explanation

Endoglin is a Type I transmembrane protein which is highly expressed on human vascular endothelial cells. Up-regulation of endoglin expression has been demonstrated in tumor vasculature and proliferating cells, suggesting that it is a proliferation-associated endothelial cell marker.^{2,3} Endoglin exists as an O- and N-glycosylated homodimer of 95 kD subunits and was designated CD105 at the Fifth International Workshop on Leucocyte Differentiation Antigens (Boston).⁴

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Anti-endoglin, SN6h, is a mouse monoclonal antibody supplied in liquid form as tissue culture supernatant (containing fetal bovine serum) dialyzed against 0.05 mol/L Tris-HCl, pH 7.2 and 0.015 mol/L sodium azide. Contains stabilizing protein.

Clone: SN6h¹ Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

Immunogen

Purified CP160 from cell membrane glycoproteins of fresh non-T/non-B acute lymphoblastic leukemia cells¹

Specificity

Endoglin has been demonstrated to be a component of the transforming growth factor (TGF) β receptor system in human umbilical vein endothelial cells and binds TGF β 1 and β 3 with high affinity, but not TGF β 2.^{4,5} Monoclonal mouse anti-endoglin, SN6h, was originally described as reacting with GP160, a human leukemia-associated cell surface glycoprotein,^{1,6} which was later identified as endoglin.⁷

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions.

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Staining procedure

Follow the recommended procedure for the detection system selected.

Specimen preparation

Paraffin Sections

Anti-endoglin, SN6h, can be used on formalin-fixed, paraffin-embedded tissue sections. In order for this antibody to perform optimally on paraffin-embedded tissues, a high sensitivity detection system is required, such as CSA (code K1500) or LSAB+ HRP (code K0679). Monoclonal anti-endoglin, SN6h, may be used at a dilution of 1:2000 in the CSA detection system, determined on formalin-fixed, paraffin-embedded tissue. When using the CSA detection system, pretreatment of tissue with proteolytic enzymes is not recommended.

Anti-endoglin, SN6h, may also be used at a dilution of 1:5 to 1:10 with other high sensitivity detection systems such as the LSAB+ HRP method. Pre-treatment of formalin-fixed, paraffin-embedded tissue sections with proteolytic enzymes should be performed prior to staining. When using anti-endoglin with this detection system, one may observe a relative decrease in sensitivity compared to results obtained with the CSA system. These recommended dilutions are guidelines only, optimal dilutions should be determined by the individual laboratory.

Cryostat Sections and Cell Smears

Anti-endoglin can also be used to label cryostat sections or cell smears with the LSAB method.

Staining interpretation

The cellular staining pattern for anti-endoglin is cytoplasmic.

Performance characteristics

Normal Cells

Endoglin has been found on endothelial cells of capillaries, arterioles, small arteries, venules and high endothelial venules in a variety of tissues. Positive tissues include liver, kidney, skin, urethra, thyroid, lymph node, tonsil, umbilical cord, ovarian tube, thymus, spleen, and lung. In tonsil, lymph node and thymus, all endothelial cells, including those of high endothelial venules, express endoglin.^{4,8,9} Endoglin expression has also been demonstrated on immature proerythroblasts, syncytiotrophoblasts, activated macrophages and stromal fibroblasts from adult and fetal lung and skin.^{4,9,10}

Tumor Cells

Endoglin expression has been demonstrated on non-T/non-B and pre-B acute lymphoblastic leukemia (ALL) and acute myelocytic and myelo-monocytic leukemia cells.^{6,8} Hairy Cell leukemic cells have also been found to express endoglin.¹¹

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Références

Literatur

1. Haruta Y, et al. Multiple epitopes of GP160, a novel human leukemia-associated cell surface glycoprotein. *Fed Proc* 1987; 46(3):1056
2. Burrows FJ, et al. Up-regulation of endoglin on vascular endothelial cells in human solid tumors: Implications for diagnosis and therapy. *Clin Canc Res* 1995; 1:1623
3. Matsuzaki H, et al. Effect of induced transformation of human leukemia cells on the expression of GP160, a novel human leukemia-associated cell surface glycoprotein. *Fed Proc* 1987; 46(3):1056

4. Letarte M, et al. CD105 (endoglin) cluster report. In Schlossman SF, et al (eds). Leucocyte Typing V White Cell Differentiation Antigens. Oxford University Press (Oxford) 1995; 1756
5. Cheifetz S, et al. Endoglin is a component of the transforming growth factor-beta receptor system in human endothelial cells. J Biol Chem 1992; 267(27):19027
6. Haruta Y and Seon BK. Distinct human leukemia-associated cell surface glycoprotein GP160 defined by monoclonal antibody SN6. Proc Natl Acad Sci USA 1986; 83:7898
7. Seon BK, et al. Multiple epitopes of human endoglin, transforming growth factor-β receptor III, detected by SN6 series monoclonal antibodies. Proc Amer Assoc Canc Res 1995; 36:23
8. Gougos A and Letarte M. Identification of a human endothelial cell antigen with monoclonal antibody 44G4 produced against a pre-B leukemic cell line. J Immunol 1988; 141(6):1925
9. Buhning H, et al. Endoglin is expressed on a subpopulation of immature erythroid cells of normal human bone marrow. Leukemia 1991;5(10):841
10. Lastres P, et al. Regulated expression on human macrophages of endoglin, an Arg-Gly-Asp-containing surface antigen. Eur J Immunol 1992; 22:393
11. Haruta Y, et al. New tumor marker for Hairy Cell leukemia. Proc Amer Assoc Canc Res 1988; 29:390

REF	Catalogue number Référence du catalogue Bestellnummer		Temperature limitation Limites de température Zulässiger Temperaturbereich		Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i>
	Manufacturer Fabricant Hersteller	LOT	Batch code Code du lot Chargenbezeichnung		Use by Utiliser jusque Verwendbar bis
EC REP	Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU		IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA

Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

EC REP

Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark

Tel +45 4485 9500
Fax +45 4485 9595

www.dako.com

PT0039/Rev C