FLEX
Monoclonal Mouse
Anti-Cytokeratin 17
 Clone E3
Ready-to-Use
 (Dako Autostainer/Autostainer Plus)

Code IS620

ENGLISH	
Intended use	For in vitro diagnostic use. FLEX Monoclonal Mouse Anti-Cytokeratin 17, Clone E3, Ready-to-Use, (Dako Autostainer/Autostainer Plus), is intended for use in immunohistochemistry together with Dako Autostainer/Autostainer Plus instruments. This antibody is useful for the identification of squamous cell carcinomas (1, 2) and for the identification of pancreatic adenocarcinomas (3) and breast carcinomas (4). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Summary and explanation	Cytokeratin 17 belongs to the intermediate filaments, which create a cytoskeleton in almost all eukaryotic cells. In contrast to other intermediate filaments, cytokeratins (CKs) are made up of a highly complex multigene family of polypeptides with molecular masses ranging from 40 to 68 kDa. CKs are generally held to belong to the most fundamental markers of epithelial differentiation, and until now, 20 distinct CK polypeptides have been revealed in various human epithelia (5). The CKs can be divided into an acidic type A (class I) and a neutral-basic type B (class II) subfamily. CK 17, a 46 kDa protein, belongs to the acidic type A subfamily, and its distribution is confined to basal and myoepithelial cells of a group of complex epithelia and transitional epithelium as well as to the outer root sheet of hair folliculi (6). Refer to Dako's <i>General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.
Reagent provided	Ready-to-use monoclonal mouse antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L sodium azide. <u>Clone:</u> E3 (6). <u>Isotype:</u> IgG2b, kappa.
Immunogen	Cytoskeletal fraction of rat colon epithelium (6).
Specificity	In one- and two-dimensional SDS-PAGE and immunoblotting of Triton X-100-extracted cytoskeleton preparations of several cell cultures and tissues, the antibody labeled a band corresponding to human cytokeratin 17 (6).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation including materials required but not supplied	The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm. Pre-treatment with heat-induced epitope retrieval (HIER) is required. Optimal results are obtained by pretreating tissues using EnVision™ FLEX Target Retrieval Solution, High pH (10x), (Dako Autostainer/Autostainer Plus) (Code K8010/K8014). <u>Deparaffinized sections:</u> Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link (Code PT100/PT101). For details, please refer to the PT Link User Guide. Follow the pre-treatment procedure outlined in the package insert for EnVision™ FLEX Target Retrieval Solution, High pH (10x), (Dako Autostainer/Autostainer Plus) (Code K8010/K8014). The following parameters should be used for PT Link: Pre-heat temperature: 65 °C; epitope retrieval temperature and time: 97 °C for 20 (±1) minutes; cool down to 65 °C. Rinse sections with diluted room temperature EnVision™ FLEX Wash Buffer (10x), (Dako Autostainer/Autostainer Plus) (Code K8010). <u>Paraffin-embedded sections:</u> As alternative specimen preparation, both deparaffinization and epitope retrieval can be performed in the PT Link using a modified procedure. See the PT Link User Guide for instructions. After the staining procedure has been completed, the sections must be dehydrated, cleared and mounted using permanent mounting medium. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of Dako Silanized Slides (Code S3003) is recommended.
Staining procedure including materials required but not supplied	The recommended visualization system is EnVision™ FLEX, High pH, (Dako Autostainer/Autostainer Plus) (Code K8010). The staining steps and incubation times are pre-programmed into the software of Dako Autostainer/Autostainer Plus instruments, using the following protocols: Template protocol: FLEXRTU2 (200 µL dispense volume) or FLEXRTU3 (300 µL dispense volume) Autoprogram: CK17 (without counterstaining) or CK17H (with counterstaining) The Auxiliary step should be set to "rinse buffer" in staining runs with ≤10 slides. For staining runs with >10 slides the Auxiliary step should be set to "none". This ascertains comparable wash times. All incubation steps should be performed at room temperature. For details, please refer to the Operator's Manual for the dedicated instrument. If the protocols are not available on the used Dako Autostainer instrument, please contact Dako Technical Services. Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory. If the evaluating pathologist should desire a different staining intensity, a Dako Application Specialist/Technical Service Specialist can be contacted for information on re-programming of the protocol. Verify that the performance of the adjusted protocol is still valid by evaluating that the staining pattern is identical to the staining pattern described in "Performance characteristics".

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References/ Références/ Literatur	
<ol style="list-style-type: none">Fisseler-Eckhoff A, Erfkamp S, Müller K-M. Cytokeratin expression in preneoplastic lesions and early squamous cell carcinoma of the bronchi. Path Res Pract 1996; 92:552-559. Smedts F, Ramaekers F, Link M, Lauerova L, Troyanovsky S, Schijf C, Vooijs G P. Detection of keratin subtypes in routinely processed cervical tissue: implications for tumour classification and the study of cervix cancer aetiology. Virchows Archiv 1994;425:145-155. Chu P G, Schwarz R E, Lau S K, Yen Y, Weiss L M. Immunohistochemical staining in the diagnosis of pancreatobiliary and Ampulla of Vater adenocarcinoma. Am J Surg Pathol 2005;29:359-367. van de Rijn M, Perou C M, Tibshirani R, Haas P, Kallioniemi O, Kononen J, Torhost J, et al. Expression of cytokeratins 17 and 5 indentifies a group of breast carcinomas with poor clinical outcome. Am J Pathol 2002;161:1991-1996. Moll R. Cytokeratins as markers of differentiation in the diagnosis of epithelial tumors. In: Hermann, Harris, editors. Subcellular Biochemistry. Volume 31, New York: Plenum Press; 1998. p. 205-62. Troyanovsky SM, Guelstein VI, Tchipyшева TA, Krutovskikh VA, Bannikov GA. Patterns of expression of keratin 17 in human epithelia: dependency on cell position. J Cell Sci 1989;93:419-26. Takahashi H, Shikata N, Senzaki H, Shintaku M, Tsubura A. Immunohistochemical staining patterns of keratins in normal oesophageal epithelium and carcinoma of the oesophagus. Histopathology 1995;26:4 Cohen-Kerem R, Madah W, Sabo E, Rahat M A, Greenberg E, Elmalah I. Cytokeratin-17 as a potential marker for squamous cell carcinoma of the larynx. Ann Otol Laryngol 2004;113:821-827.	

Explanation of symbols/ Explication des symboles/ Erläuterung der Symbole			
REF	Catalogue number Référence catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Use by Utiliser avant Verwendbar bis
IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	 Contains sufficient for <n> tests Contenu suffisant pour <n> tests Inhalt ausreichend für <n> Tests	 Manufacturer Fabricant Hersteller
	Consult instructions for use Voir les instructions d'utilisation Gebrauchsanweisung beachten	LOT Batch code Numéro de lot Chargenbezeichnung	

(115290-001)			
Dako Denmark A/S Produktionsvej 42 DK-2600 Glostrup Denmark Tel. +45 44 85 95 00 Fax +45 44 85 95 95 CVR No. 33 21 13 17			

IS620/EFG/LPF/25.10.07 p. 4/4