**Monoclonal Mouse
Anti-Human
Mammaglobin
Clone 304-1A5**

Code M3625

ENGLISH

Intended use

For In Vitro Diagnostic Use

This antibody is intended for laboratory use to identify qualitatively by light microscopy mammaglobin expressing cells in normal and neoplastic tissues using immunohistochemical (IHC) test methods. Positive results aid in the classification of breast cancer. The clinical interpretation of any positive staining or its absence should be complemented by morphological and histological studies with proper controls. Evaluations should be made within the context of the patient's clinical history and other diagnostic tests by a qualified individual.

Synonyms

Antibody synonym: MGA

Antigen synonym: Mam-1A5

Summary and explanation

Mammaglobin is a 93 amino acid protein with homology to other secretoglobin-uteroglobin family members that was originally identified in studies of differentially expressed cDNAs from breast carcinoma cell lines.¹ Mammaglobin expression has also been reported in a series of malignant and nonmalignant tissues derived from the female genital tract by RT-PCR.² Expression of the mammaglobin gene and protein have been demonstrated in neoplastic breast tissues and in a subset of normal breast epithelial cells within the acini of type I and type II lobules and within the columnar cells of terminal ducts.³ Mammaglobin expression has been found in well, moderately, and poorly differentiated invasive ductal carcinomas by IHC and in metastatic and occult breast carcinoma by RT-PCR.^{3,4} Mammaglobin expression is highly restricted to breast tissue and not detected in the majority of normal tissue types evaluated, with the exception of skin.¹

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Code M3625

Monoclonal Mouse antibody provided in liquid form as tissue culture supernatant in 0.05 mol/L Tris-HCl, pH 7.2 and 0.015 mol/L sodium azide. This product contains stabilizing protein.

Clone: 304-1A5¹ Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

M3625 may be used at a dilution of 1:100 when performing IHC using the EnVision™+ detection system. These are guidelines only. Optimal antibody concentrations may vary depending on specimen and preparation method, and should be determined individually in each laboratory.

Immunogen

Purified mammaglobin recombinant protein, containing amino acids 15-93.¹

Specificity

Monoclonal Mouse anti-Mammaglobin is reported to react with a conformational epitope, as indicated by its lack of reactivity in ELISA with overlapping peptides corresponding to the mammaglobin sequence. Clone 304-1A5 reacts in Western blots with recombinant mammaglobin and with lysates of breast carcinoma cell line MDA-MB-415 and primary breast tumor cell line 450-17. No reactivity is observed in primary breast cancer cell line 450-88. Although the predicted molecular mass for mammaglobin is 9 kD, anti-Mammaglobin detected bands of 14 kD and 21 kD by Western blot. The 9 kD form of mammaglobin was observed by Western blot after enzymatic deglycosylation of MDA-MB-415 indicating glycosylation of the mammaglobin protein.¹

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions. Suggested diluent for IHC procedures:

Antibody Diluent (code S0809)

The following negative control is recommended for IHC procedures:

Mouse IgG₁ (code X0931)

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.⁵
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin Sections

Anti-Mammaglobin can be used on formalin-fixed, paraffin-embedded tissue sections. Pretreatment of tissue with proteolytic enzymes is not recommended.

The deparaffinized tissue sections must be treated with heat prior to the IHC staining procedure. Heat-induced epitope retrieval (HIER) involves immersion of tissue sections in a pre-heated buffer solution and maintaining heat in a water bath (95–99 °C). Alternative heat sources may be used for HIER upon validation against the recommended procedure. Use a 40-minute heating protocol for HIER performed at 95–99 °C; after thermal treatment, allow the jar with buffer and slides to cool for 20 minutes at room temperature. Rinse well with buffer or deionized water following HIER. For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended. Target Retrieval Solution pH 9.0 (code S2368) or 10x Concentrate (code S2367) is recommended.

Cryostat Sections And Cell Smears

Anti-Mammaglobin can be used for labeling acetone-fixed cryostat sections or fixed cell smears.

Staining procedure

Follow the procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-Mammaglobin is cytoplasmic and/or membranous.

Performance characteristics

Normal Tissues¹

Tissue Type (# tested)	Positively Staining Tissue Elements
Breast (10)	5/10 ductal epithelial cells
Colon (5)	0/5
Kidney (6)	0/6
Liver (3)	0/3
Lung (4)	0/4
Brain (3)	0/3
Heart (2)	0/2
Bone Marrow (1)	0/1
Adrenal (3)	0/3
Cervix (2)	0/2
Skin (7)	7/7 epithelial cells lining the apocrine and eccrine sweat glands
Duodenum (2)	0/2
Gall bladder (1)	0/1
Ileum (2)	0/2
Ovary (4)	0/4
Pancreas (2)	0/2
Parotid gland (1)	0/1
Prostate (4)	0/4
Spleen (3)	0/3
Testis (1)	0/1
Skeletal muscle (1)	0/1
Endometrium (3)	0/3

Abnormal Tissues^{1,6,7}

Tumor Type (# tested)	Positively Staining Tumors
Breast cancer (35)	23/35
Breast infiltrating ductal carcinoma (IDC) (96)	78/96
Endometrial carcinoma (4)	0/4
Epithelial mesothelioma (EM) (16)	1/16
Gastrointestinal carcinoma (4)	0/4
Genitourinary carcinoma (2)	0/2
Lung carcinoma (32)	0/32
Ovarian carcinoma (9)	0/9
Ovarian serous papillary carcinoma (OSPC) (21)	0/21
Pulmonary adenocarcinoma (PAC) (23)	6/23
Uterine endometrioid carcinoma (UEC) (18)	6/18

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Bibliographie

Literaturangaben

1. Fanger GR, Houghton RL, Retter MW, Hendrickson RC, Babcook J, Dillon DC, Durham MD, Reynolds LD, Johnson JC, Carter D, Fleming TP, Roche PC, Persing DH, Reed SG. Detection of mammaglobin in the sera of patients with breast cancer. *Tumour Biol* 2002;23(4):212-21
2. Grünewald K, Haun M, Fiegl M, Urbanek M, Müller-Holzner E, Massoner A, Riha K, Propst A, Marth C, Gastl G. Mammaglobin expression in gynecologic malignancies and malignant effusions detected by nested reverse transcriptase-polymerase chain reaction. *Lab Invest* 2002;82(9):1147-53
3. Watson MA, Dintzis S, Darrow CM, Voss LE, DiPersio J, Jensen R, Fleming TP. Mammaglobin expression in primary, metastatic, and occult breast cancer. *Cancer Res* 1999;59:3028-31
4. Min CJ, Tafra L, Verbanac KM. Identification of superior markers for polymerase chain reaction detection of breast cancer metastases in sentinel lymph nodes. *Cancer Res* 1998;58(20):4581-4
5. National Committee for Clinical Laboratory Standards. Protection of laboratory workers from instrument biohazards and infectious diseases transmitted by blood, body fluids, and tissue; approved guideline. Villanova, PA. 1997; Order code M29-A
6. Ciampa A, Fanger G, Khan A, Rock KL, Xu B. Mammaglobin and CRxA-01 in pleural effusion cytology: potential utility of distinguishing metastatic breast carcinomas from other cytokeratin 7-positive/cytokeratin 20-negative carcinomas. *Cancer* 2004;102(6):368-72
7. Stachurski D, Fanger G, Xu B, Edminston K, Rock KL, Khan A. CRxA-01 and mammaglobin expression in breast cancer: potential role in differentiating breast cancer from other CK7+/CK20- tumors. (205) *Mod Path* 2003;16(1):48A

REF	Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i>
	Manufacturer Fabricant Hersteller	LOT	Batch code Code du lot Chargenbezeichnung
EC REP	Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA

Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

EC REP

Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark

Tel +45 4485 9500
Fax +45 4485 9595

www.dako.com

PT0039/Rev C

Edition 07/07