**Monoclonal Mouse
Anti-Human COX-2
Clone CX-294**

**ENGLISH
Code M3617**

Intended use

For In Vitro Diagnostic Use.

This antibody is intended for laboratory use to identify qualitatively by light microscopy COX-2 expressing cells in normal and neoplastic tissues using immunohistochemical (IHC) test methods. The clinical interpretation of any positive staining or its absence should be complemented by morphological and histological studies with proper controls. Evaluations should be made within the context of the patient's clinical history and other diagnostic tests by a qualified individual.

Synonym

Cyclooxygenase-2

Summary and explanation

The cyclooxygenase (COX) enzymes, COX-1 and COX-2, are critical in the biosynthesis of prostaglandins from arachidonic acid. The human COX-1 protein is constitutively expressed in most tissues and functions in the maintenance of tissue homeostasis. Nonsteroidal anti-inflammatory drugs (NSAIDs) are effective in decreasing the inflammatory response by blocking both COX enzymes. COX-2 is a 70 kD inducible enzyme that is responsible for prostaglandin synthesis at the site of inflammation.³ COX-2 expression has also been linked to carcinogenesis, and specific COX-2 inhibitors have been shown to have antitumor effects.^{4,5}

Expression of COX-2, as determined by immunohistochemistry, has been reported in a variety of normal and neoplastic tissues. Among neoplastic tissues, heterogeneous cytoplasmic expression of COX-2 was demonstrated in colorectal adenocarcinoma, while normal adjacent tissue showed weak staining, and normal colon was negative.^{6,7} COX-2 protein was found to be expressed in a subset of breast adenocarcinomas and adjacent ductal carcinoma in situ. Using RT-PCR, COX-2 mRNA from breast tumors exhibited an increase above mRNA expression in paired normal specimens.⁸ Increased COX-2 expression was observed in lung adenocarcinomas when compared to normal bronchial epithelial cell expression.^{9,10} A majority of esophageal squamous cell carcinomas and squamous adenocarcinomas demonstrated heterogeneous cytoplasmic COX-2 expression, whereas weak and no expression was found in normal esophageal squamous epithelium and Barrett's mucosa, respectively.¹¹ COX-2 protein expression has been detected in both normal neurons and astrocytic gliomas.¹² The majority of primary malignant melanomas were found to express COX-2, while benign nevi and adjacent normal epithelium were negative.¹³ COX-2 has also been demonstrated in: gastric adenocarcinoma,¹⁴ head and neck squamous cell carcinoma,¹⁵ pancreatic carcinoma,^{16,17} malignant pheochromocytoma,¹⁸ testicular cancer,¹⁹ and squamous cell carcinoma of the urinary bladder.²⁰ COX-2 was also demonstrated in subsets of: ovarian carcinoma,²¹ prostate carcinoma,²² invasive transitional cell carcinoma of the urinary bladder,²³ renal cell carcinoma,²⁴ and hepatocellular carcinoma.²⁵

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Monoclonal Mouse antibody provided in liquid form as tissue culture supernatant in 0.05 mol/L Tris-HCl, pH 7.2 and 0.015 mol/L sodium azide. This product contains stabilizing protein.

Clone: CX-294^{1,2} Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

M3617 may be used at a dilution of 1:100 when performing IHC using the EnVision+™ detection system. These are guidelines only. Optimal antibody concentrations may vary depending on specimen and preparation method, and should be determined individually in each laboratory.

Immunogen

Synthetic peptide corresponding to human COX-2 amino acids 580-598^{1,2}

Specificity

Clone CX-294 is specific to the peptide used for immunization. In Western blotting and immunoprecipitation experiments CX-294 was shown to identify the inducible human COX-2 using IL-1 α or PMA stimulated human umbilical vein endothelial cells (HUVEC); peptide blocking of the antibody with the COX-2 immunogen eliminated all reactivity.^{1,2}

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions. Required diluent(s) for IHC procedures:

Antibody Diluent (code S0809)

The following negative control is recommended for IHC procedures:

Mouse IgG₁ (code X0931)

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin Sections

Anti-COX-2 can be used on formalin-fixed, paraffin-embedded tissue sections. Pretreatment of tissue with proteolytic enzymes is not recommended.

The deparaffinized tissue sections must be treated with heat prior to the IHC staining procedure. Heat-induced epitope retrieval (HIER) involves immersion of tissue sections in a pre-heated buffer solution and maintaining heat in a water bath or steamer (95-99 °C). Use a 20-minute heating protocol for HIER performed at 95-99 °C; after thermal treatment, allow the jar with buffer and slides to cool for 20 minutes at room temperature. Rinse well with buffer or deionized water. For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended. Target Retrieval Solution pH 9.0 (code S2367 concentrated 10x) is recommended.

Cryostat Sections And Cell Smears

Anti-COX-2 can be used for labelling acetone-fixed cryostat sections or fixed cell smears.

Staining procedure

Follow the recommended procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-COX-2 is cytoplasmic.

Performance characteristics

Normal Tissues²⁶

Tissue Type (# tested)	Positively Staining Tissue Elements
Adrenal (3)	3/3 cytoplasm of adrenal cortical cells
Bone Marrow (3)	1/3 cytoplasm of granulocyte precursors
Brain, cerebellum (3)	0/3
Brain, cerebrum (3)	3/3 faint cytoplasmic staining of neurons
Breast (3)	3/3 cytoplasm of breast epithelial cells
Cervix (3)	0/3
Colon (2)	1/2 cytoplasm of rare glandular cells
Endometrium (3)	3/3 cytoplasm of endometrial glandular epithelial cells
Esophagus (3)	3/3 cytoplasm of skeletal muscle
Heart (3)	3/3 faint cytoplasmic staining of myocytes
Kidney (3)	3/3 cytoplasm of renal tubule cells 2/3 faint staining of glomerular mesangium 1/3 faint cytoplasmic staining of glomeruli
Liver (3)	3/3 cytoplasm of hepatocytes and bile duct epithelial cells 1/3 cytoplasm of Kupffer cells
Lung (3)	3/3 membrane and cytoplasm of the intra-alveolar membrane 2/3 cytoplasm of bronchial epithelial cells
Mesothelial Cells (3)	0/3
Nerve, peripheral (3)	0/3
Ovary (3)	1/3 cytoplasm of cells of inclusion cysts
Pancreas (3)	3/3 cytoplasm of acinar cells and ducts 2/3 cytoplasm of islet cells
Parathyroid (3)	3/3 cytoplasm of chief cells
Pituitary (3)	3/3 cytoplasm of pituicytes in adenohypophysis
Prostate (3)	3/3 cytoplasm and some membrane of prostatic glandular epithelial cells
Salivary Gland (3)	3/3 cytoplasm and some membrane of ductal epithelial cells 2/3 cytoplasm of acini
Skeletal Muscle (3)	1/3 cytoplasm of perimysium
Skin (3)	3/3 cytoplasm of sebaceous gland
Small Intestine (3)	0/3
Spleen (3)	0/3
Stomach (3)	3/3 cytoplasm of chief and parietal cells
Testis (3)	0/3
Thymus (3)	3/3 cytoplasm of epithelial cells

Thyroid (3)	1/3 cytoplasm of follicular epithelial cells
Tonsil (3)	0/3

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Bibliographie

Literatur

- Creminon C, Habib A, Maclouf J, Pradelles P, Grassi J, Frobert Y. Differential measurement of constitutive (COX-1) and inducible (COX-2) cyclooxygenase expression in human umbilical vein endothelial cells using specific immunometric enzyme immunoassays. *Biochim Biophys Acta* 1995;1254(3):341-8
- Habib A, Creminon C, Frobert Y, Grassi J, Pradelles P, Maclouf J. Demonstration of an inducible cyclooxygenase in human endothelial cells using antibodies raised against the carboxyl-terminal region of the cyclooxygenase-2. *J Biol Chem* 1993;268(31):23448-54
- Hla T, Neilson K. Human cyclooxygenase-2 cDNA. *Proc Natl Acad Sci USA* 1992;89(16):7384-8
- Stratton MS, Alberts DS. Current application of selective COX-2 inhibitors in cancer prevention and treatment. *Oncology (Williston Park)*. 2002;16(5 Suppl 4):37-51
- Masferrer JL, Leahy KM, Koki AT, Zweifel BS, Settle SL, Woerner BM, Edwards DA, Flickinger AG, Moore RJ, Seibert K. Antiangiogenic and antitumor activities of cyclooxygenase-2 inhibitors. *Cancer Res* 2000;60(5):1306-11
- Sano H, Kawahito Y, Wilder RL, Hashiramoto A, Mukai S, Asai K, Kimura S, Kato H, Kondo M, Hla T. Expression of cyclooxygenase-1 and -2 in human colorectal cancer. *Cancer Res* 1995;55(17):3785-9
- Sheehan KM, Sheahan K, O'Donoghue DP, MacSweeney F, Conroy RM, Fitzgerald DJ, Murray FE. The relationship between cyclooxygenase-2 expression and colorectal cancer. *JAMA* 1999;282(13):1254-7
- Half E, Tang XM, Gwyn K, Sahin A, Wathen K, Sinicrope FA. Cyclooxygenase-2 expression in human breast cancers and adjacent ductal carcinoma in situ. *Cancer Res* 2002;62(6):1676-81
- Achiwa H, Yatabe Y, Hida T, Kuroishi T, Kozaki K, Nakamura S, Ogawa M, Sugiura T, Mitsudomi T, Takahashi T. Prognostic significance of elevated cyclooxygenase 2 expression in primary, resected lung adenocarcinomas. *Clin Cancer Res* 1999 May;5(5):1001-5
- Hosomi Y, Yokose T, Hirose Y, Nakajima R, Nagai K, Nishiwaki Y, Ochiai A. Increased cyclooxygenase 2 (COX-2) expression occurs frequently in precursor lesions of human adenocarcinoma of the lung. *Lung Cancer* 2000 Nov;30(2):73-81
- Zimmermann KC, Sarbia M, Weber AA, Borchard F, Gabbert HE, Schror K. Cyclooxygenase-2 expression in human esophageal carcinoma. *Cancer Res* 1999;59(1):198-204
- Shono T, Tofilon PJ, Bruner JM, Owolabi O, Lang FF. Cyclooxygenase-2 expression in human gliomas: prognostic significance and molecular correlations. *Cancer Res* 2001;61(11):4375-81
- Denkert C, Kobel M, Berger S, Siegert A, Leclere A, Trefzer U, Hauptmann S. Expression of cyclooxygenase 2 in human malignant melanoma. *Cancer Res* 2001;61(1):303-8
- Lim HY, Joo HJ, Choi JH, Yi JW, Yang MS, Cho DY, Kim HS, Nam DK, Lee KB, Kim HC. Increased expression of cyclooxygenase-2 protein in human gastric carcinoma. *Clin Cancer Res* 2000;6(2):519-25
- Chan G, Boyle JO, Yang EK, Zhang F, Sacks PG, Shah JP, Edelstein D, Soslow RA, Koki AT, Woerner BM, Masferrer JL, Dannenberg AJ. Cyclooxygenase-2 expression is up-regulated in squamous cell carcinoma of the head and neck. *Cancer Res* 1999;59(5):991-4
- Okami J, Yamamoto H, Fujiwara Y, Tsujie M, Kondo M, Noura S, Oshima S, Nagano H, Dono K, Umeshita K, Ishikawa O, Sakon M, Matsuura N, Nakamori S, Monden M. Overexpression of cyclooxygenase-2 in carcinoma of the pancreas. *Clin Cancer Res* 1999;5(8):2018-24
- Tucker ON, Dannenberg AJ, Yang EK, Zhang F, Teng L, Daly JM, Soslow RA, Masferrer JL, Woerner BM, Koki AT, Fahey TJ 3rd. Cyclooxygenase-2 expression is up-regulated in human pancreatic cancer. *Cancer Res* 1999;59(5):987-90
- Salmenkivi K, Haglund C, Ristimäki A, Arola J, Heikkilä P. Increased expression of cyclooxygenase-2 in malignant pheochromocytomas. *J Clin Endocrinol Metab* 2001;86(11):5615-9
- Hase T, Yoshimura R, Matsuyama M, Kawahito Y, Wada S, Tsuchida K, Sano H, Nakatani T. Cyclooxygenase-1 and -2 in human testicular tumours. *Eur J Cancer* 2003;39(14):2043-9
- Shirahama T, Sakakura C. Overexpression of cyclooxygenase-2 in squamous cell carcinoma of the urinary bladder. *Clin Cancer Res* 2001;7(3):558-61
- Denkert C, Kobel M, Pest S, Koch I, Berger S, Schwabe M, Siegert A, Reles A, Klosterhalfen B, Hauptmann S. Expression of cyclooxygenase 2 is an independent prognostic factor in human ovarian carcinoma. *Am J Pathol* 2002;160(3):893-903
- Kirschenbaum A, Klausner AP, Lee R, Unger P, Yao S, Liu XH, Levine AC. Expression of cyclooxygenase-1 and cyclooxygenase-2 in the human prostate. *Urology* 2000;56(4):671-6
- Mohammed SI, Knapp DW, Bostwick DG, Foster RS, Khan KN, Masferrer JL, Woerner BM, Snyder PW, Koki AT. Expression of cyclooxygenase-2 (COX-2) in human invasive transitional cell carcinoma (TCC) of the urinary bladder. *Cancer Res* 1999;59(22):5647-50

24. Miyata Y, Koga S, Kanda S, Nishikido M, Hayashi T, Kanetake H. Expression of cyclooxygenase-2 in renal cell carcinoma: correlation with tumor cell proliferation, apoptosis, angiogenesis, expression of matrix metalloproteinase-2, and survival. Clin Cancer Res 2003;9(5):1741-9
25. Bae SH, Jung ES, Park YM, Kim BS, Kim BK, Kim DG, Ryu WS. Expression of cyclooxygenase-2 (COX-2) in hepatocellular carcinoma and growth inhibition of hepatoma cell lines by a COX-2 inhibitor, NS-398. Clin Cancer Res 2001;7(5):1410-8
26. IHC003 Report On File

REF	Catalogue number Référence du catalogue Bestellnummer		Temperature limitation Limites de température Zulässiger Temperaturbereich		Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i>
	Manufacturer Fabricant Hersteller	LOT	Batch code Code du lot Chargenbezeichnung		Use by Utiliser jusque Verwendbar bis
EC REP	Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU			IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA

Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

EC REP

Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark

Tel +45 4485 9500
Fax +45 4485 9595

www.dako.com

PT0039/Rev C

Edition 06/07