**FLEX
Polyclonal Rabbit
Anti-Human
Thyroglobulin
Ready-to-Use
(Dako Autostainer/Autostainer Plus)**

Code IS509

ENGLISH

Intended use	For in vitro diagnostic use. FLEX Polyclonal Rabbit Anti-Human Thyroglobulin, Ready-to-Use, (Dako Autostainer/Autostainer Plus), is intended for use in immunohistochemistry together with Dako Autostainer/Autostainer Plus instruments. This antibody is useful for the identification of thyroglobulin in thyroid tissue and is a useful tool for the identification of well differentiated thyroid carcinomas (1, 2). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Summary and explanation	Thyroglobulin (Tg) is a glycoprotein with a predominant form as a 660 kDa homodimer. Tg, the precursor of thyroid hormones, is synthesized by thyrocytes and transported to the apical surface where it is secreted into the lumen of thyroid follicles and stored as the major component of colloid (> 95%). A minor proportion of Tg is found as 330 kDa monomers or as tetramers. Reduction or degradation of 660 kDa or 330 kDa Tg molecules can lead to the formation of smaller polypeptides, some of which are present in trace amounts in the colloid. At the cell-colloid interface, post-translational modifications of Tg occur, which are characterized by coupling of tyrosyl residues with iodide, leading to the formation of thyroid hormone residues within the Tg molecule. Hormone release generally requires uptake of Tg from the colloid by thyrocytes and proteolytic cleavage along the lysosomal pathway (3). Refer to Dako's <i>General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.
Reagent provided	Ready-to-use polyclonal rabbit antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L NaN ₃ .
Immunogen	Human thyroglobulin isolated from human thyroid glands.
Specificity	The antibody labels human Tg, traces of contaminating antibodies have been removed by solid-phase absorption with human plasma proteins. In crossed immunoelectrophoresis using 12.5 µL concentrated antibody per cm ² gel area against 2 µL of human thyroid gland extract, only one precipitate corresponding to Tg appears. With 2 µL of human plasma, no precipitate appears. Staining: Coomassie Brilliant Blue. In human thyroid cell cultures, the antibody labels intracellular Tg (4).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation including materials required but not supplied	The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm. Pre-treatment with heat-induced epitope retrieval (HIER) is required using Dako PT Link (Code PT100/PT101). For details, please refer to the PT Link User Guide. Optimal results are obtained by pretreating tissues using EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). <u>Paraffin-embedded sections:</u> Pre-treatment of formalin-fixed, paraffin-embedded tissue sections is recommended using the 3-in-1 specimen preparation procedure for Dako PT Link. Follow the pre-treatment procedure outlined in the package insert for EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). Note: After staining the sections must be dehydrated, cleared and mounted using permanent mounting medium. <u>Deparaffinized sections:</u> Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link and following the same procedure as described for paraffin-embedded sections. After staining the slides should be mounted using aqueous or permanent mounting medium. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended.
Staining procedure including materials required but not supplied	The recommended visualization system is EnVision™ FLEX, High pH (Dako Autostainer/Autostainer Plus) (Code K8010). The staining steps and incubation times are pre-programmed into the software of Dako Autostainer/Autostainer Plus instruments, using the following protocols: Template protocol: FLEXRTU2 (200 µL dispense volume) or FLEXRTU3 (300 µL dispense volume) Autoprogram: Thyrog (without counterstaining) or ThyrogH (with counterstaining) The Auxiliary step should be set to "rinse buffer" in staining runs with ≤10 slides. For staining runs with >10 slides the Auxiliary step should be set to "none". This ascertains comparable wash times. All incubation steps should be performed at room temperature. For details, please refer to the Operator's Manual for the dedicated instrument. If the protocols are not available on the used Dako Autostainer platform, please contact Dako Technical Services. Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory. If the evaluating pathologist should desire a different staining intensity, a Dako Application Specialist/Technical Service Specialist can be contacted

(115240-002)

IS509/EFG/MNI/2009.12.04 p. 1/4

for information on re-programming of the protocol. Verify that the performance of the adjusted protocol is still valid by evaluating that the staining pattern is identical to the staining pattern described in "Performance characteristics".

Counterstaining in hematoxylin is recommended using EnVision™ FLEX Hematoxylin (Dako Autostainer/Autostainer Plus) (K8018). Non-aqueous, permanent mounting medium is recommended.

Positive and negative controls should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include thyroid and the cells/structures should display reaction patterns as described for this tissue in "Performance characteristics" in all positive specimens. The recommended negative control reagent is FLEX Negative Control, Rabbit (Dako Autostainer/Autostainer Plus) (Code IS600).

Staining interpretation	The cellular staining pattern is cytoplasmic and lumen of thyroid follicles and the apical surface of thyrocytes.
Performance characteristics	<u>Normal tissues:</u> The antibody labels Tg in thyroid follicles and at the apical surface of thyrocytes. Thyroglobulin is found in thyroid follicular cells and to a lesser extent, in interstitial tissue and thyroid stroma due to normal leakage. <u>Abnormal tissues:</u> In hyperplastic glands, thyroglobulin has been detected in follicles and the glandular cells of smaller follicles. This antibody stains papillary, follicular and multi-hormone producing medullary carcinomas (5-7). Thyroglobulin has also been detected in follicular clear cell carcinomas whereas other clear cell neoplasms are not labeled (5, 8). Some anaplastic spindle and giant cell carcinoma cells are also reactive. Sarcomas or metastatic carcinomas are negative (5).

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