**Monoclonal Mouse
Anti-Human
Cytokeratin 10**
Clone DE-K10
Code No./ Code/ Code-Nr. M 7002
Edition/ Ausgabe 24.03.03

ENGLISH

Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Cytokeratin 10, Clone DE-K10, is intended for use in immunocytochemistry. The antibody labels all suprabasal cells in the epidermis and is a useful tool for the identification of more differentiated squamous cell carcinomas (1, 2). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Introduction	The cytokeratins (CKs) belong to the intermediate filaments, which create a cytoskeleton in almost all eukaryotic cells. In contrast to other intermediate filaments, CKs are made up of a highly complex multigene family of polypeptides with molecular masses ranging from 40 to 68 kDa. CKs are generally held to belong to the most fundamental markers of epithelial differentiation, and until now, 20 distinct CK polypeptides have been revealed in various human epithelia (3). The CKs can be divided into an acidic type A (class I) and a neutral-basic type B (class II) subfamily. CK 10 is an intermediate-sized, acidic type I cytokeratin, with a molecular mass of 56.5 kDa, expressed in all suprabasal cells simultaneously with CK 1 but is absent in basal cells. Together, CK 1 and CK 10 represent some of the first markers of epidermal differentiation (1, 4). CK 10 has been implicated as marker for more differentiated cells in squamous cell carcinomas (1, 2).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> DE-K10 (1). <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> See label on vial.
Immunogen	Cytoskeletal preparation extracted from human ectocervical epithelium (1).
Specificity	In Western blotting of cytoskeletal preparations from human epidermis, the antibody labels a single 56.5 kDa band corresponding to CK 10 (1). In two-dimensional immunoblotting of cytoskeletal proteins from human skin, the antibody labels a single dot corresponding to CK 10 (1).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is not required, but optimal results are obtained with Dako Target Retrieval Solution, High pH, code No. S 3308, Target Retrieval Solution, code No. S 1700, 10 mmol/L citrate buffer, pH 6.0 or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone-fixed, frozen sections (1).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human Cytokeratin 10, code No. M 7002, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human epidermis and using 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. <u>Visualization:</u> DAKO LSAB™+HRP kit, code No. K 0679, and DAKO EnVision™+HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.
Product-specific limitations	The antibody labels all the cells in the sebaceous glands in frozen sections, but only some cells in the basal layer of the acini on formalin-fixed, paraffin-embedded sections (1).
Performance characteristics	Cells labelled by the antibody display a cytoplasmic staining pattern (1). <u>Normal tissues:</u> The antibody labels all suprabasal layers of the epidermis, thymic Hassal's bodies and a variable amount of cells in some non-cornyfing stratified epithelia including vagina, ectocervix and tongue (1). Some cells are also labelled in hair follicles, sebaceous glands, and in the inner layer of the large secretory ducts of sweat glands (1). All simple and glandular epithelia and transitional epithelium of urinary bladder and urethra are not labelled by the antibody (1). <u>Abnormal tissues:</u> The antibody labelled 21/26 (81%) squamous cell carcinomas of the vulva, mostly in more differentiated parts (1). Various adenocarcinomas were not labelled by the antibody (1). Among 268 cervical biopsies, including 55 normal, 20 low grade squamous intraepithelial lesions (LG-SIL), 8 high grade squamous intraepithelial lesions (HG-SIL) and 29 invasive carcinomas, the median percentage of immunostaining for CK 10 was 43% in normal

tissu, 41% in LG-SIL, 17% in HG-SIL and 1% in invasive carcinomas. Using the antibody together with labelling of CK 17 (Dako Mouse Monoclonal Anti-Cytokeratin 17, code No. M 7046) a distinction of cancer specimens were facilitated (2).

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
1. Ivanyi D, Asink A, Groeneveld E, Hageman PC, Mooi WJ, Heintz APM. New monoclonal antibodies recognizing epidermal differentiation-associated keratins in formalin-fixed, paraffin-embedded tissue. Keratin 10 expression in carcinoma of the vulva. J Pathol 1989;159:7-12.
2. Maddox P, Sasieni iP, Szarewski A, Anderson M, Hanby A. Differential expression of keratins 10, 17, and 19 in normal cervical epithelium, cervical intraepithelial neoplasia, and cervical carcinoma. J Clin Pathol 1999;52:41-6.
3. Moll R. Cytokeratins as markers of differentiation in the diagnosis of epithelial tumors. In: Hermann, Harris, editors. Subcellular Biochemistry. Volume 31, New York: Plenum Press;
4. Moll R, Franke WW, Schiller DL, Geiger B, Krepler R. The catalog of human cytokeratins: patterns of expression in normal epithelia, tumors and cultured cells. Cell 1982;31:11-24.
5. Moll R. Cytokeratins as markers of differentiation in the diagnosis of epithelial tumors. Subcellular Biochem; Intermediate Filaments. Hermann, Harris, editors. Plenum Press, New York; 1998; Volume 31. p. 205-60.

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