Monoclonal Mouse Anti-Human CD35
Clone Ber-MAC-DRC
Code No./ Code/ Code-Nr. M 0846
Edition/ Edition/ Ausgabe 09.07.03

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human CD35, Clone Ber-MAC-DRC, is intended for use in immunocytochemistry. The antibody labels cells expressing CD35. The antibody may be a useful tool for characterization of histiocytic/dendritic cell neoplasms (1) and follicular dendritic cell sarcomas (2, 3). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Synonym for antigen	C3b receptor (C3bR), C4bR, Complement Receptor type 1 (CR1), Immune Adherence Receptor (4, 5).
Introduction	CD35 is a type 1 single-chain cell surface membrane glycoprotein with molecular mass polymorphism due to differences in the number of short consensus repeats (SCRs) in the extracellular region and variation in glycosylation (190 kDa, 23 SCRs; 220 kDa, 30 SCRs; 250 kDa, 37 SCRs; 280 kDa, 44 SCRs). CD35 is expressed on erythrocytes, monocytes, neutrophils, eosinophils, B lymphocytes, and 10-15% of T lymphocytes. The protein is also found on follicular glomerular podocytes, dendritic cells and some astrocytes. CD35 is a receptor for C3b and C4b bound to immune complexes and functions as a cofactor for the specific proteolytic cleavage of C3b and C4b by the plasma serine protease, factor I. CD35 also accelerates the decay of both C3 and C5 convertases. Thus, by facilitating C3b and C4b cleavage by factor I and accelerating the decay of the C3 and C5 convertases, CD35 limits complement activation and produces ligands for other complement receptors.
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L Na ₃ N. <u>Clone:</u> Ber-MAC-DRC. <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> See label on vial.
Immunogen	Suspended human tonsillar cells.
Specificity	In dot blotting, the antibody labels a recombinant form of CD35 lacking the transmembrane and cytoplasmic domains (5). In Western blotting (non-reducing) of human peripheral blood cell extract containing erythrocytes, granulocytes and mononuclear cells, the antibody labels two bands of 190 kDa and 220 kDa corresponding to 220 kDa and 250 kDa in reduced form, respectively (5). As demonstrated by immunocytochemistry, the antibody cross-reacts with the CD35-equivalent protein in cynomolgus monkey (6).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700. The following solutions, DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, and 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0, or pre-treatment of tissues with proteinase K were found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections (7, 8).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human CD35, code No. M 0846, may be used at a dilution range of 1:25 -1:50 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer (2).
Performance characteristics	Cells labelled by the antibody display a staining pattern confined to the cell membrane. <u>Normal tissues:</u> The antibody labels CD35 in renal glomeruli and in tonsils. In secondary lymphoid follicles the antibody labels both the light, dark and mantle zones (7). <u>Abnormal tissues:</u> In 3 cases of follicular dendritic cell sarcoma of the head and neck region, the antibody labelled tumour cells (3). However, a reticulum cell sarcoma with mixed follicular dendritic cell and fibroblastic reticular cell features was not labelled by the antibody (2). In 15 patients with primary Sjögren syndrome, the antibody labelled follicular dendritic cells in labial salivary glands (8).

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

(104284-002)

M 0846/EFG/KLI/09.07.03 p. 3/4

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- Pileri SA, Grogan TM, Harris NL, Banks P, Campo E, Chan JKC, et al. Tumours of histiocytes and accessory dendritic cells: an immunohistochemical approach to classification from the International Lymphoma Study Group based on 61 cases. *Histopathology* 2002;41:1-29.
- Jones D, Amin M, Ordonez NG, Glassman AB, Hayes KJ, Medeiros LJ. Reticulum cell sarcoma of lymph node with mixed dendritic and fibroblastic features. *Mod Pathol* 2001;14:1059-67.
- Biddle DA, Ro JY, Yoon GS, Yong Y-MH, Ayala AG, Ordonez NG, et al. Extranodal follicular dendritic cell sarcoma of the head and neck region: three new cases, with a review of the literature. *Mod Pathol* 2002;15:50-8.
- Nickells M, Atkinson JP. CD35. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. *Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference*; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p.780-81.
- Nickells MW, Marsh Jr. HC, Atkinson JP. CD35 Workshop Panel Report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference*; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 984-988.
- Gommerman JL, Mackay F, Donskoy E, Meier W, Martin P, Browning JL. Manipulation of lymphoid microenvironments in nonhuman primates by an inhibitor of the lymphotoxin pathway. *J Clin Invest* 2002;110:1359-69.
- Yamakawa M, Imai Y. Complement activation in the follicular light zone of human lymphoid tissues. *Immunology* 1992;76:378-84.
- Aziz KE, McCluskey PJ, Wakefield D. Characterisation of follicular dendritic cells in labial salivary glands of patients with primary Sjogren syndrome: comparison with tonsillar lymphoid follicles. *Ann Rheum Dis* 1997;56:140-3.

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