**Monoclonal Mouse
Anti-Human
Muscle Actin
Clone HHF35**

ENGLISH
Code M0635

Intended use

For In Vitro Diagnostic Use.

This antibody is intended for laboratory use to identify qualitatively by light microscopy an epitope present on muscle actin in normal and neoplastic tissues using immunohistochemical (IHC) test methods. Clone HHF35 has been demonstrated to be a reliable marker for soft tissue tumors with muscle differentiation, i.e. leiomyomas (LM), leiomyosarcomas (LMS), and rhabdomyosarcomas (RMS).² The clinical interpretation of any positive staining or its absence should be complemented by morphological and histological studies with proper controls. Evaluations should be made within the context of the patient's clinical history and other diagnostic tests by a qualified individual.

Summary and explanation

Actin, a highly conserved, ubiquitous cytoskeletal protein of muscle and nonmuscle cells, exists in three isotypes (α , β , γ) that differ by their amino acid sequences and isoelectric points. The monoclonal mouse anti-human Muscle Actin, clone HHF35 was made by immunizing mice with a polypeptide fraction of human myocardium from a case of idiopathic hypertrophic subaortic stenosis.

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Monoclonal Mouse antibody provided in liquid form as tissue culture supernatant in 0.05 mol/L Tris-HCl, pH 7.2 and 0.015 mol/L sodium azide. This product contains stabilizing protein.

Clone: HHF35 Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

M0635 may be used at a dilution of 1:50 when performing IHC using the LSAB2 detection system. These are guidelines only. Optimal antibody concentrations may vary depending on specimen and preparation method, and should be determined by each individual laboratory.

Immunogen

SDS extracted protein fraction of human myocardium

Specificity

Actin does not react with the α -actin of non-muscle (endothelial cells) sources.¹ Gel electrophoresis and immunoblots show the specificity of HHF35 to be for the α - and γ -actin isotypes of skeletal, cardiac and smooth muscle.¹

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions. Suggested diluent for IHC procedures: Dilution of this antibody in a buffer containing 0.08 mol/L EDTA is recommended to reduce nonspecific background staining.³ The following negative control is recommended for IHC procedures:
Mouse IgG₁ (code X0931).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.^{4,5}
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin Sections

Anti-Muscle Actin, HHF35 can be used on formalin-fixed, paraffin-embedded tissue sections. Pretreatment of tissue with proteolytic enzymes is not required.

Cryostat Sections and Cell Smears

Anti-Muscle Actin, HHF35 can be used for labelling acetone-fixed cryostat sections or fixed cell smears.

Staining procedure

Follow the procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-muscle actin is cytoplasmic.

Product specific limitations

1. Schmidt, et al. and others found that the addition of EDTA to an HHF35 primary antibody diluent reduced nonspecific staining^{2,3} and also decreased the chances of false-positive staining of neuroblastomas, retinoblastomas, and Ewing Sarcomas while maintaining adequate sensitivity for myogenic tumors.³
2. Miettinen found mild enzyme predigestion (pepsin, pronase, trypsin) to improve staining quality of formalin-fixed, paraffin-embedded tissue;⁷ however, Dako does not recommend tissue pretreatment.
3. Only rarely was immunoreactivity with HHF35 observed in isolated spindle cells of the liver, lymph nodes, kidney, pancreas, and the adrenal gland.⁸
4. Neoplastic cells of some pleomorphic undifferentiated sarcomas (malignant fibrous histiocytomas, MFH) have been reported positive, localized only to the smooth muscle cells and pericytes of blood vessels.^{2,7}

Performance characteristics

Normal Tissues

In normal tissue, HHF35 demonstrates cytoplasmic staining of striated fibers of skeletal muscles, the smooth muscles of arteries, veins and pericytes of smaller arteries, the tunica muscularis of the GI tract, the myometrium of the uterus, prostatic stroma, the capsule cells of several parenchymal organs, including liver, kidney, lymph nodes and spleen, and the myoepithelial layers of the mammary ducts and glands, and the eccrine sweat, bronchial and salivary glands.^{1,2,7-9} Other non-muscle cells are non-reactive, including vascular endothelial cells, epithelial cells, lymphoid cells, macrophages, connective tissue, and neural cells.^{1,2,8,9}

Abnormal Tissues

In pathological tissues, HHF35 was demonstrated to be a reliable marker for soft tissue tumors with muscle differentiation, i.e. leiomyomas (LM), leiomyosarcomas (LMS) and rhabdomyosarcomas (RMS), for which it displayed a higher degree of sensitivity than desmin antibodies.² This was confirmed by Schmidt, et al.³ who found 29/30 RMS, including embryonal, alveolar, botryoid and pleomorphic subtypes, regardless of the degree of differentiation, to be HHF35 positive. A study comprising 285 well characterized soft tissue tumors⁷ found 17/17 RMS, 31/32 LMS, 23/23 LM and 3/5 pleomorphic liposarcomas to be immunoreactive with HHF35. The majority of glomus tumors also reacted with HHF35.^{7,10} Desmoid tumors showed occasional positive cells in 9/15 cases.⁷ Similar results were reported by others⁸ who found 34/35 RMS, 11/22 LMS, 5/6 LM and 4/4 rhabdomyomas to be HHF35 reactive. The myofibroblasts of some lesions, including reactive tissue, healing wounds and atherosclerotic plaques also stained with HHF35 in the majority of cases.^{1,2,8,11} HHF35 was also used successfully for the differentiation of noninvasive (consistently actin positive) from invasive breast tumors (actin negative).⁹ Non-muscle sarcomas and neoplastic cells of carcinomas, melanomas, and lymphomas are non-reactive.^{2,7}

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Références

Literatur

1. Tsukada T, Tippens D, Gordon D, Ross R, Gown AM. HHF35, a muscle-actin specific monoclonal antibody. I. Immunocytochemical and biochemical characterization. *Amer J Pathol* 1987;126(1):51-60
2. Tsukada T, McNutt MA, Ross R, Gown AM. HHF35, A muscle-actin-specific monoclonal antibody. II. Reactivity in normal, reactive, and neoplastic human tissues. *Amer J Pathol* 1987;127(2):389-402
3. Schmidt RA, Cone R, Haas JC, Gown AM. Diagnosis of rhabdomyosarcomas with HHF35, a monoclonal antibody directed against muscle actins. *Amer J Pathol* 1988;131(1):19-28
4. Department of Health, Education and Welfare, National Institute for Occupational Safety and Health, Rockville, MD. "Procedures for the decontamination of plumbing systems containing copper and/or lead azides." 1976
5. Center for Disease Control Manual Guide - Safety Management, No. CDC-22, Atlanta, GA. "Decontamination of laboratory sink drains to remove azide salts." April 30, 1976
6. Clinical Laboratory Improvement Amendments of 1988: Final Rule, 57CFR7163. February 28, 1992
7. Miettinen M. Antibody specific to muscle actins in the diagnosis and classification of soft tissue tumors. *Amer J Pathol* 1988;130(1):205-15
8. Rangdaeng S, Truong LD. Comparative immunohistochemical staining for desmin and muscle-specific actin. *Amer J Clin Pathol* 1991;96(1):32-45
9. Gottlieb C, Raju U, Greenawald KA. Myoepithelial cells in the differential diagnosis of complex benign and malignant breast lesions: An immunohistochemical study. *Mod Pathol* 1990;3(2):135-40
10. Porter PL, Bigler SA, McNutt M, Gown AM. The immunophenotype of hemangiopericytomas and glomus tumors, with special reference to muscle protein expression: An immunohistochemical study and review of the literature. *Mod Pathol* 1991;4(1):46-52

11. Gown AM, Tsukada T, Ross R. Human atherosclerosis II. Immunocytochemical analysis of the cellular composition of human atherosclerotic lesions. Amer J Pathol 1986;125(1):191-207

 Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i>
 Manufacturer Fabricant Hersteller	 Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis
 Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU		 In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA

Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

 **EC REP**

Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark

Tel +45 4485 9500
Fax +45 4485 9595

www.dako.com

PT0039/Rev C

Edition 05/07