**Monoclonal Mouse Anti-Human CD7**  
Clone CBC.37  
**Code No./ Code/ Code-Nr. M 7255**  
Edition/ Ausgabe 19.12.02

ENGLISH	
<b>Intended use</b>	For in vitro diagnostic use.  Monoclonal Mouse Anti-Human CD7, Clone CBC.37, is intended for use in immunocytochemistry. CD7 is expressed by the majority of peripheral blood T cells, NK cells, and all thymocytes, and the antibody is a useful tool in the diagnosis of T-cell lymphoma (1). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
<b>Introduction</b>	CD7 is a 40 kDa transmembrane, single-chain glycoprotein belonging to the immunoglobulin superfamily (2). As one of the earliest surface antigens on T- and NK-cell lineages, the presence or absence of CD7 is a useful marker for the classification of T- and NK-cell malignancies. CD7 expression is found on most immature T-cell malignant cells, i.e. acute lymphocytic leukaemia (ALL) and T-cell lymphoblastic lymphoma cells. CD7 is also found in NK-cell leukaemias, including ALL syndromes that are CD16+, CD56+, and the syndrome of large granular lymphocytosis with malignant CD7+ NK cells (3).  The CD7 molecule is absent in a subset of human T cells under certain physiologic and pathologic conditions. Thus, an increase in CD7-negative T cells can be seen in some inflammatory skin diseases, and CD7 is absent in a substantial part of mature T-cell neoplasms, especially of primary cutaneous subtype (1).
<b>Reagent provided</b>	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN <sub>3</sub> . <u>Clone:</u> CBC.37 (1). <u>Isotype:</u> IgG2b, kappa. <u>Mouse IgG concentration:</u> see label on vial.
<b>Immunogen</b>	CEM T-cell line (ATCC CCL-199), a T-lymphoblastoid cell line, established from a patient with acute lymphoblastic leukaemia (1).
<b>Specificity</b>	SDS-PAGE analysis of immunoprecipitates formed between CEM cell lysate and the antibody shows that the antibody precipitates a 40 kDa polypeptide corresponding to CD7 (1).  In flow cytometry the antibody strongly labels the cell surface of Jurkat and CEM T-cell lines, whereas the Raji B-cell line is negative. Additionally the antibody blocks the reaction of the CD7 antibody, clone DK24 (1).
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
<b>Specimen preparation</b>	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin or Duboscq-Brasil (1). Pre-treatment of tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700. However, DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0 and pre-treatment of tissues with proteinase K were found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections (1).
<b>Staining procedure</b>	<u>Dilution:</u> Monoclonal Mouse Anti-Human CD7, code No. M 7255, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG2b, code No. X 0944, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.
<b>Performance characteristics</b>	Cells labelled by the antibody display staining confined to the cell membrane. <u>Normal tissues:</u> In thymus, the antibody strongly labels lymphoid cells in the medulla and cortex. In lymph node and tonsil, interfollicular T-cell areas are strongly labelled. No labelling is observed in other normal human tissues in which only small T cells are strongly labelled (1). <u>Abnormal tissues:</u> Among 110 T-cell lymphomas examined, all T-lymphoblastic lymphomas were strongly positive with the antibody (n=15), while only 25/95 of peripheral T-cell neoplasms were positive. Of note was the diagnostic value of the CD3+, CD7- T-cell phenotype in the diagnosis of lymph node involvement by cutaneous T-cell lymphoma. In Hodgkin's disease (n=15) of different categories, neoplastic cells were clearly negative and Reed Sternberg cells were often surrounded by strongly labelled small T lymphocytes. All but one of 35 B-cell

lymphomas were negative with the antibody, which labelled only small reactive T cells. Nonlymphoid tumours (n=69) were invariably negative, except for the expected labelling of small reactive T cells (1).

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