**Monoclonal Mouse
Anti-Human CD8**
Clone C8/144B
Code No./ Code/ Code-Nr. M 7103
Edition/ Ausgabe 17.12.02

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human CD8, Clone C8/144B, is intended for use in immunocytochemistry. The antibody labels cytotoxic/suppressor T cells and is a useful tool for the identification of these cells and their neoplastic counterparts (1), and for the differentiation, for example, between mycosis fungoides and cutaneous inflammatory processes (2). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Introduction	CD8 is a transmembrane glycoprotein with a molecular mass of 68 kDa. It is expressed as a disulfide-linked heterodimer comprising a 32-34 kDa α - and a 30-32 kDa β -chain, or as a homodimer comprising two α -chains. Both CD8 α and CD8 β have a typical immunoglobulin variable region-like domain in an N-terminal extracellular portion that makes them members of the immunoglobulin gene superfamily. CD8 is expressed mostly as the $\alpha\beta$ heterodimer by a majority of thymocytes, and by class I major histocompatibility complex restricted, mature, suppressor/cytotoxic T cells. A proportion of $\gamma\delta$ T cells and NK cells express the CD8 $\alpha\alpha$ homodimer (3).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> C8/144B (1). <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial.
Immunogen	Synthetic peptide corresponding to the 13 C-terminal amino acids of cytoplasmic domain of human CD8 α coupled to thyroglobulin (1).
Specificity	SDS-PAGE analysis of immunoprecipitates formed between lysates of ¹²⁵ I-labelled human T lymphoblasts and the antibody shows reaction primarily with a 32 kDa polypeptide corresponding to CD8 α (1).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin or Bouin's fixative (1). Pre-treatment of tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, or 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone-fixed, frozen sections or cell preparations (1, 2).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human CD8, code No. M 7103, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval solution, High pH, code No. S 3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.
Product-specific limitations	The antibody labels native cytokeratin 15 in human hair follicle bulge keratinocytes and a discrete area of the outer root sheet basal layer near the attachment of the arrector pili muscle without labelling the remaining hair follicle (4).
(104424-001)	M 7103/EFG/CE/17.12.02 p. 1/4
DakoCytomation Denmark A/S · Produktionsvej 42 · DK-2600 Glostrup · Denmark · Tel. +45 44 85 95 00 · Fax +45 44 85 95 95 · CVR No. 33 21 13 17	

Performance characteristics Cells labelled by the antibody display staining confined to the surface membrane.
Normal tissues: The antibody labels cytotoxic/suppressor T cells in T-cell areas of spleen and tonsil, and splenic sinusoidal lining cells (1). In colorectal biopsies, the antibody labels cytotoxic/suppressor T cells in rectum, sigmoid-, descending-, transverse-, and ascending colon, and caecum (5).
Abnormal tissues: In infectious colitis, Crohn's disease, and ulcerative colitis, the antibody demonstrated a significantly higher density of cytotoxic/suppressor T cells at affected sites of the colon than in normal controls (5). In mycosis fungoides, the CD8 antibody together with anti-CD4 demonstrated an elevated CD4:CD8 ratio in the epidermis compared with the ratio in inflammatory lesions (2).

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
- Mason DY, Cordell JL, Gaulard P, Tse AGD, Brown MH. Immunohistological detection of human cytotoxic/suppressor T cells using antibodies to a CD8 peptide sequence. J Clin Pathol 1992;45:1084-8.
- Nuckols JD, Shea CR, Horenstein MG, Burchette JL, Prieto VG. Quantitation of intraepidermal T-cell subsets in formalin-fixed, paraffin-embedded tissue helps in the diagnosis of mycosis fungoides. J Cutan Pathol 1999;26:169-75.
- Nakauchi H. TC9. CD8 Workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 65-7.
- Lyle S, Christofidou-Solomidou M, Liu Y, Elder DE, Albelda S, Cotsarelis G. Human hair follicle bulge cells are biochemically distinct and possess an epithelial stem cell phenotype. J Invest Dermatol Symp Proc 1999;4:296-301.
- Yamagata K, Tanaka M, Kudo H. A quantitative immunohistochemical evaluation of inflammatory cells at the affected and unaffected sites of inflammatory bowel disease. J Gastroenterol 1998;13:801-8.

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