**Abnormal tissues:** The antibody labelled 24/26 leiomyomas, 6/7 atypical leiomyomas and 21/25 leiomyosarcomas of the uterus, as well as 13/13 extrauterine nongastrointestinal spindle leiomyosarcomas (2). Moreover, the antibody labelled a variable amount of cells in 8/8 pseudosarcomatous myofibroblastic tumours of the urinary bladder in children (7). In pleomorphic adenomas, the antibody labelled tumour epithelial cells (myoepithelial cells) in 19/20 cases (3). In frozen tissues, the antibody, in addition to the labelling of 5/5 leiomyomas and 6/7 leiomyosarcomas, also labelled 4/22 malignant fibrous histiocytomas and 1/2 rhabdomyosarcomas. 6/6 malignant schwannomas were negative, as also 13/13 other soft tissue tumours, including 1 fibrosarcoma, 6 liposarcomas, 1 angio-sarcoma, 1 capillary haemangioma, 1 Triton tumour, and 3 synovial sarcomas (1).

## ENGLISH

**Intended use**  
For in vitro diagnostic use.  
Monoclonal Mouse Anti-Human Smooth Muscle Actin, Clone 1A4, is intended for use in immunohistochemistry. The antibody labels smooth muscle cells, myofibroblasts and myoepithelial cells, and is a useful tool for the identification of leiomyomas, leiomyosarcomas (1, 2), and pleomorphic adenomas (3). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

**Summary and explanation**  
Cyttoplasmic actins, which belong to the microfilament system of cytoskeleton proteins, are some of the most conserved eukaryotic proteins being expressed in mammals and birds. The actin protein consists of six isoforms, varying in their amino acid sequence, but all having the same molecular mass of 42 kDa. The isoforms show more than 90% overall sequence homology, but only 50-60% homology in their 18 N-terminal residues. The N-terminal region appears to be a major antigenic region (4). There are different  $\alpha$  isoforms specific for muscle tissues, i.e. skeletal muscle  $\alpha$ , cardiac muscle  $\alpha$ , and smooth muscle  $\alpha$ , respectively (1). The  $\beta$ - and  $\gamma$ -actins may be present in muscle cells as well as most other cell types in the body, including non-muscle cells (5).

**Reagent provided**  
Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>.  
**Clone:** 1A4. The 1A4 clone is identical to the anti-asm-1 described in (4). **Isotype:** IgG2a, kappa.  
**Mouse IgG concentration:** see label on vial.

The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

**Immunogen**  
N-terminal synthetic decapeptide of  $\alpha$ -smooth muscle actin coupled to keyhole limpet haemocyanin (KLH) (4).

**Specificity**  
In Western blotting and SDS-PAGE immunoblotting of the  $\alpha$ -smooth muscle isoform of actin, the antibody labels a band corresponding to  $\alpha$ -smooth muscle actin (4).  
As demonstrated by Western blotting and/or immunocytochemistry, the antibody cross-reacts with the  $\alpha$ -smooth muscle actin-equivalent protein in chicken, cow and rat (4).

**Precautions**  
1. For professional users.  
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.  
3. As with any product derived from biological sources, proper handling procedures should be used.  
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.  
5. Unused solution should be disposed of according to local, State and Federal regulations.

**Storage**  
Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Services.

**Specimen preparation**  
**Paraffin sections:** The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is recommended. Optimal results are obtained with 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with Dako Target Retrieval Solution, High pH, code S3308, or 10 mmol/L citrate buffer, pH 6.0. However, Dako Target Retrieval Solution, code S1700 was found inefficient. Pre-treatment of tissues with proteinase K was found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.  
**Frozen sections and cell preparations:** The antibody can be used for labelling acetone-fixed, frozen sections (1).

**Staining procedure**  
**Dilution:** Monoclonal Mouse Anti-Human Smooth Muscle Actin, code M0851, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of normal human colon and using 20 minutes heat-induced epitope retrieval in 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG2a, code X0943, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code S0809. Positive and negative controls should be run simultaneously with patient specimen.

**Visualization:** Dako LSAB<sup>TM</sup>+/HRP kit, code K0679, and Dako EnVision<sup>TM</sup>+/HRP kits, codes K4004 and K4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code K0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.

**Automation:** The antibody is well-suited for immunohistochemical staining using automated platforms, such as the Dako Autostainer.

**Performance characteristics**  
Cells labelled by the antibody display a cytoplasmic staining pattern.

**Normal tissues:** The antibody labels smooth muscle cells in blood vessels and, additionally, salivary ducts and myoepithelial cells around acini in salivary glands (3). Smooth muscle cells in 35/36 normal uterine myometria were also positively labelled (2). Further, a temporal labelling of perisinusoidal liver cells has been observed (6). In frozen tissues, the antibody labels myofibroblasts and myoepithelial cells around acini and ducts of the breast, whereas epithelia (adeno, squamous), lymphocytes, cardiac- and skeletal muscle cells, endothelial cells, fat cells, Schwann cells and fibroblast are negative (1).

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
- Roholl PJM, Elbers HRJ, Prinsen I, Claessens JAJ, van Unnik JAM. Distribution of actin isoforms in sarcomas: an immunohistochemical study. Hum Pathol 1990;21:1269-74.
- Rizeq MN, van de Rijn M, Hendrickson MR, Rouse RV. A comparative immunohistochemical study of uterine smooth muscle neoplasms with emphasis on the epithelioid variant. Hum Pathol 1994;25:671-7.
- Brennan PA, Umar T, Zaki GA, Langdon JD, Spedding A, Buckley J, et al. Are myoepithelial cells responsible for the widespread expression of inducible nitric oxide synthase in pleomorphic adenoma? An immunohistochemical study. J Oral Pathol Med 2000;29:279-83.
- Skalli O, Ropraz P, Trzeciak A, Benzonana G, Gillessen D, Gabbian G. A monoclonal antibody against α-smooth muscle actin: a new probe for smooth muscle differentiation. J Cell Biol 1986;103:2787-96.
- Alberts B, Bray D, Lewis J, Raff M, Roberts K, Watson JD. The cytoskeleton. In: Alberts B, Bray D, Lewis J, Raff M, Roberts K, Watson JD, editors. Molecular biology of the cell. 2nd ed. New York and London: Garland Publishing; 1989. p. 613-629.
- Schmitt-Gräff A, Krüger S, Bochard F, Gabbiani G, Denk H. Modulation of alpha smooth muscle actin and desmin expression in perisinusoidal cells of normal and diseased livers. Am J Pathol 1991;138:1233-42.
- Hojo H, Newton WA, Hamoudi AB, Qualman SJ, Wakasa H, Suzuki S, et al. Pseudosarcomatous myofibroblastic tumor of the urinary bladder in children: a study of 11 cases with review of the literature. Am J Surg Pathol 1995;19:1224-36.

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