Monoclonal Mouse Anti-Human Carcinoembryonic Antigen

Clone II-7
Code No./ Code/ Code-Nr. M 7072
 Edition/ Ausgabe 20.01.03

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Carcinoembryonic Antigen (CEA), Clone II-7, is intended for use in immunocytochemistry. The antibody labels the CEA-positive glycocalyx surface of gastrointestinal cells and is a useful tool for the identification of colon carcinomas (1). In addition, the characteristic CEA expression of secretory meningiomas, is labelled by the antibody (2). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Introduction	CEA is an oncofetal heterogeneous family of related glycoproteins (3) designated CD66 (4). The family consists of more than 17 expressible closely related genes, belonging to the immunoglobulin superfamily (4). CEA is normally present during fetal life, but occurs at low concentrations in adults, and circulates in high concentrations in patients with certain malignancies, particularly epithelial tumours (3). The protein part of CEA consists of a single polypeptide chain, containing a 107 amino acid NH-terminal domain and three highly homologous domains of 178 amino acids. CEA is secreted into the glycocalyx surface of gastrointestinal cells and hence it is especially associated with carcinomas from the gastrointestinal tract (5) – e.g. adenocarcinomas of colon, stomach and pancreas. In addition, CEA is found in adenocarcinomas of lung and breast (6).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> II-7 (3). <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial.
Immunogen	Purified human CEA from liver metastasis of lung adenocarcinoma (3).
Specificity	In ELISA, the antibody reacts with CEA preparations from both colo-rectal carcinomas, lung adenocarcinomas and the colon carcinoma cell line LS 174T, whereas the antibody does not cross-react with isolated CEA-related antigens, e.g. NCA-55, NCA-95, NCA-160 and BGP I (3). Based on ELISA, using reduced, carboxymethylated and Smith degraded CEA, the antibody is predicted to react with a conformation dependent epitope, in which the carbohydrate structures of the CEA glycoprotein are involved (3). Using solid-phase competitive immunoassays, monoclonal antibodies to CEA have been classified into five essentially non-interacting epitope groups called GOLD 1-5. Monoclonal antibodies with a high degree of CEA specificity generally belong to epitope groups GOLD 1 or 3. Monoclonal Mouse Anti-Human Carcinoembryonic Antigen (CEA), Clone II-7, belongs to epitope group GOLD 1 (7).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found less efficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections and cell preparations (1, 8).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human Carcinoembryonic Antigen (CEA), code No. M 7072, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of normal human colon or colon carcinoma, and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+ /HRP kit, code No. K 0679, and DAKO EnVision™+ /HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.
Performance characteristics	Cells labelled by the antibody display a cytoplasmic staining pattern. In tissue sections of normal colon, CEA is mainly localized at the apical border of the epithelial cells. In tissue sections of colon carcinoma, CEA is mainly localized at the apical border border of glandular structures, whereas cytoplasmic labelling predominates in more solid parts of the tumour. <u>Normal tissues:</u> The antibody labels normal colon. In frozen sections, the antibody labelled normal colon mucosa and gastric foveolae. In addition, the antibody labelled mesenchymal cells in 2/3 stomach-, 1/9 colon- and 1/3 spleen sections (1). No labelling was observed in pancreas duct- and acinar cells, liver bile canaliculi, liver Kupffer cells, spleen monocytes and spleen polymorphnuclear cells (8), the
(103474-001)	M 7072/EFG/HEW/20.01.03 p. 1/4

parenchyma and mesenchymal cells in 1 duodenum-, 2 ileum-, 3 liver-, 2 pancreas-, 1 lung-, 1 thymus- and 8 salivary gland sections nor in the parenchyma of 2 stomach-, 1 placenta-, 1 skin-, 1 tonsil- and 2 spleen sections (1).

Abnormal tissues: The antibody labelled 6/6 colon carcinomas, 2/2 gall bladder carcinomas and 2/3 carcinomas of the Vater-ampulla. In addition, the antibody labelled 31/31 secretory meningiomas (2). In frozen sections, the antibody labelled 8/9 colon carcinomas, 1/2 gastric adenocarcinomas and 1/1 small cell bronchial carcinoma (1). No labelling was observed in 1 pleomorphic adenoma of parotid, 1 bronchial squamous cell carcinoma, 1 bronchial adenocarcinoma, 2 laryngeal carcinomas, 2 bladder carcinomas, 3 breast carcinomas, 9 carcinomas of the uterine cervix, 3 malignant melanomas, 2 malignant non-Hodkin's lymphomas and 2 thymomas (1).

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