**Monoclonal Mouse  
Anti-Human CD15  
Clone Carb-3**

---

**ENGLISH**  
Code M3631**Intended use**

For In Vitro Diagnostic Use.

Monoclonal Mouse Anti-Human Clone Carb-3 is intended for laboratory use in immunohistochemistry. Antibodies to CD15 may be useful for the identification of Hodgkin's disease. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

**Synonyms for antigen**

Lewis X (Le<sup>x</sup>) antigen (1,2), 3-fucosyl-N-acetyllactosamine [3-FL] (3), X-hapten (2,4), stage-specific embryonic antigen-1 [SSEA-1](1,2), lacto-N-fucopentaose III [LNFP III] (2).

**Summary and explanation**

CD15 is expressed on Reed-Sternberg cells of Hodgkin's disease and by various other cell types including myeloid cells and epithelial cells (2,4,5). Antibodies to CD15 recognize a pentasaccharide sequence occurring in lacto-N-fucopentaose III ceramide (also referred to as X hapten of Le<sup>x</sup>) found in higher glycolipids and glycoproteins (1,4). A review by Arber et al. has reported that antibodies to CD15 demonstrate positive staining in 87% of Hodgkin's disease including nodular sclerosing, mixed cellularity, and lymphocyte depletion, whereas the lymphocyte predominant variant exhibits a lower rate of positivity (37%). Among non-Hodgkin's lymphoma, 13% express CD15 including 4.1% B-cell, 21% T-cell, and 17% null-cell. CD15 expression has also been demonstrated in acute myeloid leukemia (65%) and chronic myelogenous leukemia (96% chronic phase and 54% blast phase). A relatively low level of CD15 expression has been reported in acute lymphoblastic leukemia (5.7% overall) with positivity observed in 7.7% common or precursor B-cell, 0% B-cell, 7.7% T-cell and 17.3% null-cell. Carcinomas derived from various organs have also been shown to be CD15 positive (56%) including adenocarcinomas, squamous cell carcinomas and undifferentiated large and small cell carcinomas (4).

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

**Reagent provided**

Monoclonal Mouse antibody provided in liquid form as cell culture supernatant (containing fetal bovine serum) dialyzed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 0.015 mol/L sodium azide. This product contains stabilizing protein.

Clone: Carb-3      Isotype: IgM  
Mouse IgM concentration mg/L: See label on vial.

**Immunogen**

Immuno-affinity purified X-Hapten oligopeptide fragment from myelomonocytic leukemia cells

**Specificity**

In Western blotting of HL60 cell lysates the antibody labels a major band of ~220 kD corresponding to the expected molecular weight of CD15 (6).

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused solution should be disposed of according to local, State and Federal regulations.

**Storage**

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

### Specimen preparation including materials required but not supplied

#### Paraffin sections:

The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin.

Pre-treatment of deparaffinized tissue sections with heat-induced epitope retrieval (HIER) is required. Optimal results are obtained by pretreating tissues with HIER using Dako Target Retrieval Solution (Codes S1700/S1699).

Waterbath HIER: 20 minutes at 95–99 °C. Preheat target retrieval solution prior to immersing slides.

After thermal treatment, allow the jar with buffer and slides to cool for 20 minutes at room temperature. Rinse sections gently with buffer or deionized water. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of Dako Silanized Slides (Code S3003) is recommended.

#### Frozen sections and cell smears:

The antibody can be used for labeling acetone-fixed, frozen sections or fixed cell smears.

### Staining procedure including materials required but not supplied

**Dilution:** Monoclonal Mouse Anti-Human CD15, Code M3631, may be used at a dilution of 1:50 when applied on pretreated, formalin-fixed, paraffin-embedded sections using a 30 minute incubation at room temperature. It is recommended that the antibody is diluted in Dako Antibody Diluent (Code S0809). These are guidelines only. Optimal conditions may vary depending on specimen and preparation method, and should be determined individually by each laboratory. The recommended negative control reagent is Dako Negative Control, Mouse IgM (Code X0942), diluted to the same mouse IgM concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately prior to use. Positive and negative controls should be run simultaneously with patient specimens.

**Visualization:** Dako EnVision™+ System/HRP, Dual Link Rabbit/Mouse (Code K4061) is recommended. Follow the procedure enclosed with the selected visualization system.

**Automation:** The antibody is well-suited for immunohistochemical staining using automated platforms, such as the Dako Autostainer.

### Staining interpretation

The cellular staining pattern is cytoplasmic and/or membranous (4, 7).

### Performance characteristics

#### Normal tissues (7):

Tissue Type (#tested)	Positively Staining Tissue Elements
Adrenal (3)	3/3 Endocrine cells of adrenal cortex (100%), cytoplasmic
Bone Marrow (3)	2/3 Myelocytes (100%), cytoplasmic 1/3 All cells (100%), cytoplasmic
Breast (2)	2/2 Ductal epithelium (25-100%), cytoplasmic
Brain/Cerebellum (3)	3/3 White matter, grey matter (100%), diffuse
Brain/Cerebrum (3)	3/3 White matter, grey matter (100%), diffuse
Cervix (2)	0/2
Colon (2)	1/2 Mucosal glandular epithelium (50%), punctuate, apical cytoplasmic 1/2 Mucosal epithelium (20%), apical cytoplasmic
Esophagus (2)	2/2 Non-basal squamous epithelium (100%), cytoplasmic
Heart (3)	0/3
Kidney (3)	3/3 Cortical tubular epithelium (100%), cytoplasmic 2/3 Cortical and medullary endothelium (100%), cytoplasmic
Liver (3)	1/3 Sinusoidal endothelium (40%), cytoplasmic
Lung (3)	0/3
Mesothelial Cells (3)	2/3 Mesothelial cells (100%), cytoplasmic
Nerve (3)	0/3
Ovary (3)	0/3
Pancreas (3)	2/3 Acinar cell (25-50%), cytoplasmic
Parathyroid (3)	0/3
Pituitary (3)	3/3 Endocrine cells of adenohypophysis (75-100%), cytoplasmic 2/3 Neuropil of neurohypophysis (75-100%), diffuse
Prostate (2)	1/2 Glandular epithelium (50%), cytoplasmic
Salivary gland (3)	1/3 Ductal epithelium (100%), apical cytoplasmic 1/3 Glandular epithelium (20%), cytoplasmic
Skeletal muscle (3)	0/3
Skin (3)	0/3
Small intestine (3)	1/3 Mucosal epithelium (20%), apical cytoplasmic
Spleen (3)	3/3 Myeloid cells (100%), cytoplasmic
Stomach (3)	3/3 Gastric mucosa (100%), cytoplasmic

Testis (3)	0/3
Thymus (3)	2/3 Hassel's corpuscles (100%), non-cellular
Thyroid (2)	0/2
Tonsil (3)	2/3 Overlying squamous epithelium (50-100%), cytoplasmic
Uterus (2)	2/2 Glandular epithelium (15-75%), cytoplasmic

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
**References/ Bibliographie/ Literaturangaben**

- Ball ED. CD15 cluster workshop report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, Ritz J, Shaw S, Silverstein R, Springer T, Tedder TF, Todd RF (eds). Leucocyte Typing V. White cell differentiation antigens. Oxford: Oxford Univ Press 1993;1:790-805
- Kannagi R. CD15 Workshop Panel report. In: Kishimoto T, Kikutani H, von dem Borne AEGK, Goyert SM, Mason DY, Miyasaka M, Moretta L, Okumura K, Shaw S, Springer TA, Sugamura K, Zola H (eds). Leucocyte Typing VI. White Cell Differentiation Antigens. Garland Publishing Inc. New York, 1998;348-51
- Gocht A, Struckhoff G, Löhler J. Invited review. CD15-containing glycoconjugates in the central nervous system. *Histol Histopathol* 1996;11:1007-28
- Arber DA, Weiss LM. CD15: a review. *Applied Immunohistochem* 1993;1(1):17-30
- Kornstein MJ, Bonner H, Gee B, Cohen R, Brooks JJ. Leu M1 and S100 in Hodgkin's disease and non-Hodgkin's lymphomas. *AJCP* 1986;85(4):433-7
- M3631 WB011607 Report On File
- M3631 IHC003 Report On File

 Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i>
 Manufacturer Fabricant Hersteller	 Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis
 Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU		 In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum

PT0039/Rev C

Dako North America, Inc.  
6392 Via Real  
Carpinteria, California 93013 USA

Tel 805 566 6655  
Fax 805 566 6688  
Technical Support 800 424 0021  
Customer Service 800 235 5763

**EC REP**

Dako Denmark A/S  
Produktionsvej 42  
DK-2600 Glostrup Denmark

Tel +45 4485 9500  
Fax +45 4485 9595

www.dako.com

Edition 02/07