Monoclonal Mouse Anti-Human Thyroid Peroxidase
Clone MoAb47
Code No./ Code/ Code-Nr. M 7257
Edition/ Ausgabe 19.12.02

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Thyroid Peroxidase (TPO), Clone MoAb47, is intended for use in immunocytochemistry. The antibody labels TPO and may be a useful aid in the differentiation between benign and malignant thyroid tumours (1-3). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Introduction	Thyroid peroxidase (TPO) is a transmembrane protein of 107 kDa containing a heme prosthetic group and it is present as a dimer on the apical surface of the thyroid follicular cells (4, 5). TPO is the primary enzyme involved in thyroid hormone synthesis, catalysing iodide oxidation, iodination of tyrosine residues, and coupling of iodotyrosines to generate the iodothyronines T ₃ and T ₄ (5). The TPO gene consists of 17 exons, and is located on the short arm of chromosome 2 (2). Malignant thyroid tumours exhibit an anomaly in TPO resulting in a lower affinity for Anti-TPO, MoAb47. Decreased TPO immunoreactivity is an early event in follicular tumourigenesis, taking place before development of invasiveness in parallel with an acceleration of cell growth and appearance of cell atypia (2).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> MoAb47 (6). <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial.
Immunogen	Human TPO (hTPO) purified from thyroid microsomes by immunoaffinity chromatography (6).
Specificity	As demonstrated by ELISA and inhibition experiments, the antibody labels purified hTPO (6). The epitope recognized by the antibody was found to be within nucleotides 2224-2260 in exon 12 of the TPO gene (2). Pooled sera from patients with Grave's disease and Hashimoto's thyroiditis at a dilution of 1:20 hinder the binding of the antibody to hTPO coated onto ELISA wells, indicating that the patient sera share some epitope reactivity with the antibody (6).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, or 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections (1).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human Thyroid Peroxidase (TPO), code No. M 7257, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of human thyroid gland and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.
Performance characteristics	Normal thyroid follicular cells labelled by the antibody show a slightly granular staining of the cytoplasm. In the perinuclear area, a well-defined ring can be observed in most cells (1). <u>Normal tissues:</u> The antibody labels thyroid follicular cells. Connective tissue, nerve sheets, vessels and pancreas are negative with the antibody (2). Additionally, no labelling is observed in tonsil, liver, colon, cerebellum, kidney and prostate. <u>Abnormal tissues:</u> Frozen sections of fine needle aspiration biopsies from 150 thyroid nodules were labelled with the antibody. In 113/125 of the benign nodules 80-100% cells were positive, whereas all samples from malignant tumours yielded less than 80% positive cells. Most papillary cancers were completely negative (1). In another study (3), formalin-fixed, paraffin-embedded fine needle aspiration biopsies from 124 solitary cold thyroid nodules were labelled with the antibody. Samples were considered benign if 80% or more of the epithelial-looking
(103479-001)	M 7257/EFG/CE/19.12.02 p. 1/4
DakoCytomation Denmark A/S · Produktionsvej 42 · DK-2600 Glostrup · Denmark · Tel. +45 44 85 95 00 · Fax +45 44 85 95 95 · CVR No. 33 21 13 17	

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cells were positive, and malignant if more than 20% of these cells were negative. Tissue sections from the nodules obtained during a subsequent operation and labelled with the antibody or stained with haematoxylin-eosin, respectively, were used as controls. Also in this study, the fine needle aspiration biopsies with less than 80% positive cells were subsequently proven to be malignant. A universal and reliable pattern with more than 80% positive cells was found in all 97 subsequently proven benign lesions, with the exception of one out of 26 follicular adenomas. This gave the method a 100% sensitivity and a 99% specificity for the diagnosis of malignancy on optimally sampled and immunostained fine needle aspiration biopsies. No labelling has been observed in carcinomas of breast and colon, and in carcinoid and metastatic melanomas.

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

References/ Références/ Literatur	<ol style="list-style-type: none">De Micco C, Zoro P, Garcia S, Skoog L, Tani EM, Carayon P, et al. Thyroid peroxidase immunodetection as a tool to assist diagnosis of thyroid nodules on fine-needle aspiration biopsy. Eur J Endocrinol 1994;131:474-9. De Micco C, Kopp F, Vassko V, Grino M. In situ hybridization and immunohistochemistry study of thyroid peroxidase expression in thyroid tumors. Thyroid 2000;10:109-15. Christensen L, Blichert-Toft M, Brandt M, Lange M, Sneppen SB, Ravnsbæk J, et al. Thyroperoxidase (TPO) immunostaining of the solitary cold thyroid nodule. Clin Endocrinol 2000;53:161-9. Rapoport B, McLachlan SM. Thyroid autoimmunity. J Clin Invest 2001;108:1253-9. McLachlan SM, Rapoport B. The molecular biology of thyroid peroxidase: cloning, expression and role as autoantigen in autoimmune thyroid disease. Endocr Rev 1992;13:192-206. Ruf J, Toubert M-E, Czarnocka B, Durand-Gorde J-M, Ferrand M, Carayon P. Relationship between immunological structure and biochemical properties of human thyroid peroxidase. Endocrinology 1989;125:1211-8.	
Explanation of symbols/ Légende des symboles/ Erläuterung der Symbole		
REF <p>Catalogue number Référence du catalogue Bestellnummer</p>	 <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p>	 <p>Manufacturer Fabricant Hersteller</p>
IVD <p>In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum</p>	LOT <p>Batch code Code du Lot Chargenbezeichnung</p>	
 <p>Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten</p>	 <p>Use by Utiliser jusque Verwendbar bis</p>	

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