**Monoclonal Mouse
Anti-Human
MCM3 Protein
Clone 101
Code No./ Code/ Code-Nr. M 7263**

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human MCM3 Protein, Clone 101, is intended for use in immunocytochemistry. The antibody labels minichromosome maintenance 3 (MCM3) protein in the nucleus of proliferating cells and in the nucleus of cells that have ceased to proliferate, but are not terminally differentiated (1). Antibodies against MCM3 protein may be useful for the characterization of different tumour cell populations (1). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Introduction	The exact duplication of a genome once per cell division is required of every proliferating cell. In eukaryotic cells this is achieved by temporally separating the assembly of a pre-replication complex (pre-RC) from initiation of DNA synthesis at replication origins. A key component of the pre-RC is the hexameric MCM complex composed of the highly conserved proteins MCM2, MCM3, MCM4, MCM5, MCM6 and MCM7 (2, 3). This MCM complex is also a DNA helicase. To ensure that initiation of DNA synthesis occurs only once per cell division, the recruitment of the MCM complex is temporally separated from its activation. In G ₁ phase, the MCM complex is recruited to the pre-RC in an inactive globular conformation. In S phase, the MCM complex becomes phosphorylated, which leads to a conformational change that converts the inactive MCM complex into an enzymatically active helicase that becomes topologically linked to DNA and disassociates from the pre-RC. This again leads to the melting of dsDNA and recruitment of proteins involved in DNA synthesis (2). Like Ki-67 antigen, MCM3 protein is expressed in proliferating cells (1, 4). However, MCM3 protein disappears more slowly after initiation of differentiation than the Ki-67 antigen, and antibodies against MCM3 protein label more cells than antibodies against Ki-67 antigen (1, 4-6). This observation, correlated with the fact that MCM3 protein is not expressed when a marker of terminal differentiation, such as p27, is expressed, indicates that antibodies against MCM3 protein label cells in the early stages of cell cycle exit, i.e. cells that have not terminally differentiated (1). Antibodies against MCM3 protein in combination with a panel of other antibodies may identify dysplastic cells that remain in cycle owing to deregulation of normal controls over cell proliferation as a step in their progression from normal to neoplastic cells (6).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> 101 (1). <u>Iso</u> <u>type:</u> IgG1, kappa. <u>Mouse IgG concentration mg/L:</u> See label on vial.
Immunogen	Recombinant human MCM3 protein (1).
Specificity	In Western blotting of lysates derived from human tonsil, as well as from cell lines of human (HeLa, S3, IM-9), mouse (Swiss 3T3) or hamster (CHO-K1) origin, the antibody labels a ~100 kDa band corresponding to MCM3 protein. An additional band of ~85 kDa was observed in varying amounts. The nature of this band is not entirely clear, but it could be proteolytically degraded MCM3 protein (1). In two-colour immunofluorescence of human tonsil, the antibody labels cells that express Ki-67 antigen, but not cells expressing p27. Additionally, the antibody labels cells in the intermediate layer of the oral mucosa epithelium that are negative for Ki-67 antigen (1). As demonstrated by immunocytochemistry, the antibody cross-reacts with the MCM3-equivalent protein in mouse (1, 7), cow, dog, horse, rat, and swine.
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, and 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human MCM3 Protein, code No. M 7263, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.
Performance characteristics	Cells labelled by the antibody display a nuclear staining pattern. <u>Normal tissues:</u> In human tonsil, the antibody predominantly labels cells in proliferative regions, i.e. cells of the dark zone within germinal centres, and cells near the basal layer of the normal mucosa (1).
(109524-001)	M 7263/EFG/KLI/23.04.04 p. 1/4

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

- Endl E, Kausch I, Baack M, Knippers R, Gerdes J, Scholzen T. The expression of Ki-67, MCM3, and p27 defines distinct subsets of proliferating, resting, and differentiated cells. J Pathol 2001;195:457-62.
- Lei M, Tye BK. Initiating DNA synthesis: from recruiting to activating the MCM complex. J Cell Sci 2001;114:1447-54.
- Tye BK, Sawyer S. The hexameric eukaryotic MCM helicase: building symmetry from nonidentical parts. J Biol Chem 2000;275:34833-6.
- Musahl C, Holthoff HP, Lesch R, Knippers R. Stability of the replicative Mcm3 protein in proliferating and differentiating human cells. Exp Cell Res 1998;241:260-4.
- Ishimi Y, Okayasu I, Kato C, Kwon HJ, Kimura H, Yamada K, et al. Enhanced expression of MCM proteins in cancer cells derived from uterine cervix. Eur J Biochem 2003;270:1089-1101.
- Freeman A, Morris LS, Mills AD, Stoeber K, Laskey RA, Williams GH, et al. Minichromosome maintenance proteins as biological markers of dysplasia and malignancy. Clin Cancer Res 1999;5:2121-32.
- Birner P, Ritzl M, Musahl C, Knippers R, Gerdes J, Voigtländer T, et al. Immunohistochemical detection of cell growth fraction in formalin-fixed and paraffin-embedded murine tissue. Am J Pathol 2001;158:1991-6.

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