**Monoclonal Mouse
Anti-Human Chromogranin A**
Clone DAK-A3
Code No./ Code/ Code-Nr. M 0869
Edition/ Ausgabe 19.12.02

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Chromogranin A, Clone DAK-A3, is intended for use in immunocytochemistry. The antibody labels chromogranin A present in secretory granules of endocrine cells. Chromogranin A is a marker for neuroendocrine differentiation and antibodies to chromogranin A are a useful tool for the identification of neuroendocrine-derived tumours. Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Introduction	Chromogranin A is a protein prohormone with a molecular mass of 48 kDa. The gene encoding the 439 amino acids protein is localized on chromosome 14 (1-3). Potentially biologically active peptides derived from chromogranin A are; vasostatins (aa 1-17/113), chromostatin (aa 124-143), chromacins (aa 173-194), pancreastatin (aa 243-294), WE-14 (aa 316-329), catestatin (aa 344-364), parastatin (aa 347-419) and GE-25 (aa 367-391) (1). Chromogranin A was first isolated from chromaffin cells in the adrenal medulla but was later found to be present in secretory granules in cells of most neuroendocrine organs. Since then it has been used as an important broadspectrum histochemical marker for the identification of neuroendocrine cells and tumours (1, 4, 5).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 50 mmol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> DAK-A3. <u>Iso</u> type: IgG2b, kappa. <u>Mouse IgG concentration:</u> see label on vial.
Immunogen	The C-terminal 20 kDa fragment of human chromogranin A corresponding to amino acids 210-439. The fragment has been isolated from the urine of patients with carcinoid syndrome, and was pure as shown by N-terminal amino acid analysis.
Specificity	In ELISA, the antibody reacts with the pure chromogranin A fragment used for immunization. In immunocytochemistry the antibody shows a staining pattern essentially as that of an established polyclonal antibody to chromogranin A, DakoCytomation code No. A 0430, when tested on 84 different tumours and 30 different normal tissues. The labelling corresponds to the known distribution of chromogranin A. In addition, labelling of melanin and lipofuscin pigment was observed with both antibodies as also with an unrelated monoclonal antibody to chromogranin A, Clone LK2H10, after demasking of antigens with citrate buffer.
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, or 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections:</u> The antibody can be used for labelling acetone-fixed frozen sections.
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human Chromogranin A, code No. M 0869, may be used at a dilution range of 1:100-1:200 when applied on formalin-fixed, paraffin-embedded sections of human colon and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval solution, High pH, code No. S 3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG2b, code No. X 0944, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen.

Visualization: DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. Follow the procedure enclosed with the selected visualization kit.

Automation: The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.

Product-specific limitations	The antibody labels melanin and lipofuscin pigment after heat induced epitope retrieval in 10 mmol/L citrate buffer, pH 6.0.
Performance characteristics	Cells labelled by the antibody display a cytoplasmic, granular staining pattern. <u>Abnormal tissues:</u> The antibody labelled 13/18 of neuroendocrine tumours of various types. None of 39 carcinomas, 7 sarcomas, 4 melanomas, 6 lymphomas, 2 seminomas, 2 mesotheliomas, 1 each of meningioma, astrocytoma and Warthin's tumour were labelled.

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- Portela-Gomes GM. Chromogranin A immunoreactivity in neuroendocrine cells in the human gastrointestinal tract and pancreas. Adv Exp Med Biol 2000;482:193-203.
- Konecki DS, Benedum UM, Gerdes H-H, Huttner WB. The primary structure of human chromogranin A and pancreastatin. J Biol Chem 1987;262:17026-30.
- Murray SS, Deaven LL, Burton DW, O'Connor DT, Mellon PL, Deftos LJ. The gene for human chromogranin A (CgA) is located on chromosome 14. Biochem Biophys Res Comm 1987;42:141-6.
- Lloyd RV, Cano M, Rosa P, Hille A, Huttner WB. Distribution of chromogranin A and secretogranin I (chromogranin B) in neuroendocrine cells and tumors. Am J Pathol 1988;130:296-304.
- Degorce F. Assessment of chromogranin A using two-site immunoassay. Selection of a monoclonal antibody pair unaffected by human chromogranin A processing. Adv Exp Med Biol 2000;482:339-50.

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